




ELECTROFORCE BIODYNAMIC TEST INSTRUMENTS



Welcome to transformational biomedical technology, engineered to advance tissue growth

ElectroForce BioDynamic Series	4
BioDynamic Series Technology	6
Perfusion Bioreactors	8
•3DCulturePro Bioreactor	8
•BioDynamic Bioreactor	10
Single Specimen Mechanical Stimulation Bioreactors	12
•Single Axis	12
•Multi-Axis	14
Multi-Specimen Mechanical Stimulation Bioreactors	16
Other Bioreactor Solutions	18
Accessories	20
Software	21
Support	22
BioDynamic Flow Loop Diagrams	24
Mechanical Stimulation Bioreactor Specifications	26

Lowest Controllable Force

Multi-Specimen

Made for an Incubator

The **MOST VERSATILE**
bioreactor system
outside of the body

Load any Tissue

Widest Range of Control

Autoclavable

Design simplicity provides

unmatched
performance
& **reproducible**
test results

ElectroForce BioDynamic test instruments provide long-term tissue engineering solutions within a sterile cell culture environment. With a full range of capabilities, choose the configuration that is right for you. From introductory bioreactors to the most versatile mechanical stimulation bioreactors on the market.

Apply physiologically-relevant loading conditions

Transfer mechanical stimuli precisely while perfusing media and maintaining sterility

Engineered for reliability and durability

Designed for long-term use and experimental repeatability

Simultaneously stimulate and characterize any type of tissue

Widest range of force and displacement control to meet the increasing mechanical needs of your tissue


ElectroForce BioDynamic Series

Blood Vessels
Ligaments
Bone
Cartilage
Biomaterials
Stem Cells
Scaffolds
Biocompatible
Mechanotransduction
Regeneration
Differentiation
Proliferation
Migration
Alignment

Perfusion
Bioreactors

Single
Specimen
Mechanical
Stimulation
Bioreactors


Multi-Specimen
Mechanical
Stimulation
Bioreactors



3DCulturePro
6 Chamber System



BioDynamic Pulsatile
Test Instrument



BioDynamic 5110
Axial Test Instrument



BioDynamic 5270
Axial/Pulsatile Test Instrument



Industry **LEADER** in biomedical testing applications **for over 20 years**

Superior chamber architecture

ElectroForce bioreactors are engineered so tests can be set up quickly and easily, while still giving users the greatest versatility to run limitless experiments.

Scalable configurations meet growing research needs

Chambers are optimized to conserve space while allowing for the addition of multi-specimen fixtures, multiple bioreactors, and added loading capability.

Unmatched waveform control and fidelity


Patented ElectroForce frictionless motor design facilitates precise force, displacement and pressure control with unrivaled responsiveness.

Explicitly engineered for 3D cell culture in an incubator

ElectroForce maintenance-free motors withstand years of use in the challenging conditions of environmental chambers.


The industry's only 10 year motor warranty

We give you confidence that your system will continue to perform as your research evolves.




Introducing a
complete portfolio
of tissue culture solutions

- **Perfusion Bioreactors**
- **Single Specimen Mechanical Stimulation Bioreactors**
- **Multi-Specimen Mechanical Stimulation Bioreactors**



3D perfusion culture made **SIMPLE**




Set up a test in minutes

- Easy to use, tool-less chamber design
- Integrated media reservoir
- Portable and compact

Adaptable design

- Includes fixtures for a variety of sample types
- Can be positioned in 3 orientations
- Accommodates multiple flow loops

Keep your cells
viable during a multi-month
tissue culture experiment
so you are measuring
**the cellular response that
actually matters***



6 Week Flow Culture of Smooth Muscle Cells

* SV Biechler. 2015. Perfusion flow keeps cells viable in long-term 3D culture. TA ElectroForce Application Note ESG001.


Perfusion Bioreactors

BIODYNAMIC BIOREACTOR


All components in contact with fluid are designed to last many autoclave cycles

Transparent viewing windows allow you to image your sample during stimulation

Chamber fillers minimize fluid volume and reduce consumable costs



Apply fluid shear stress to
stimulate cellular remodeling
of the **extracellular matrix (ECM)**
and **enhance** the mechanical
and biological properties of
vascular **tissue***



*F Boccafoschi, M Basetti, C Mosca, D Mantovani, and M Cannas. 2012. The role of shear stress on mechanically stimulated engineered vascular substitutes: influence on mechanical and biological properties. Journal of Tissue Engineering and Regenerative Medicine 6(1):60-67

Single Specimen Mechanical Stimulation Bioreactors

SINGLE AXIS



BioDynamic 5110
Axial Test Instrument



BioDynamic Pulsatile Test Instrument


Optimized design gives you the freedom to focus on your research

- Perform simultaneous periodic loading experiments with supplemental chambers
- Rigid shafts are hollow to accommodate dynamic flow as well as precise loading of stiff samples
- Chamber and flow loop can be autoclaved, assembled, and sealed before transferring to the test frame

Versatile

- Loading a different sample type is as easy as changing fixtures
- Flow ports let you choose the flow loop that is right for your test

Increase airway strength
and **promote** physiological
cellular alignment
by applying **dynamic**
circumferential
strain*



*CE Ghezzi, B Marelli, MB Donelli, A Alessandrino, G Freddi, and SN Nazhat. 2014. The role of physiological mechanical cues on mesenchymal stem cell differentiation in an airway tract-like dense collagen-silk fibroin construct. *Biomaterials* 35(24):6236-6247.

Single Specimen Mechanical Stimulation Bioreactors

MULTI-AXIS



BioDynamic 5175
Axial/Torsion/Pulsatile
Test Instrument*



* Axial/Torsion and Axial/Pulsatile configurations also available

The **ONLY**
tri-axis stimulation bioreactor system
on the market.



Multi-Specimen Mechanical Stimulation Bioreactors

Patented flexure design ensures equivalent loading to all samples

- Controlled deformation of all samples with a shared actuator
- Each chamber has its own force sensor for independent sample characterization

A high throughput, yet highly flexible solution

- Multiple chambers enable reproducible results with statistical relevance
- Individual or shared flow loop options so you can customize the number of experimental variables




BioDynamic 5270*
4 Chamber Axial/Pulsatile Test Instrument



* Axial or Pulsatile configurations also available

Apply periodic **tensile** loading to **increase** expression of tenogenic **differentiation** markers in tissue engineered tendon*



* JG Barber, AM Handorf, TJ Allee, and WJ Li. 2013. Braided nanofibrous scaffold for tendon and ligament tissue engineering. Tissue Engineering Part A 19(11-12):1265-1274.

Need a system that goes above and beyond tissue engineering?

Add BioDynamic chambers to standard load frame products and experience the full breadth of ElectroForce testing capabilities.




**ElectroForce 3200 Test Instrument
with BioDynamic Bioreactor**

**ElectroForce 5500
Test Instrument
with Off-Axis
Pulsatile Bellows**



**ElectroForce Planar Biaxial Test Instrument with Planar
Biaxial BioDynamic Bioreactor**



Apply cyclic **compression**
to promote cartilage
regeneration *



**ElectroForce 5500 Test Instrument
with 24-Well Plate Fixture**

*N Peake, N Su, M Ramachandran, P Achan, DM Salter, DL Bader, AJ Moyes, AJ Hobbs, TT Chowdhury. 2013. Natriuretic peptide receptors regulate cytoprotective effects in a human ex vivo 3D/bioreactor model. Arthritis Research and Therapy 15(4):R76.



4-Well Platens

ElectroForce® BioDynamic test instruments can be outfitted with a variety of specimen fixtures, upgrade options and software modules to make your test yield the most biologically-relevant results.

Grips and Platens

- DMA Grips
- Tissue Grips
- BioDynamic® Tensile Grips
- BioDynamic Compression Platens
 - Porous (40 µm and 100 µm)
 - Nonporous
- Porous Membrane Platens
- 4-Well Nonporous Platens

Fixtures and Chambers

- Barbed Fittings
- 3 and 4 Point Bend Fixture
- Multi-Specimen Fixtures
- Additional Bioreactor Chambers
- MRI-Compatible BioDynamic Chambers

Sensors

- Force/Torque
- Displacement/Rotation
- Strain
- Pressure
- Digital Video Extensometer
- Laser Micrometer



MRI-Compatible BioDynamic Chamber

The Most Flexible Control System Available

Single comprehensive package provides an intuitive user interface, closed-loop waveform controls, and data acquisition:

- Powerful waveform generation tools to quickly create standard waveforms for basic stimulation, complex waveforms with block grouping, or user-imported non-standard waveforms
- Integrated data acquisition algorithms so a variety of data collection techniques can be utilized
- Advanced controls including multi-channel synchronization of phase and amplitude, and cross-channel compensation
- Calculated channels to provide real-time mathematical calculations of measured values
- Additional options include:
 - External Waveform Input
 - Dynamic Link Libraries
 - Dynamic Mechanical Analysis




Industry-Leading Sales & Support

TA Instruments' leadership position results from the fact that we offer the best overall product in terms of technology, performance, quality, and customer support. While each is important, our demonstrated commitment to after-sales support is a primary reason for the continued loyalty of our customers. To provide this level of support, TA Instruments has assembled the largest worldwide team of field technical and service professionals in the industry. Others promise good service. Talk to our customers and learn how TA Instruments consistently delivers on our promise to provide exceptional service.

With direct support staff in **23 countries** and **5 continents**, TA Instruments can extend its exceptional support to you, wherever you are.




BioDynamic Flow Loop Diagrams




3DCulturePro

Flow around and through tubular sample




BioDynamic 5110/5115* /5210

Flow around tensile sample



BioDynamic 5110/5115* /5210

Flow through compressive sample



BioDynamic 5170/5175* /5270

Flow through pulsatile sample

Mechanical Stimulation Bioreactor Specifications

	BioDynamic Pulsatile Test Instrument	BioDynamic 5110/ BioDynamic 5115*	BioDynamic 5170/ BioDynamic 5175*
No. of Chambers	1	1	1
Max. Force	•	± 200 N	± 200 N
Max. Displacement	•	± 6.5 mm	± 6.5 mm
Min. Displacement Increment	•	0.001 mm	0.001 mm
Max. Frequency	•	20 Hz	20 Hz
Max. Deformation Rate	•	740 mm/s	740 mm/s
Max. Pressure	2000 mm Hg	•	2000 mm Hg
Max. Pulse Volume	8.8 mL	•	8.8 mL
Max. Frequency	5 Hz	•	5 Hz
Max. Torque	•	± 2.8 Nm (*BioDynamic 5115)	± 2.8 Nm (*BioDynamic 5115)
Max. Rotation	•	± 30° (*BioDynamic 5115)	± 30° (*BioDynamic 5115)
Max. Frequency	•	20 Hz (*BioDynamic 5115)	20 Hz (*BioDynamic 5115)
Pump Type	Gear	Peristaltic	Gear
Flow Range	17-1760 mL/min	0.1 – 280 mL/min	17-1760 mL/min

- Not Available
- Upgrade Available

	BioDynamic 4 Chamber Pulsatile Test Instrument	BioDynamic 5210 Test Instrument	BioDynamic 5270 Test Instrument
No. of Chambers	4	4	4
Max. Force	•	± 200 N (± 50 per Chamber)	± 200 N (± 50 per Chamber)
Max. Displacement	•	± 6.5 mm	± 6.5 mm
Min. Displacement Increment	•	0.001 mm	0.001 mm
Max. Frequency	•	15 Hz	15 Hz
Max. Deformation Rate	•	350 mm/s	350 mm/s
Max. Pressure	2000 mm Hg (500 mm Hg per Chamber)	•	2000 mm Hg (500 mm Hg per Chamber)
Max. Pulse Volume	6.0 mL	•	6.0 mL
Max. Frequency	5 Hz	•	5 Hz
Max. Torque	—	—	—
Max. Rotation	—	—	—
Max. Frequency	—	—	—
Pump Type	Gear	Peristaltic	Gear
Flow Range	17-1760 mL/min	0.36 – 36 mL/min	17-1760 mL/min



AMERICAS

- New Castle, DE USA
- Lindon, UT USA
- Saugus, MA USA
- Eden Prairie, MN USA
- Chicago, IL USA
- Montreal, Canada
- Toronto, Canada
- Mexico City, Mexico
- São Paulo, Brazil

EUROPE

- Hüllhorst, Germany
- Eschborn, Germany
- Wetzlar, Germany
- Elstree, United Kingdom
- Brussels, Belgium
- Etten-Leur, Netherlands
- Paris, France
- Barcelona, Spain
- Milano, Italy
- Warsaw, Poland
- Prague, Czech Republic
- Sollentuna, Sweden
- Copenhagen, Denmark

ASIA & AUSTRALIA

- Shanghai, China
- Beijing, China
- Tokyo, Japan
- Seoul, South Korea
- Taipei, Taiwan
- Guangzhou, China
- Petaling Jaya, Malaysia
- Singapore
- Bangalore, India
- Sydney, Australia



tainstruments.com