



Systemic Rheology

revised by A. Franck, TA Instruments Germany

keywords: systemic rheology, transient/steady mechanical testing, dynamic mechanical testing

THE RHEO-REACTOR

Systemic rheology is the determination of rheological behavior under near real process conditions

INTRODUCTION

The elaboration of complex fluids with special properties and functionality is based on formulation engineering, a multi-level science, encompassing process engineering, physico-chemistry, chemistry, rheology and interface sciences among others. Some objectives of formulation engineering are the quantitative analysis of the different steps of the manufacturing process, the mastering of those steps in order to relate to formulations established in the factory laboratory, the knowledge and monitoring of the rheology of the complex fluids in the course of its formulation.

An original tool integrating the knowledge of process engineering and rheology and which has

been named “systemic rheology” has developed by Gemico in France /1/. To this effect, traditional geometries of the rheometers have been changed into stirred micro-reactors so that within certain steps of a process could be reproduced to some extent if not copied exactly. The test assembly constitutes the rheo-reactor operating in batch or semi-batch mode and which can perform mixing procedures, rheological monitoring in situ, i.e. without sampling and even rheological characterisations which are often difficult to realise with complex fluids. Some of these are actually extremely sensitive to their thermo-mechanical history and therefore cannot be tested in an application relevant manner outside of the process itself. Also complex materials may be of heterogeneous nature with a micro-structure of the same size as the gap of traditional geometries (cone-plate, parallel plates, Couette etc.). Conditioning and analysis of this kind of samples often meet with major experimental difficulties such as structural alterations of the material, slipping, sedimentation or characterization during an uncontrolled pre-strain state.

The configuration of the rheo-reactor itself eliminates all these problems and at the same time it fulfills some of the objectives of formulation engineering /2/. The purpose of this paper is to show that it is possible to characterise a complex fluid from a rheological point of view within the configuration of the rheo-reactor.

METHOD

Working with the rheo-reactor is based on two steps: first, the traditional geometries are replaced by stirred micro-vases and, secondly, the assembly vase/stirring rod must be gauged to allow the




Figure 1. : Helicoidal ribbon and anchor




Figure 2.: Rate sweep

determination of the characteristic material functions of the respective material; the aim of the rheometer thus remains preserved. The first step is quite simply and consists of using a universal test fixture holder and replacing the test geometry by adapting a vase and stirring rod to the end of the axes (Fig. 1). This procedure can be performed on all rotational rheometers such as ARES, AR2000 and their predecessors RDAII, RFSII, AR1000.

As to the second step, advantage is taken from the option “special geometry” available in the TAOrchestrator software. Two constants were




Figure 3.: Frequency sweep

introduced namely $(K\tau)$, which relates the torque M and the stress τ and $(K\gamma)$, relating the angular deformation θ and the deformation γ :

$$\tau = K\tau M \text{ and } \gamma = K\gamma \theta$$

These two constants can be determined by a gauging procedure allowing to extract the material

functions from instrument data [3] obtained with the special geometry. It consists in calculating the stresses and the deformations in a Couette, equivalent to a stirring system.

RESULTS

The graphs in fig. 2 to 4 show the compared test results measured in steady and oscillatory modes. The following geometries have been used:

- double wall Couette (DC), (\varnothing Bob = 47 mm, \varnothing Cup = 50 mm),
- helicoidal ribbon (HR), (\varnothing Ribbon = 31 mm, step = 30 mm, \varnothing Cup = 40 mm), and
- ancre (AN), (\varnothing Ancre 47 mm, \varnothing Cup 50mm).

The sample used was a preparation for salad dressing (Maille). All measurements were con-




Figure 4.: Strain sweep

ducted at room temperature. The repeatability of the tests is comparable to that of measurements with conventional geometries.

Figure 2 shows the flow curves obtained from tests performed with the three geometries, DWC, HR, AC. The data obtained from the DWC cover a wider range than the data obtained from HR and AC. Good agreement between all three geometries exists after shifting the HR and AC data to superpose the DWC data.

Figure 3 shows dynamic frequency data. HR and AC deviate slightly at higher frequency. Complex flow behaviour around the the ribbon and the ancre as well as inertia effects are probably the reason for these deviations.

Figure 4 shows a dynamic strain sweep. In the linear range, very good agreement can be seen between the three geometries. As expected from Figure 1, the sensitivity of the HR and the AC is slightly inferior to the Couette geometry at low strains. Minor differences can also be seen at the on-set on non-linear behavior. This is due to the difference in the flow pattern (simple flow in a Couette, complex flow around the ribbon and ancre), which causes the sample to behave differently. These differences, when more significant, can be the reason for the failure of using conventional rheometer data to model real processes.

REFERENCES

1. L. Choplin and P. Marchal, Proceedings XIIth Int. Congr. on Rheology, 665-666 (1996)
2. L. Choplin and P. Marchal, Proceedings Ist Int. Symp. on Food Rheology , 40-44 (1997)
3. A. Ait-Kadi, P. Marchal, A.S. Chrissemant, M. Bousmina and L. Choplin, American Society of Mechanical Engineers, Dallas (Texas - USA), Nov. 16-21, S.G. , FED Vol. 243, MD Vol 78, 159-170 (1997)