

RHEOLOGY SOLUTIONS

EVALUATION OF DENTURE ADHESIVES

PROBLEM

Suppliers of denture adhesives need a rapid method for comparing different product formulations for a variety of properties including their ability to coat a complex surface such as the denture cavity.

SOLUTION

Controlled stress rheology provides several measurements which can be used to compare these materials. Figure 1, for example, shows the creep results for three different formulations. Creep measurements are based on the application of a small constant stress to the material and monitoring the resultant strain with time. If the applied stress is sufficiently small, the material deforms within its linear viscoelastic limits, and the shape of the strain curve indicates the relative viscous/elastic behavior. In this case, adhesives B and C are more compliant than adhesive A and

hence should deform more easily to fill the denture cavity. Subsequent test panels confirmed this projection with B and C providing acceptable denture seals while A was unacceptable, indicating that creep curves, once a threshold compliance is established with suitable known materials, can be used to rapidly assess coating acceptability of new formulations.

All three materials show good recovery in the relaxation made after stress is removed, implying that the pastes will “bounce back” from applied stress either during dispensing from a container or during placement in the mouth.

If it is desirable to quantitatively compare formulations and model their behavior, this can also be accomplished from these creep curves. Figure 2 shows the quantitative parameters determined for Adhesive C when fitted to a simple Voigt model using the *Rheology Solutions* software.




Figure 2. ANALYSIS RESULTS (Voigt Model)

Analysis Results - Voigt

File name and step

ADHESIVE C

compliance	3.1580 X 10 ⁻⁴ m ² /N
newtonian viscosity	515900 Pa·s
newtonian shear rate	2.907 X 10 ⁻⁵ 1/s
compliance	5.0051X 10 ⁻⁴ m ² /N
1st voigt viscosity	84290 Pa·s
1st voigt time	42.19 s
compliance	2.6296 X 10 ⁻⁴ m ² /N
2nd voigt viscosity	17160 Pa·s
2nd voigt time	4.512 s
compliance	0 m ² /N
3rd voigt viscosity	0 Pa·s
3rd voigt time	0 s
compliance	0 m ² /N
4th voigt viscosity	0 Pa·s
4th voigt time	0 s
standard error	6.582
Batch	2134.0

OK

ACKNOWLEDGMENT: This brief is based on work by P. Whittingstall, TA Instruments Applications Laboratory (US).

For more information or to place an order, contact:

TA Instruments, Inc., 109 Lukens Drive, New Castle, DE 19720, Telephone: (302) 427-4000, Fax: (302) 427-4001

TA Instruments S.A.R.L., Paris, France, Telephone: 33-01-30489460, Fax: 33-01-30489451

TA Instruments N.V./S.A., Gent, Belgium, Telephone: 32-9-220-79-89, Fax: 32-9-220-83-21

TA Instruments GmbH, Alzenau, Germany, Telephone: 49-6023-30044, Fax: 49-6023-30823

TA Instruments, Ltd., Leatherhead, England, Telephone: 44-1-372-360363, Fax: 44-1-372-360135

TA Instruments Japan K.K., Tokyo, Japan, Telephone: 813-5434-2771, Fax: 813-5434-2770

Internet: <http://www.tainst.com>

TA Instruments
Thermal Analysis & Rheology
A SUBSIDIARY OF WATERS CORPORATION