



## Analysis of a Gelatin Capsule Using Controlled Humidity Dynamic Mechanical Analysis

### ABSTRACT

This paper discusses the dynamic mechanical analysis of a gelatin capsule material under varying conditions of temperature and relative humidity.

### INTRODUCTION

Gelatin is a naturally occurring material which is used in a wide array of applications, most typically in edible/foods, pharmaceutical, photographic and technical products. The use of gelatin in the manufacture of various pharmaceutical dosage forms dates back to the early 19th century and possibly earlier<sup>1</sup>. When stored in an ambient, low-humidity environment gelatin is remarkably stable. However, when combined with water, gelatin forms a semi-solid colloid gel and this can profoundly affect its mechanical properties. Particularly as it relates to storage and handling conditions, it is important to understand and quantify these effects *in situ*.

The TA Instruments DMA-RH Accessory allows the mechanical properties of a sample to be analyzed under constant and/or varying conditions of both relative humidity and temperature. It is designed for use with the Q800 Dynamic Mechanical Analyzer. The DMA-RH accessory is an integrated unit and contains the following components:



**Figure 1: The TA Instruments Q800 Dynamic Mechanical Analyzer and DMA-RH Accessory**

1. The sample chamber mounts to the DMA in place of the standard furnace and encloses the sample. Peltier elements in the chamber precisely control the temperature to within  $\pm 0.1^\circ\text{C}$ . The sample chamber accommodates standard DMA clamps including tension,

cantilever, and 3-point bending, and can be easily removed for rapid conversion back to the standard DMA furnace.

2. The DMA-RH Accessory contains the humidifier and electronics which continuously monitor and control temperature and humidity of the sample chamber. The DMA Q800 and the DMA-RH Accessory are fully software-integrated.
3. A heated vapor transfer line is maintained above the dew point temperature of the humidified gas in order to avoid condensation and provide accurate results.




Figure 2: Sample Chamber of the DMA-RH Accessory


The DMA-RH accessory allows for the control of temperature over the range 5-120°C, and humidity over the range 5-95% RH. As such, it is well-suited to investigate gelatin and other polymeric materials which exhibit structural and mechanical instability within this temperature range.

## RESULTS & DISCUSSION


A sample of dimension 9.2 mm x 3.2 mm x 0.12 mm was cut from the sidewall of a two-piece hard gelatin capsule. It is important to minimize the thickness of the sample to facilitate efficient water transfer throughout the polymer matrix. The sample was analyzed at a frequency of 1 Hz and amplitude of 20  $\mu\text{m}$  (~0.2% strain) over the humidity range 5-95% RH at a rate<sup>i</sup> of 0.25% RH/min. The experiment was performed at both 25°C and 40°C. The resultant storage modulus and tan delta are shown in Figures 3 and 4 on the next page. The data demonstrate that the gelatin material undergoes a multi-step transition as the humidity is increased. At 25°C, an initial loss in modulus is seen around 25% RH, with a more dramatic loss of structure near 80% RH. At the elevated temperature of 40°C, the initial loss remains near 25% RH, with the subsequent transitions shifting to *ca.* 70% RH. This suggests that even at ambient temperature, the storage conditions for gelatin capsules should be maintained below 70% RH to ensure mechanical integrity of the material.

---

<sup>i</sup> The RH ramp was composed of a series of short isohume steps, 0.25%RH / 1 min hold.



**Figure 3: DMA-RH Analysis of Gelatin**



**Figure 4: DMA-RH Analysis of Gelatin (Expanded Scale)**

**CONCLUSIONS**

The data presented illustrate the dramatic effect of water on the mechanical properties of gelatin. The TA Instruments Q800 DMA and DMA-RH Accessory provide the ideal platform for this investigation.

**REFERENCES**

1. [http://www.gelatin-gmia.com/html/rawmaterials\\_app.html](http://www.gelatin-gmia.com/html/rawmaterials_app.html) (Gelatin Manufacturer's Institute of America)

## **TA INSTRUMENTS**

### **United States**

109 Lukens Drive, New Castle, DE 19720 • Phone: 1-302-427-4000 • E-mail: [info@tainstruments.com](mailto:info@tainstruments.com)

### **Canada**

Phone: 1-905-309-5387 • E-mail: [shunt@tainstruments.com](mailto:shunt@tainstruments.com).

### **Mexico**

Phone: 52-55-5200-1860 • E-mail: [mdominguez@tainstruments.com](mailto:mdominguez@tainstruments.com)

### **Spain**

Phone: 34-93-600-9300 • E-mail: [spain@tainstruments.com](mailto:spain@tainstruments.com)

### **United Kingdom**

Phone: 44-1-293-658-900 • E-mail: [uk@tainstruments.com](mailto:uk@tainstruments.com)

### **Belgium/Luxembourg**

Phone: 32-2-706-0080 • E-mail: [belgium@tainstruments.com](mailto:belgium@tainstruments.com)

### **Netherlands**

Phone: 31-76-508-7270 • E-mail: [netherlands@tainstruments.com](mailto:netherlands@tainstruments.com)

### **Germany**

Phone: 49-6196-400-7060 • E-mail: [germany@tainstruments.com](mailto:germany@tainstruments.com)

### **France**

Phone: 33-1-304-89460 • E-mail: [france@tainstruments.com](mailto:france@tainstruments.com)

### **Italy**

Phone: 39-02-2742-11 • E-mail: [italia@tainstruments.com](mailto:italia@tainstruments.com)

### **Sweden/Norway**

Phone: 46-8-555-11-521 • E-mail: [sweden@tainstruments.com](mailto:sweden@tainstruments.com)

### **Japan**

Phone: 813-5479-8418 • E-mail: [j-marketing@tainstruments.com](mailto:j-marketing@tainstruments.com)

### **Australia**

Phone: 613-9553-0813 • E-mail: [sshamis@tainstruments.com](mailto:sshamis@tainstruments.com)

### **India**

Phone: 91-80-2839-8963 • E-mail: [india@tainstrument.com](mailto:india@tainstrument.com)

### **China**

Phone: 8610-8586-8899 • E-mail: [info@tainstruments.com.cn](mailto:info@tainstruments.com.cn)

### **Taiwan**

Phone: 886-2-2563-8880 • E-mail: [skuo@tainstruments.com](mailto:skuo@tainstruments.com)

### **Korea**

Phone: 82.2.3415.1500 • E-mail: [dhchee@tainstruments.com](mailto:dhchee@tainstruments.com)

To contact your local TA Instruments representative visit our website at [www.tainstruments.com](http://www.tainstruments.com)

© Copyright 2007 TA Instruments