

Thermomechanical analysis (TMA) is a technique widely used to determine transition temperatures and expansion/contraction properties of materials, including the coefficient of thermal expansion (α). When determining α , however, there are three factors which need to be considered to obtain reproducible quantitative results. These are TMA sensitivity and baseline drift, sample thickness, and the magnitude of α .

The TA Instruments TMA has a maximum sensitivity of 15 nm, a displacement resolution of 0.5 μm , and a baseline drift of 1 μm (typically, 0.5 μm between 0 to 500 $^{\circ}\text{C}$). Therefore, the recommended measured sample dimensional change over this range for a determinations is >5 μm (roughly 10x the baseline drift). Knowing this, the minimum sample thickness required for testing can be determined from:

$$\alpha = \frac{\Delta L}{\Delta T} * \frac{1}{L_o} \rightarrow L_o = \frac{\Delta L}{\Delta T} * \frac{1}{\alpha} \quad (1)$$

where:

- ΔL = measured dimension change
- L_o = original sample thickness
- ΔT = temperature range of measurement

For example, a material with an expected α of 50 $\mu\text{m}/\text{m}^{\circ}\text{C}$ would yield a recommended minimum sample thickness for testing of:

$$L_o = \left(\frac{5 \mu\text{m}}{500 \text{ }^{\circ}\text{C}} \right) * \left(\frac{1}{50 \mu\text{m}/\text{m }^{\circ}\text{C}} \right) = 200 \mu\text{m} \quad (2)$$

This minimum sample thickness will obviously change with a and temperature range. Clearly, it is better to use samples much thicker than this minimum if possible, as this will further improve signal-to-noise ratio.

ACKNOWLEDGMENTS

This applications note was submitted by Jon Foreman of the Applications Laboratory (US).

For more information or to request a product quote, please visit www.tainstruments.com/ to locate your local sales office information.