

## THERMAL APPLICATIONS NOTE

### DETERMINING MINIMUM USABLE SAMPLE THICKNESS IN TMA

Thermomechanical analysis (TMA) is a technique widely used to determine transition temperatures and expansion/contraction properties of materials, including the coefficient of thermal expansion ( $\alpha$ ). When determining  $\alpha$ , however, there are three factors which need to be considered to obtain reproducible quantitative results. These are TMA sensitivity and baseline drift, sample thickness, and the magnitude of  $\alpha$ .

The TA Instruments TMA 2940 has a maximum sensitivity of 0.1  $\mu\text{m}$ , a maximum resolution of 3nm, and a baseline drift (0 to 500°C) of about 0.5 $\mu\text{m}$ . Therefore, the recommended measured sample dimensional change over this range for  $\alpha$  determinations is >5 $\mu\text{m}$  (roughly 10x the baseline drift). Knowing this, the minimum sample thickness required for testing can be determined from:

$$\alpha = \frac{\Delta L}{\Delta T} * \frac{1}{L_o} \Rightarrow L_o = \frac{\Delta L}{\Delta T} * \frac{1}{\alpha} \quad [1]$$

where:  $\Delta L$  = measured dimension change  
 $L_o$  = original sample thickness  
 $\Delta T$  = temperature range of measurement

For example, a material with an expected  $\alpha$  of 50 $\mu\text{m}/\text{m}^\circ\text{C}$  would yield a recommended minimum sample thickness for testing of:

$$L_o = \left( \frac{5\mu\text{m}}{500^\circ\text{C}} \right) * \left( \frac{1}{50\mu\text{m} / \text{m}^\circ\text{C}} \right) = 200\mu\text{m} \quad [2]$$

This minimum sample thickness will obviously change with  $\alpha$  and temperature range. Clearly, it is better to use samples much thicker than this minimum if possible, as this will further improve signal-to-noise ratio.

Acknowledgment: This applications note was submitted by Jon Foreman of the Applications Laboratory (US).

For more information or to place an order, contact:

**TA Instruments, Inc.**, 109 Lukens Drive, New Castle, DE 19720, Telephone: (302) 427-4000, Fax: (302) 427-4001

**TA Instruments S.A.R.L.**, Paris, France, Telephone: 33-01-30489460, Fax: 33-01-30489451

**TA Instruments N.V./S.A.**, Gent, Belgium, Telephone: 32-9-220-79-89, Fax: 32-9-220-83-21

**TA Instruments GmbH**, Alzenau, Germany, Telephone: 49-6023-30044, Fax: 49-6023-30823

**TA Instruments, Ltd.**, Leatherhead, England, Telephone: 44-1-372-360363, Fax: 44-1-372-360135

**TA Instruments Japan K.K.**, Tokyo, Japan, Telephone: 813-5434-2771, Fax: 813-5434-2770

Internet: <http://www.tainst.com>