

## THERMAL APPLICATIONS NOTE

### Sapphire Specific Heat Capacity Literature Values

Specific heat capacity measurements by DSC require the use of a well-characterized reference material (usually sapphire) to obtain results. The table below summarizes specific heat capacity data (literature values) for sapphire as a function of temperature. These values can be used in the required calculations.

Temperature °C	K	Specific Heat J/g°C	Temperature °C	K	Specific Heat J/g°C
-183.15	90	0.0949	126.85	400	0.9423
-173.15	100	0.1261	136.85	410	0.9545
-163.15	110	0.1603	146.85	420	0.9660
-153.15	120	0.1968	156.85	430	0.9770
-143.15	130	0.2349	166.85	440	0.9875
-133.15	140	0.2739	176.85	450	0.9975
-123.15	150	0.3134	186.85	460	1.0070
-113.15	170	0.3526	196.85	470	1.0161
-103.15	170	0.3913	206.85	480	1.0247
-93.15	180	0.4291	216.85	490	1.0330
-83.15	190	0.4659	226.85	500	1.0409
-73.15	200	0.5014	236.85	510	1.0484
-63.15	210	0.5356	246.85	520	1.0557
-53.15	220	0.5684	256.85	530	1.0627
-43.15	230	0.5996	266.85	540	1.0692
-33.15	240	0.6294	276.85	550	1.0756
-23.15	250	0.6579	286.85	560	1.0817
-13.15	260	0.6848	296.85	570	1.0876
-3.15	270	0.7103	306.85	580	1.0932
0.00	273.15	0.7180	316.85	590	1.0987
6.85	280	0.7343	326.85	600	1.1038
16.85	290	0.7572	336.85	610	1.1089
26.85	300	0.7788	346.85	620	1.1137
36.85	310	0.7994	356.85	630	1.1183
46.85	320	0.8188	366.85	640	1.1228
56.85	330	0.8373	376.85	650	1.1271
66.85	340	0.8548	386.85	660	1.1313
76.85	350	0.8713	396.85	670	1.1353
86.85	360	0.8871	406.85	680	1.1393
96.85	370	0.9020	416.85	690	1.1431
106.85	380	0.9161	426.85	700	1.1467
116.85	390	0.9296	446.85	720	1.1538

Temperature °C	K	Specific Heat J/g°C	Temperature °C	K	Specific Heat J/g°C
466.85	740	1.1604	806.85	1080	1.2383
486.85	760	1.1667	826.85	1100	1.2417
506.85	780	1.1726	846.85	1120	1.2451
526.85	800	1.1783	866.85	1140	1.2484
546.85	820	1.1837	886.85	1160	1.2516
566.85	840	1.1888	906.85	1180	1.2548
586.85	860	1.1937	926.85	1200	1.2578
606.85	880	1.1985	976.85	1250	1.2653
626.85	900	1.2030	1026.85	1300	1.2724
646.85	920	1.2074	1076.85	1350	1.2792
666.85	940	1.2117	1126.85	1400	1.2856
686.85	960	1.2159	1176.85	1450	1.2917
706.85	980	1.2198	1226.85	1500	1.2975
726.85	1000	1.2237	1276.85	1550	1.3028
746.85	1020	1.2275	1326.85	1600	1.3079
766.85	1040	1.2312	1376.85	1650	1.3128
786.85	1060	1.2348			

Table taken from D. A. Ditmars, et. al., *J. Res. Nat. Bur. Stand.*, 87, (2), 159-163 (1982).

The values in the table were determined by Ginnings and Furukawa of the National Bureau of Standards on aluminum oxide in the form of synthetic sapphire (corundum). The sapphire pieces passed a #10 sieve but were retained by a #40 sieve, and had 99.98 to 99.99 % purity by weight. Heat capacity values below the experimental range were obtained by extrapolation of the Debye equation fitted to the experimental value at the lowest temperature. The units J/g-mole are absolute joules per degree per gram mole (molecular weight - 101.9613) at a constant pressure of 1 atmosphere (100 Pa).

---

For more information or to place an order, contact:

**TA Instruments, Inc.**  
109 Lukens Drive  
New Castle, DE 19720  
Telephone: (302) 427-4000  
Fax: (302) 427-4001

**TA Instruments N.V./S.A.**  
Gent, Belgium  
Telephone: 32-9-220-79-89  
Fax: 32-9-220-83-21

**TA Instruments, Ltd.**  
Leatherhead, England  
Telephone: 44-1-372-360363  
Fax: 44-1-372-360135

**TA Instruments S.A.R.L.**  
Paris, France  
Telephone: 33-01-30489460  
Fax: 33-01-30489451

**TA Instruments GmbH**  
Alzenau, Germany  
Telephone: 49-6023-30044  
Fax: 49-6023-30823

**TA Instruments Japan K.K.**  
Tokyo, Japan  
Telephone: 813-5434-2771  
Fax: 813-5434-2770

---

Internet: <http://www.tainst.com>

TN-8A

**TA Instruments**  
Thermal Analysis & Rheology  
A SUBSIDIARY OF WATERS CORPORATION