

THERMAL SOLUTIONS

STORAGE EFFECTS ON T_g FOR EPOXY MOLDING COMPOUND*

PROBLEM

A manufacturer of epoxy molding compound received complaints from one of its customers that the batches (drums) of material supplied were pre-cured to different levels. This conclusion was based on the fact that the glass transition temperature was lower for some batches than others. Since all batches (T_g) were prepared using the same process, the differences seen could not be explained by the manufacturer.

SOLUTION

By carefully monitoring the procedure the customer followed when using the epoxy, the manufacturer noticed that drums of the material were not closed tightly in cases where only a partial batch was initially used. Evaluation of the epoxy

after varying exposure to air by TMA (Figure 1) subsequently showed that the glass transition shifted to lower temperature with longer exposure. This shift, which was as much as 8°C after 93 hours of exposure, was traced to moisture in the air. The effects of moisture were not reversible since post exposure baking did not restore a higher T_g.

TMA was ideal for evaluating these materials because of its high sensitivity to the T_g based on the associated change in expansion coefficient. Although several thermal analysis techniques can be used to detect the T_g, the relative sensitivity of each technique depends on the amount of change in the measured parameter at the T_g. Table 1 summarizes these factors.




TABLE 1 COMPARISON OF THERMAL TECHNIQUES FOR T _g MEASUREMENT		
TECHNIQUE	MEASURED PARAMETER	RELATIVE SENSITIVITY (10 is most sensitive)
DSC	Heat Capacity	1
TMA	Coefficient of expansion	5
DMA	Modulus, damping	10
DEA	Permittivity, loss factor	10

*ACKNOWLEDGEMENT: This work was performed by W. Lee, EXTRAD Instrumentation Pte. (Singapore)

For more information or to place an order, contact:

TA Instruments, Inc., 109 Lukens Drive, New Castle, DE 19720, Telephone: (302) 427-4000, Fax: (302) 427-4001

TA Instruments S.A.R.L., Paris, France, Telephone: 33-01-30489460, Fax: 33-01-30489451

TA Instruments N.V./S.A., Gent, Belgium, Telephone: 32-9-220-79-89, Fax: 32-9-220-83-21

TA Instruments GmbH, Alzenau, Germany, Telephone: 49-6023-30044, Fax: 49-6023-30823

TA Instruments, Ltd., Leatherhead, England, Telephone: 44-1-372-360363, Fax: 44-1-372-360135

TA Instruments Japan K.K., Tokyo, Japan, Telephone: 813-5434-2771, Fax: 813-5434-2770

Internet: <http://www.tainst.com>