

Thermal Analysis & Rheology

THERMAL SOLUTIONS

Determination of the Dimensional Stability of a Thin PET Film

PROBLEM

Polyethylene Terephthalate (PET) is often used as a thin film substrate for magnetic recording tapes. These PET films are prepared by rapid extrusion so that high biaxial orientation and high dimensional stability are obtained. Processors of these films are interested in quickly evaluating their dimensional stability (expansion) properties as well as the presence of a glass transition. The latter is undesirable in the temperature range of usage because the dimensional stability usually decreases substantially at the glass transition temperature (Tg.)

SOLUTION

Figure 1 shows the thermomechanical analysis (TMA) for three commercial tape materials. TMA is an ideal technique

for evaluating thin films since it directly measures dimensional change with temperature. Furthermore, since the coefficient of expansion changes dramatically at the Tg, TMA is more sensitive for detecting this and related transitions than differential scanning calorimetry (DSC). The sensitivity of TMA is particularly valuable in oriented films which usually have high crystallinity which further limits detection by techniques such as DSC.

In this study, the results obtained indicated that the Tg's for all three films evaluated were above the highest temperature of actual end-use and that the dimensional stabilities up to the Tg were acceptable for tape applications. Above Tg, tape C had the highest dimensional stability. These TMA tests were run using small samples (4-5 mm) in about one hour, making this an easy method for evaluating the effects of changes in process conditions.



For more information or to place an order, contact:

TA Instruments, Inc., 109 Lukens Drive, New Castle, DE 19720, Telephone: (302) 427-4000, Fax: (302) 427-4001
TA Instruments S.A.R.L., Paris, France, Telephone: 33-01-30489460, Fax: 33-01-30489451
TA Instruments N.V./S.A., Gent, Belgium, Telephone: 32-9-220-79-89, Fax: 32-9-220-83-21
TA Instruments GmbH, Alzenau, Germany, Telephone: 49-6023-30044, Fax: 49-6023-30823
TA Instruments, Ltd., Leatherhead, England, Telephone: 44-1-372-360363, Fax: 44-1-372-360135
TA Instruments Japan K.K., Tokyo, Japan, Telephone: 813-5434-2771, Fax: 813-5434-2770

Internet: http://www.tainst.com