




# **Thermal Properties Calculations Using Finite-Difference Method**

Akhan Tleoubaev and  
Andrzej Brzezinski

**LASERCOMP, Inc., Saugus, Massachusetts**

Presented at the 5<sup>th</sup> Annual Vacuum Insulation Association  
Symposium (VIA-2002) May 22-23, 2002, Atlanta, Georgia

# FOX Heat Flow Meter Instruments ASTM Standard C 518



- Two isothermal plates - hot and cold - at temperatures  $T_H$  and  $T_C$
- Heat Flow Meters:  $Q_H$  ,  $Q_C$
- Flat sample of thickness  $DX$
- Thermal conductivity  $\lambda$
- Thermal diffusivity  $a$  - now can be determined as well

# Regular Formula for the Heat Flow Meter Method

Steady-State - after reaching full thermal equilibrium  
- when Fourier number  $Fo = at/DX^2 \gg 1$

- $\lambda = DX/(T_H - T_C) * [S_H(T_H) * Q_H + S_C(T_C) * Q_C] / 2$
- average of the two heat flow meters
- $S_H$  and  $S_C$  are calibration factors

# Thermal Conductivity Equation (Non-Steady Transient State)

$$\partial^2 T(x,t) / \partial x^2 = (1/a) \partial T(x,t) / \partial t$$

- $a$  - thermal diffusivity [ $\text{m}^2/\text{s}$ ]
- Boundary conditions (B.C.):
- For hot plate ( $x=0$ )  $T(0,t)=T_H(t)$
- For cold plate ( $x=DX$ )  $T(DX,t)=T_C(t)$
- Initial conditions (I.C.) ( $t=0, 0 < x < DX$ ):  
 $T(x,0)$

# One-dimensional temperature field evolution with time: $T(x,t)$

- Initially, sample has some initial temperature distribution  $T(x,0)$  (usually, room temperature).
- After placing flat sample between the instrument's plates, its inner temperature distribution  $T(x,t)$  starts to change according to the thermal conductivity equation, eventually reaching the final thermal equilibrium condition.
- This transient process contains information about both thermal conductivity and diffusivity.

# Thermal Conductivity Equation using Finite-Difference Method:

- $[T(x+\delta x, t) - 2T(x, t) + T(x-\delta x, t)] / (\delta x)^2 \cong \cong (1/a) [\underline{T(x, t+\delta t)} - T(x, t)] / \delta t$
- Next moment temperatures  $T(x, t+\delta t)$  can be calculated using previous moment temperatures  $T(x+\delta x, t)$ ,  $T(x, t)$ , and  $T(x-\delta x, t)$
- Boundary conditions:
- $T(0, t) = T_H(t); \quad T(DX, t) = T_C(t)$
- $q_H(t) = \lambda [T(0, t) - T(\delta x, t)] / \delta x = S_H * Q_H(t)$
- $q_C(t) = \lambda [T(DX, t) - T(DX - \delta x, t)] / \delta x = S_C * Q_C(t)$

# Two pairs of the heat flow arrays:

- 1) Experimental:  $Q_H(t)$  and  $Q_C(t)$
- 2) Calculated:  $Q_{Hcalc}(t)$  and  $Q_{Ccalc}(t)$  which are calculated using two input parameters - values of thermal conductivity  $\lambda$  and thermal diffusivity  $\mathbf{a}$ . Their correct values can be found using Least-Squares Method (weighted):
- $F(\lambda, \mathbf{a}) = \Sigma \{ [Q_H - Q_{Hcalc}]^2 + [Q_C - Q_{Ccalc}]^2 \} m^2$   
sum of squares of differences  $\Rightarrow$  minimum

# $F(\lambda, a) \Rightarrow \text{minimum}$


- Any analytical function has minimum if all (two in our case) its partial derivatives are equal to zero:
- $f_{\lambda} = \partial F(\lambda, a) / \partial \lambda = 0$
- $f_a = \partial F(\lambda, a) / \partial a = 0$
- Newton's method is most useful and convenient to find the best pair of  $\lambda$  and  $a$  values using their initial guess values and iteration calculation procedure.



# Newton's method:

- $f_{\lambda} + (\partial f_{\lambda} / \partial \lambda)[\lambda^{(j+1)} - \lambda^{(j)}] + (\partial f_{\lambda} / \partial a)[a^{(j+1)} - a^{(j)}] = 0$
- $f_a + (\partial f_a / \partial \lambda)[\lambda^{(j+1)} - \lambda^{(j)}] + (\partial f_a / \partial a)[a^{(j+1)} - a^{(j)}] = 0$
- System of two linear equations with two unknowns:  $\lambda^{(j+1)}$  and  $a^{(j+1)}$ :
  - 1)  $(\partial f_{\lambda} / \partial \lambda)\lambda^{(j+1)} + (\partial f_{\lambda} / \partial a)a^{(j+1)} =$   
 $= (\partial f_{\lambda} / \partial \lambda)\lambda^{(j)} + (\partial f_{\lambda} / \partial a)a^{(j)} - f_{\lambda};$
  - 2)  $(\partial f_a / \partial \lambda)\lambda^{(j+1)} + (\partial f_a / \partial a)a^{(j+1)} =$   
 $= (\partial f_a / \partial \lambda)\lambda^{(j)} + (\partial f_a / \partial a)a^{(j)} - f_a;$

# Heat Flow Meters Signals - calculated and experimental (in microvolts) versus time (in seconds) - 1"-thick "good" vacuum panel




Thermal conductivity  
 $\lambda = 5.79 \text{ mW}/(\text{m} \cdot \text{K})$

Thermal diffusivity  
 $a = 4.3 \cdot 10^{-8} \text{ m}^2/\text{s}$

Volumetric specific heat  
 $C_p \rho = 135 \text{ kJ}/(\text{m}^3 \cdot \text{K})$

Thermal effusivity  
 $\varepsilon = 27.9 \text{ W sec}^{1/2}/(\text{m}^2 \cdot \text{K})$

# Heat Flow Meters Signals - calculated and experimental (in microvolts) versus time (in seconds) - 1"-thick "bad" vacuum panel




Thermal conductivity  
 $\lambda = 33.3 \text{ mW}/(\text{m} \cdot \text{K})$

Thermal diffusivity  
 $a = 2.9 \cdot 10^{-7} \text{ m}^2/\text{s}$

Volumetric specific heat  
 $C_p \rho = 115 \text{ kJ}/(\text{m}^3 \cdot \text{K})$

Thermal effusivity  
 $\varepsilon = 61.8 \text{ W} \cdot \text{sec}^{1/2}/(\text{m}^2 \cdot \text{K})$

# Heat Flow Meters Signals - calculated and experimental (in microvolts) versus time (in seconds) - 2"-thick "good" vacuum panel



Thermal conductivity  
 $\lambda=5.23 \text{ mW}/(\text{m}\cdot\text{K})$

Thermal diffusivity  
 $a=4.3\cdot 10^{-8} \text{ m}^2/\text{s}$

Volumetric specific heat  
 $C_p\rho=122 \text{ kJ}/(\text{m}^3\text{K})$

Thermal effusivity  $\varepsilon=25.2$   
 $\text{W sec}^{1/2}/(\text{m}^2\text{K})$

# New mathematical algorithm of calculations was developed

- Two thermal properties - thermal conductivity  $\lambda$  and thermal diffusivity  $\mathbf{a}$  can be calculated long before reaching full thermal equilibrium.
- Also two more thermal properties can be calculated - volumetric specific heat  $\mathbf{C}_p\rho = \lambda/\mathbf{a}$  and thermal effusivity  $\varepsilon = \lambda/\sqrt{\mathbf{a}}$

# Prospective

- The new algorithm will be used in the LaserComp's "WinTherm" software for FOX family of the Heat Flow Meter instruments.
- Tests duration can be made many times shorter using the new algorithm. This will be especially efficient in case of vacuum super-insulation panels and thick samples of thermal insulation materials.