

# TAM IV Calorimeter



## Site Preparation Guide

# Table of Contents


Table of Contents .....	2
Ideal Setup .....	3
System Components.....	4
Instrument Measurements .....	5
Utility Requirements.....	6-7
Laboratory .....	6
Power .....	6
Gas .....	7
Computer Requirements .....	8-9
Accessories .....	10-11
Air Dryer.....	10
Uninterrupted Power Source (UPS) .....	11
Site Preparation Checklist .....	12
TA Instrument Offices.....	13

# Ideal Setup



## IDEAL PLACEMENT AND BENCH MEASUREMENTS

Select a location with adequate floor space for the instrument and a bench for the computer (controller) and any required accessories.



Floor space width: 1–1.2 m (3.5–4 ft)

Floor space depth: 1 m (3.5 ft)

Space on each side of the instrument: 300 mm (12 in)


Bench space for computer: 40 cm (15.7 in)

Bench space for optional accessories: 0.6–1 m (2–3 ft)

# System Components



## MAIN SYSTEM COMPONENTS



- A. Computer (Controller)
- B. Instrument
- C. Uninterruptible Power Source (Optional)
- D. Gas Dryer (Optional)
- E. Accessories (Optional)



Uninterrupted Power Source (UPS) models are updated frequently by the manufacturer and may not look exactly like the one pictured.

# Instrument Measurements



## MAIN INSTRUMENT



Height: 100 cm (40 in)

Width: 58 cm (23 in)

Depth: 76 cm (30 in)

Weight: 160 kg (350 lbs)

# Utility Requirements





## LABORATORY

Item	Requirement
Temperature	15–30°C
Relative Humidity	5–80% (non-condensing)
Temperature Stability	±1°C over 24 hours with changes in temperature being gradual shifts instead of fast changes
Instrument Location Environment	<ul style="list-style-type: none"><li>• Dust-free</li><li>• Vibration-free</li><li>• Away from exposure to direct sunlight and direct air drafts</li><li>• Pollution Degree 2 Environment</li><li>• Maximum altitude: 2000 m (6560 ft)</li></ul>



## POWER

Item	Requirement
Instrument power	100–240V, 50/60 Hz, 1400 Watts
Power line	<ul style="list-style-type: none"><li>• Grounded, single-phase line for instrument and computer, not shared with motors, heaters, or compressors</li><li>• No fluctuation between ground and neutral</li></ul>
Electrical power cord	<ul style="list-style-type: none"><li>• The plug of the cord must be rated to carry at least 125% of the product current rating.</li><li>• The cord length must be less than 4.5 m (15 ft) and must be UL or CSA approved.</li><li>• 20 A socket for the UPS with 110 VAC; Standard socket for the UPS with 220 VAC</li></ul> <p> <b>Use power cords with plugs appropriate for your circuit.</b></p>
Customer-supplied	Surge suppressor power strip 

# Utility Requirements



## GAS

Item	Requirement
Conditions	<ul style="list-style-type: none"><li>• Must be dry</li><li>• Must be free from oil, water, and dirt</li></ul>
Purging	<ul style="list-style-type: none"><li>• Nitrogen required for purging the oil reservoir. For installation only; not required for normal daily use.</li><li>• Nitrogen required for RH Perfusion accessory (2 bar/29 psig).</li><li>• Dry gas (air) for purging the cabinet when operating near/below ambient dew point. 3 L/min (higher when loading samples. Compressed air or nitrogen from tanks can be an acceptable source for low temperature operation lasting a few days. For long-term experiments, use a continuously available gas supply. Refer to <a href="#">Air Dryer Requirements</a> for specifications for on-site compressed air.</li></ul>
Regulator	<ul style="list-style-type: none"><li>• Low pressure dual stage regulator that provides 0.6–3 bar (10–50 psig)</li></ul>



# Computer Requirements



## HARDWARE REQUIREMENTS


Item	Requirement
Processor	<ul style="list-style-type: none"><li>• Intel® Core™ i5 8400 or better</li><li>• 2.8 GHz with 9 MB L2 cache</li></ul>
Memory	≥ 16 GB RAM DDR4 2666 SDRAM
Hard drive	≥ 80 GB free space
DVD (optional)	≥ 48x CD-ROM or DVD (optional for installing software)
USB Ports	<ul style="list-style-type: none"><li>• 1 USB port for the TAM IV, and at least 1 other USB port for the RH Perfusion Accessory, Titration, I/O controls, and SolCal accessories.</li><li>• 2.0</li></ul>



# Computer Requirements



## SOFTWARE REQUIREMENTS

Item	Requirement
Operating System	<u>TAM Assistant</u> <ul style="list-style-type: none"><li>• Windows 10 or 11 Pro &amp; Enterprise</li><li>• Home version not supported</li><li>• 64-bit version</li></ul>
Service Pack	Microsoft Operating System Service Pack
Windows Updates	Turn off Windows power-saving settings  <b>Updates must be controlled.</b> Refer to the document “Controlling Windows Updates” on the TA Instruments website.
Network	<i>TA Instruments is not responsible for resolving issues associated with connections to your corporate network.</i>
Conflicts	<i>TA Instruments is not responsible for resolving hardware/software conflicts created by the addition of third-party hardware or software to the computer.</i>

# Accessories



## AIR DRYER REQUIREMENTS\*

\*Optional; Only when working below the dew point with the low temperature kit.

### Requirements



- Inlet air: 100–130 psig (0.7–0.9 MPa); air temperature  $\sim 20^{\circ}\text{C}$ ; Relative Humidity of 70% or less at RT with particle size of 5 microns (0.0002 in) or less
- Air source into dryer should be oil-less compressed air



- The dryer weighs 3.2 kg (7 lbs)
- It has two mounting holes 22.3 cm (8.8 in) apart
- Must be mounted upright to the wall within 183–244 cm (6–8 ft) of the air source



Customer must provide:

- Means to connect to a 3/8" NPT male connector on the inlet hose (provided by TA Instruments) of the air dryer
- Gauge to monitor the air into the air dryer
- Water trap if there is excessive moisture levels that result in immediate condensation into water (installed by Customer's maintenance personnel)



# Accessories



## UPS MEASUREMENTS



Height: 33cm (12.8 in)

Width: 21 cm (8.4 in)

Depth: 41 cm (16 in)

Weight: 34 kg (76 lbs)



UPS models are updated frequently by the manufacturer and may not look exactly like the one pictured.



## UPS REQUIREMENTS

### Requirements



#### 120V:

- Input voltage range 90–144 VAC
- Frequency 50/60 Hz
- Nominal current 16.6 A
- Power minimum 2200 VA (2.2 kVA)






#### 230V:

- Input voltage range 160-276 VAC
- Frequency 50/60 Hz
- Nominal Current 8.7 A
- Power minimum 2000 VA (2 kVA)

# Site Preparation Checklist



## TAM IV

	<p>Sufficient floor space for instrument and bench space for computer (controller) and accessory (if needed):</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Floor space: 1–1.2 m (width) x 1 m (depth)</li> <li><input type="checkbox"/> Space on each side of instrument: 300 mm</li> <li><input type="checkbox"/> Bench space: 40 cm</li> </ul> <p>Laboratory conditions meet the following requirements:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Dust-free, vibration-free, away from exposure to direct sunlight/air drafts, in a pollution degree 2 environment</li> <li><input type="checkbox"/> Temperature is 15–30°C with a stability of <math>\pm 1^\circ\text{C}</math> over 24 hours</li> <li><input type="checkbox"/> Relative humidity of 5–80% non-condensing</li> <li><input type="checkbox"/> Maximum altitude is 2000 m (6560 ft)</li> </ul>
	<ul style="list-style-type: none"> <li><input type="checkbox"/> Power is 100–240V, 50/60 Hz, 1400 Watts</li> </ul>
	<ul style="list-style-type: none"> <li><input type="checkbox"/> Nitrogen for purging the oil reservoir</li> <li><input type="checkbox"/> Dry gas (air) for purging the cabinet (3 L/min)</li> <li><input type="checkbox"/> Low pressure dual stage regulator that provides 10–50 psig</li> </ul>
	<ul style="list-style-type: none"> <li><input type="checkbox"/> The Customer assumes responsibility for any damage that occurs when the instrument is moved by someone other than a trained TA Instruments Service Representative.</li> </ul>
	<ul style="list-style-type: none"> <li><input type="checkbox"/> Computer meets all hardware requirements</li> <li><input type="checkbox"/> Computer meets all software requirements</li> <li><input type="checkbox"/> Customer's IT personnel has provided Administrative privileges on the controller computer</li> <li><input type="checkbox"/> The Customer's IT personnel will be on site the day of installation</li> </ul>

I hereby acknowledge that all utility requirements have been met per the checklist above and that they will be ready at the agreed time of installation.

**If all utility requirements are not met at the agreed time of installation, additional charges may be incurred for a return Service trip.**

Customer \_\_\_\_\_ DD / MM / YYYY

Company \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Country \_\_\_\_\_

Please send a signed copy of the completed checklist to your local Service representative.

# TA Instruments Offices

For information on our latest products, contact information, and more, see our website at:  
<http://www.tainstruments.com>.

To find your local TA Instruments office and contact information, visit  
<http://www.tainstruments.com/contact/ta-directory/>

TA Instruments – Waters LLC  
Corporate Headquarters  
159 Lukens Drive  
New Castle, DE 19720  
USA

Telephone: 302-427-4000

Fax: 302-427-4001

Email: [info@tainstruments.com](mailto:info@tainstruments.com)