


RHEOLOGY – FROM LIQUIDS TO SOLIDS



©2018 TA Instruments


TA | TAINSTRUMENTS.COM

Topics in the Waters – TA Characterization Workshops


- Optimizing DSC and TGA Testing
- GPC Techniques and Applications
- Rheology from Liquids to Solids
- DMA: Rheology of Solid-Like Materials




Rheology: An Introduction



Rheology: The study of the flow and deformation of matter





Basic Material Behaviors



Rheological Testing - Rotational

➤ 2 Basic Rheological Methods

1. Apply Force (Torque) and measure Deformation and/or Deformation Rate (Angular Displacement, Angular Velocity) - Controlled Force, Controlled Stress
2. Control Deformation and/or Deformation Rate and measure Force needed (Controlled Displacement or Rotation, Controlled Strain or Shear Rate)



TA Instruments Rotational Rheometers

Discovery HR Rheometer



Combined Motor and Transducer
"Native Mode" = Force (Stress)

ARES-G2




Separate Motor and Transducer
"Native Mode" = Deformation/Deformation Rate
(Strain/Shear Rate)


With computer feedback, the instruments can do both,
deformation/deformation rate control and force control.



Rotational Rheometer Designs



Steady Simple Shear Flow



$$\dot{\gamma} = dv_x/dy = V_0/H \quad \text{Shear Rate, sec}^{-1}$$


$$\sigma = F/A \quad \text{Shear Stress, Pascals}$$

$$\eta = \sigma/\dot{\gamma} \quad \text{Viscosity, Pa-sec}$$

➤ These are the fundamental flow parameters. Shear rate is always a change in velocity with respect to distance.




Parallel Plate Shear Flow



Stress	σ	Force or torque
Strain	γ	Linear or angular displacement
Strain Rate	$\dot{\gamma}$	Linear or angular velocity




Representative Flow Curve



TA | TAINSTRUMENTS.COM

Simple Shear Deformation

Top Plate Displacement = X_0 ; Area = A ; Force = F



Bottom Plate Displacement = 0

$$x = (y/H) \cdot X_0$$

$$\gamma = dx_x/dy = X_0/H \quad \text{Shear Strain, unitless}$$

$$\sigma = F/A \quad \text{Shear Stress, Pascals}$$


$$G = \sigma/\gamma \quad \text{Modulus, Pa}$$

➤ These are the fundamental deformation parameters. Shear strain is always a change in displacement with respect to distance.

TA | TAINSTRUMENTS.COM

Stress Relaxation

Stress Relaxation of Soy Flour, Overlay



- Instantaneous Strain
- Note the decrease in the modulus as a function of time.



Rheological Parameters

FLUIDS TESTING

Parameter	Shear	Elongation	Units
Rate	$\dot{\gamma}$	$\dot{\epsilon}$	Seconds ⁻¹
Stress	σ	τ	Pascals
Viscosity	$\eta = \sigma/\dot{\gamma}$	$\eta_E = \tau/\dot{\epsilon}$	Pascal-seconds

SOLIDS TESTING


Parameter	Shear	Elongation	Units
Strain	γ	ϵ	Unitless
Stress	σ	τ	Pascals
Modulus	$G(t) = \sigma/\gamma$	$E = \tau/\epsilon$	Pascals

Rheological Parameters


CREEP TESTING

Parameter	Shear	Elongation	Units
Stress	σ	τ	Pascals
Strain	γ	ϵ	Unitless
Compliance	$J = \gamma/\sigma$	$D = \epsilon/\tau$	1/Pascals

Creep Testing



Geometry Options



Examples for Common Configurations

Geometry	Examples
Concentric Cylinder	Coatings Beverages Slurries (vane rotor option) Starch pasting
Cone and Plate	Low viscosity fluids Viscosity standards Sparse materials Polymer melts in steady shear
Parallel Plate	Widest range of materials Adhesives Polymer melts Hydrogels Asphalt Curing of thermosetting materials Foods Cosmetics
Torsion Rectangular	Thermoplastic solids Thermoset solids



Converting Machine to Rheological Parameters in Rotational Rheometry

$$\frac{M \times K_{\sigma}}{\Omega \times K_{\gamma}} = \frac{\sigma}{\dot{\gamma}} = \eta$$

$$\frac{M \times K_{\sigma}}{\theta \times K_{\gamma}} = \frac{\sigma}{\gamma} = G$$

Machine Parameters

M: Torque

Ω : Angular Velocity

θ : Angular Displacement

Conversion Factors

K_{σ} : Stress Conversion Factor

K_{γ} : Strain (Rate) Conversion Factor

Rheological Parameters

σ : Shear Stress (Pa)

$\dot{\gamma}$: Shear Rate (sec^{-1})

η : Viscosity (Pa-sec)

γ : Shear Strain

G : Shear Modulus (Pa)




The conversion factors, K_{σ} and K_{γ} will depend on the following:

Geometry of the system – concentric cylinder, cone and plate, parallel plate, and torsion rectangular

Dimensions – gap, cone angle, diameter, thickness, etc.



Shear Rate and Shear Stress Calculations

Conversion Factor	Geometry		
	Couette	Cone & Plate	Parallel Plates
			
K_{γ}	$R_{avg}/(R_o - R_i)$	$1/\beta$	R/h
K_{σ}	$1/(2 * \pi R_i^2 L)$	$3/(2\pi R^3)$	$2/(\pi R^3)$



Cone & Plate and Parallel Plate Geometric Considerations

Standard DHR Peltier
plate geometry diameters

20mm

40mm

60mm



Shear Stress

Shear Stress

At given torque,

Increase in diameter → Decrease in shear stress

Shear Rate

At given angular velocity

Increase in cone angle → Decrease in shear rate

Increase in gap (parallel plate) → Decrease in shear rate

So, for low viscosity fluids, use the largest diameter cone or plate.
For high viscosity fluids, use the smallest diameter cone or plate



Geometry Size Selection

• DHR

- Most common is 40-mm parallel plate; 1000 micron gap
- Use 60-mm cone and plate and parallel plate for low viscosity materials, say, up to 100 mPa-sec, but 40-mm geometries can often handle these materials too.
- 20-mm plates are often used at higher viscosities.
- 25-mm parallel plates are the preferred choice for polymer melts.
- 40-mm 2-degree is the most common cone geometry. This is often used to verify an instrument with a viscosity standard.
- 8-mm plates are often used for pressure sensitive adhesives and for asphalt around room temperature.


• ARES-G2

- The most common geometry on the ARES-G2 is the 25-mm parallel plate. Examples would be polymer melts and thermosetting materials.
- Low viscosity fluids are run with 50-mm plates or cone-and-plate.
- Again, 8-mm plates are used for adhesives and asphalt at room temperature.




Rheological Methods – Unidirectional Testing

- Flow
 - Stress/Rate Ramp
 - Stress/Rate Sweep
 - Time sweep/Peak Hold/Stress Growth
 - Temperature Ramp
- Creep (constant stress)
- Stress Relaxation (constant strain)




TA | TAINSTRUMENTS.COM

Idealized Flow Curve




TA | TAINSTRUMENTS.COM


Thixotropy: Up & Down Flow Curves



Rate Ramp




Stress Ramp with Yield



Ramp Selection

- Do a ramp as a preliminary scouting test, prior to the more fundamentally sound rate or stress sweep.
- For rate ramps, a common acceleration rate is 1 sec^{-1} per second. For example, 0 to 100 sec^{-1} in 100 seconds. This is a starting point. The operator can select a rate or a range that is more appropriate for the sample in question.
- Ramp up/Ramp down tests are common for determining thixotropy. The area between the up and down curves is often reported as a thixotropy parameter.
- There have been times when the reproducibility is better with the down curve than it is with the up curve.
- Ramps are good for characterizing materials that may slip or exude from the gap as the shear rate is increased. Often one can get to higher shear rates with ramps than with sweeps because one doesn't dwell at the high rates as long.
- Stress ramps are often used to get the yield point of a material. Sometimes these are not always clear-cut. Also, one has to be cautious when working with models. There have been instances where negative (!?) yield stresses are determined by software for the selected model.
- For stress ramps, use the Step Termination feature to prevent over-speed.

Rate Sweeps - Polymer Melt Flow Curves



Structured fluids

Structured fluids are mostly colloidal dispersions and can be classified into three categories.

- Suspension Solid particles in a Newtonian fluid
- Emulsion Fluid droplets in a Newtonian fluid
- Foam Gas in a fluid (or solid)


■ Examples:

- Paints
- Coatings
- Inks
- Personal Care Products
- Cosmetics
- Foods


■ Properties:

- Yield Stress
- Non-Newtonian Viscous Behavior
- Thixotropy
- Elasticity


Steady State Flow at 25 C - Coatings



Viscosity curve of various structured fluids



Viscous response of suspensions




Rheological parameters of interest: Yield stress, viscosity, time dependence, linear viscoelasticity

H.M. Laun *Angew. Makromol. Chem.* **335**, 124 (1984)



Normal Force Measurements with Cone & Plate




Normal Stress Difference:

- In steady flow, polymeric materials can exert a force that tries to separate the cone and the plate.
- A parameter to measure this is the Normal Stress Difference, N_1 , which equals $\sigma_{xx} - \sigma_{yy}$ from the Stress Tensor.
- $N_1 = 2F/(\pi R^2)$, where F is the measured force.
- $\Psi_1 = N_1/\dot{\gamma}^2$ This is the primary normal stress coefficient.



Effect of HDPE Variations in Blow Molding



No differences in MFI, Viscosity, or GPC!

M-2 produces heavier bottles in blow molding due to increased parison swell




Rate and Stress Sweeps

- Sweeps are preferable to ramps because the material is given a chance to equilibrate at a particular shear rate or stress.
- These are useful tests to see how materials will perform in flow, such as transporting fluid through pipes and tubing, after structure has been broken.
- The most common rate range is 0.1 to 100 $1/\text{sec}$, 5 points per logarithmic decade.
- The Steady State sensing feature is a useful tool to perform valid rate or stress sweeps in the shortest amount of time.
- For materials that exhibit flow instability, the ramp may be preferable. With sweeps, the material is exposed to high shear rates for prolonged periods, whereas, with ramps, the dwell time at a particular rate is shorter.




Creep Testing on US and UK Paints




TA | TAINSTRUMENTS.COM

Stress Relaxation




TA | TAINSTRUMENTS.COM


Flow Temperature Ramp – Printing Inks



Tack Testing



Dynamic Testing




Dynamic Rheological Parameters

Parameter	Shear	Elongation	Units
Strain	$\gamma = \gamma_0 \sin(\omega t)$	$\epsilon = \epsilon_0 \sin(\omega t)$	---
Stress	$\sigma = \sigma_0 \sin(\omega t + \delta)$	$\tau = \tau_0 \sin(\omega t + \delta)$	Pa
Storage Modulus (Elasticity)	$G' = (\sigma_0/\gamma_0) \cos \delta$	$E' = (\tau_0/\epsilon_0) \cos \delta$	Pa
Loss Modulus (Viscous Nature)	$G'' = (\sigma_0/\gamma_0) \sin \delta$	$E'' = (\tau_0/\epsilon_0) \sin \delta$	Pa
Tan δ	G''/G'	E''/E'	---
Complex Modulus	$G^* = (G'^2 + G''^2)^{0.5}$	$E^* = (E'^2 + E''^2)^{0.5}$	Pa
Complex Viscosity	$\eta^* = G^*/\omega$	$\eta_E^* = E^*/\omega$	Pa-sec

Cox-Merz Rule for Linear Polymers: $\eta^*(\omega) = \eta(\dot{\gamma})$ @ $\dot{\gamma} = \omega$


Rheological Methods – Dynamic Oscillatory Testing

- Dynamic Strain Sweep/Dynamic Stress Sweep
- Isothermal Dynamic Time Sweep
- Isothermal Dynamic Frequency Sweep
- Dynamic Temperature Ramp
- Dynamic Temperature Sweep at 1 or Multiple Frequencies.



TA | TAINSTRUMENTS.COM


Linear Viscoelasticity




TA | TAINSTRUMENTS.COM

Stability is Related to Structure in Inks

Ink Samples: Oscillation Stress Sweeps @ 6.28 rad/s




Frequency Sweep on PDMS




Surface Defects during Pipe Extrusion

Indicated by
Elasticity at
low
frequency


Caused by
Broader
MWD




Coatings Frequency Sweep



Example of Cox-Merz Rule




Dynamic testing: Dependence on M_w and MWD




When shifting along an axis of -1, all the curves can be superposed, unless the width of the MWD is not the same.

Zero shear viscosity increases with increasing MW



MWD and Dynamic Moduli



• The storage and loss modulus of a typical polymer melt cross over between 1 and 100 rad/s.


- $\omega_c \propto 1/M_w$.
- $G_c \propto 1/\text{MWD}$.

Higher crossover frequency = lower M_w
Higher crossover Modulus = narrower MWD




ETC Application: TTS on Polymer Melt


Polystyrene Frequency Sweeps from 160°C to 220°C



Idealized Flow Curve - Polymers



Tack and Peel of Adhesives




Desirable PSA characteristics

Tack: high G' at low frequency
 Peel: low G' at high frequency


Correlation of Rheological Parameters to Adhesive Properties

Property	Rheological Properties	Practical Adhesive Property
Tack	<ul style="list-style-type: none"> Low $\tan \delta$ and Low G' Low Cross-links ($G'' > G'$ @ ~1 Hz) 	High Tack
Shear Resistance	<ul style="list-style-type: none"> High G' @ < 0.1 Hz High Viscosity @ Low Shear Rates 	High Shear Resistance
Peel Strength	<ul style="list-style-type: none"> High G'' @ ~> 100 Hz 	High Peel Strength
Cohesive Strength	<ul style="list-style-type: none"> High G', low $\tan \delta$ 	High Cohesive Strength
Adhesive Strength	<ul style="list-style-type: none"> High G'', high $\tan \delta$ 	High Adhesion Strength with Surface


Dynamic Time Sweep: Material Response




Cure of a "5 minute" Epoxy



Isothermal Cure of Tire Compound Effect of Curing Temperature




Dynamic time sweep following shear




This shows how rapidly the structure is rebuilding.

Figure 7: Recovery of structure


Dynamic Temperature Ramp - Torsion



Effect of Molecular Weight




Effect of Crosslinking



Dynamic Testing

- Of all the tests performed on a rheometer, dynamic oscillatory testing is the most common. Typically, this is the most convenient way of getting a material's viscous and elastic nature.
- Dynamic Stress or Strain sweeps are useful for determining a dynamic yield point of a material and can suggest strains or stresses to use in subsequent tests, such as frequency sweeps or temperature ramps.
- Dynamic Frequency sweeps are the most useful tests for characterizing polymer melts and adhesives. They provide information on the molecular weight and molecular weight distribution of a material.
- Time temperature superposition can often be used to provide knowledge beyond the usual limits of 0.1 to 100 rad/sec.
- The rheometer can be used as a DMA to provide glass transition temperatures and thermal mechanical integrity.



TA Instruments Rotational Rheometers

Discovery HR Rheometer



Combined Motor and Transducer
"Native Mode" = Force (Stress)

ARES-G2



Separate Motor and Transducer
"Native Mode" = Deformation/Deformation Rate
(Strain/Shear Rate)

With computer feedback, the instruments can do both,
deformation/deformation rate control and force control.



DHR Instrument Specifications

Specification	HR-3	HR-2	HR-1
Bearing Type, Thrust	Magnetic	Magnetic	Magnetic
Bearing Type, Radial	Porous Carbon	Porous Carbon	Porous Carbon
Motor Design	Drag Cup	Drag Cup	Drag Cup
Minimum Torque (nN.m) Oscillation	0.5	2	10
Minimum Torque (nN.m) Steady Shear	5	10	20
Maximum Torque (mN.m)	200	200	150
Torque Resolution (nN.m)	0.05	0.1	0.1
Minimum Frequency (Hz)	1.0E-07	1.0E-07	1.0E-07
Maximum Frequency (Hz)	100	100	100
Minimum Angular Velocity (rad/s)	0	0	0
Maximum Angular Velocity (rad/s)	300	300	300
Displacement Transducer	Optical encoder	Optical encoder	Optical encoder
Optical Encoder Dual Reader	Standard	N/A	N/A
Displacement Resolution (nrad)	2	10	10
Step Time, Strain (ms)	15	15	15
Step Time, Rate (ms)	5	5	5
Normal/Axial Force Transducer	FRT	FRT	FRT
Maximum Normal Force (N)	50	50	50
Normal Force Sensitivity (N)	0.005	0.005	0.01
Normal Force Resolution (mN)	0.5	0.5	1



TA | TAINSTRUMENTS.COM

DHR Instrument Model Features

All Discovery Hybrid Rheometers Feature:	Moving from HR-1 to HR-2 Adds:	Moving from HR-2 to HR-3 Adds:
<ul style="list-style-type: none"> - Patented Ultra-low Inertia Drag-Cup Motor - Single-Thrust & Dual-Radial Bearing Design - Patented Second Generation Magnetic Bearing - Nano-Torque Motor Control - High-Resolution Optical Encoder - Superior Stress and Strain Control - Force Rebalance Normal Force (FRT) - Patented Smart Swap Geometries - True Position Sensor (Patent Pending) - Ultra-low Compliance Single-Piece Frame - Heat and vibration Isolated Electronics Design - Smart Swap™ Temperature Systems - Superior Peltier Technology - Patented Heat Spreader Technology - Patented Active Temperature Control - Color Display - Capacitive Touch Keypad - TRIOS Software - Navigator Software - Electronic Bearing Lock - NIST Traceable Torque Calibration 	<ul style="list-style-type: none"> - 5X better low torque in Oscillation - 2X better low torque in steady Shear - 25% Higher torque - 2X better NF Sensitivity - Direct Strain Oscillation - Fast data sampling - Transient Data Acquisition/LAOS - Stress Growth (Transient NF) - Access to UV Curing Options - Access to SALS Option - Access to Interfacial Options 	<ul style="list-style-type: none"> - 4X better low torque in Oscillation - 2X better low torque in steady Shear - Optical Encoder Dual Reader (pat. Pend.) - 5X better angular resolution - 3X better phase angle resolution - No encoder drift



TA | TAINSTRUMENTS.COM

ARES-G2 Specifications

Force/Torque Rebalance Transducer (Sample Stress)

Transducer Type	Force/Torque Rebalance
Transducer Torque Motor	Brushless DC
Transducer Normal/Axial Motor	Brushless DC
Minimum Transducer Torque in Oscillation	0.05 $\mu\text{N}\cdot\text{m}$
Minimum Transducer Torque in Steady Shear	0.1 $\mu\text{N}\cdot\text{m}$
Maximum Transducer Torque	200 $\text{mN}\cdot\text{m}$
Transducer Torque Resolution	1 $\text{nN}\cdot\text{m}$
Transducer Normal/Axial Force Range	0.001 to 20 N
Transducer Bearing	Groove Compensated Air

Drive Motor (Sample Deformation)

Maximum Motor Torque	800 $\text{mN}\cdot\text{m}$
Motor Design	Brushless DC
Motor Bearing	Jeweled Air, Sapphire
Displacement Control/Sensing	Optical Encoder
Strain Resolution	0.04 μrad
Min. Angular Displacement in Oscillation	1 μrad
Max. Angular Displacement in Steady Shear	Unlimited
Angular Velocity Range	1 x 10 ⁻⁴ rad/s to 300 rad/s
Angular Frequency Range	1 x 10 ⁻² rad/s to 628 rad/s
Step Change in Velocity	5 ms
Step Change in Strain	10 ms

Orthogonal Superposition and DMA modes

Motor Control	Force Rebalance Transducer
Minimum Transducer Force in Oscillation	0.001 N
Maximum Transducer Force	20 N
Minimum Displacement in Oscillation	0.5 μm
Maximum Displacement in Oscillation	50 μm
Displacement Resolution	10 nm
Axial Frequency range	1 x 10 ⁻⁴ Hz to 16 Hz

Stepper Motor


Movement/Positioning	Micro-stepping Motor/ Precision Lead Screw
Position Measurement	Linear Optical Encoder
Positioning Accuracy	0.1 micron

Temperature Systems

Smart Swap	Standard
Forced Convection Oven, FCO	-150 °C to 600 °C
FCO Camera Viewer	Optional
Advanced Peltier System, APS	-10 °C to 150 °C
Peltier Plate	-40 °C to 180 °C
Sealed Bath	-10 °C to 150 °C

JMENTS.COM

Minimum Torque Specs




DHR Accessories – Visual Display



Peltier Plate
Temperature Systems



Advanced Peltier Plate



Dual Stage Peltier Plate



Upper Heated Plate for
Peltier Plate



Peltier Concentric
Cylinders



Electrically Heated
Cylinder (EHC)



Pressure Cell




Electrically Heated
Plates




Environmental Test
Chamber



Relative Humidity
Accessory



Modular Microscope
(MMA)



Optical Plate

You can see the
updated list of
accessories on our
website,
www.tainstrument.com.



DHR Accessories – Visual Display



Small Angle Light
Scattering



Interfacial Accessories



Tribo-Rheometry
Accessory



Magneto-Rheology



Electro-Rheology



UV Curing Accessories



Dielectric Measurement



Immobilization Cell




Starch Pasting Cell




Dynamic Mechanical
Analysis




ARES-G2 Accessories




Forced Convection Oven
(FCO)



Advanced Peltier
System (APS)




Orthogonal
Superposition & 2D-
SAOS



Dynamic Mechanical
Analysis (DMA)



UV Curing Accessory



Dielectric Thermal
Analysis Accessory
(DETA)



Extensional Viscosity
Fixture (EVF)




Tribo-Rheometry
Accessory



Air Chiller System



Cone and Partitioned
Plate Accessory



Interfacial Rheology



Electrorheology (ER)
Accessory



DMA Capabilities




Motor Control	Force Rebalance Transducer
Minimum Force in Oscillation	0.1 N
Maximum Axial Force	50 N
Minimum Displacement in Oscillation	1 μm
Maximum Displacement in Oscillation	100 μm
Displacement Resolution	10 nm
Axial Frequency Range	1×10^{-5} to 16 Hz

- DHR Film/Fiber Tension Clamp Accessory kit
- DHR 3-Point Bending Clamp Accessory kit
- DHR Cantilever Bending Clamp Accessory kit

The DMA capabilities of the DHR and ARES-G2 are unique for commercial rheometers.








ARES-G2 DMA Mode

Features and Benefits

- Exclusive to the ARES-G2 rheometer
- Wide range of geometries:
 - 3-Point Bending
 - Film/Fiber Tension
 - Single and Dual Cantilever (Clamped Bending)
 - Parallel Plates Compression
- Axial Force Control tracks material stiffness and automatically adapts static load
- AutoStrain adjusts applied strain to changing sample stiffness
- Responsive FCO temperature control: $-150^{\circ}C$ to $600^{\circ}C$
- Sample visualization with FCO camera





ARES-G2 DMA Testing



DMA Specifications


	RSA G2	DMA 850	ARES G2 DMA	DHR DMA (optional)
Max Force	35N	18N	20N	50N
Min Force	0.0005N	0.0001N	0.001N	0.1N
Frequency Range	1e-5 to 628 rad/s (1.6e-6 to 100 Hz)	6.28e-3 to 1250 rad/s (0.001 to 200 Hz)	6.3e-5 to 100 rad/s (1.0e-5 to 16 Hz)	6.3e-5 to 100 rad/s (1.0e-5 to 16 Hz)
Dynamic Deformation Range	+/- 0.05 to 1,500µm	+/- 0.005 to 1e4 µm	+/- 1 to 50 µm	+/- 1 to 100 µm
Control Stress/Strain	Control Strain (SMT)	Control Stress (CMT)	Control Strain (CMT)	Control Stress (CMT)
Heating Rate	0.1°C to 60°C/min	0.1°C to 20°C/min	0.1°C to 60°C/min	0.1°C to 60°C/min
Cooling Rate	0.1°C to 60°C/min	0.1°C to 20°C/min	0.1°C to 60°C/min	0.1°C to 60°C/min

SER2 for DHR Rheometers




This is an interesting application of using the rotational rheometer to determine elongational viscosity


SER2 DHR Data



Extensional Viscosity Measurements




Fix drum connected to transducer




Rotating drum connected to the motor:

- rotates around its axis
- rotates around axis of fixed drum

LDPE (High branching)



Large Angle Oscillatory Shear: Higher Harmonics



The transition from linear to nonlinear regime for viscoelastic materials shows in a decrease of G' & G'' and an increase of the magnitude of the harmonics.

Trios Software



Trios Software

- Trios is non-restricted. Anyone can download Trios from our website at no cost and install it on any computer.
- Trios is versatile and friendly. It combines the best of Rheology Advantage (AR Rheometers) and Orchestrator (ARES) and has additional features too.
- Trios has all the necessary procedures for conventional as well as specialized testing.
- In Trios, with the DHR rheometer, you can add steps and modify steps that haven't started while the test is in progress.
- There is the Auto-tension feature, which is used in torsion rectangular testing and linear DMA testing.
- There is the Auto-strain feature, which allows the strain to change during an experiment based on the torque. This is useful for temperature ramps and curing tests.
- The graphics quality is excellent. Formatting graphs is easy with the Automatic options and very versatile with the Customized graphing.
- Data can be read in spreadsheet form. If you copy and paste into Excel, you will get the data and the headings.
- There are numerous variables to select from, but you can also construct user-defined variables.
- Trios can read data from Rheology Advantage and Orchestrator data files.



Trios Screen

The screenshot displays the Trios software interface. At the top, the title bar reads "Offline-Discovery HR-3 : TA Instruments Trios v4.3.0.38388". Below this is a menu bar with options: Experiment, Scripting, Instrument, Engineering, and Script Buttons. A toolbar contains icons for Offline, Geometry, Setup Procedure, Start Experiment, Stop Experiment, Point Display, Image View, Log, and Help.

The main window is divided into several panes. On the left is the "File Manager" pane showing a tree view with "Experiments" and "[Experiment 2]". Below this is a "Log" pane with a table for recording experiment data. The table has columns for "Information", "Code", "Log Time", and "Description".

The central pane shows the setup for "Experiment 2". It includes fields for "Sample: S600", "Geometry: 40mm 2.0° cone plate, Pelier plate Steel", and "Procedure: Field Service Test - stress ramp on S600 procedure". Below these fields is a checkbox for "1: Flow Ramp 25°C, 00:03:00h:mm:ss, 0 to final 5000uN/m".

A red box highlights the text: "The Control Panel is shown on the next slide."



You can control many things from the Control Panel section:

1. Gap
2. Temperature
3. Spindle speed
4. The laser for SALS
5. Axial force.

TA | TAINSTRUMENTS.COM

Geometry Information

Geometry: 40mm parallel plate, Peltier plate Steel

Diameter: 40 mm
 Gap: 1.0 mm
 Loading gap: 45.0 mm
 Trim gap offset: 0.05 mm

Material: Steel
 Environmental system: Peltier plate
 Minimum sample volume is 1.25664 cm³

Constants

Gap temperature compensation
 Expansion coefficient: 0.6 μm/°C
☒ Enabled

Compliance: 1.58 mrad/N.m
 Geometry inertia: 7.117 μN.m.s²
 Friction: 0.897 μN.m/(rad/s) ☒ Enabled
 Stress constant: 79577.5 Pa/N.m
 Strain constant: 20.0 1/rad
 Normal stress constant: 1591.55 Pa/N
 Fluid density constant: 8.378e-11

This information will be stored with the geometry file. It will be accessed when the magnetic reader scans the magnetic strip at the top of the geometry.

There is a friendly Geometry wizard that guides you through preparing a new geometry.

Double-pressing the Bearing Lock icon on the touch-pad and loading the geometry in the same location. enables you to retain at least most of the mapping

TA | TAINSTRUMENTS.COM

Trios Software Procedures

Oscillation
Flow
Step (Transient)
Other
Conditioning

Frequency
Temperature Ramp
Temperature Sweep
Time
Amplitude
Fast Sampling
Manual

Oscillation
Flow
Step (Transient)
Other
Conditioning

Axial
Extensional
Event
Temperature Ramp IsoForce
Temperature Ramp IsoStrain
Arbitrary Wave

Oscillation
Flow
Step (Transient)
Other
Conditioning

Temperature Ramp
Temperature Ramp
Ramp
Sweep
Peak Hold

Oscillation
Flow
Step (Transient)
Other
Conditioning

Sample
Options
End Of Test
Sample Loading

Step (Transient)
Oscillation
Flow
Step (Transient)
Other
Conditioning

Creep
Stress Relaxation
Stress Growth
Repeated Creep

These are the choices in Trios for the DHR.
There are minor differences with the
choices for the ARES-G2.



Oscillation Procedure

Procedure: Sample Frequency Sweep


1: Oscillation Frequency

Environmental Control
Temperature 25 °C ☐ Inherit Set Point
Soak Time 00:03:00 hh:mm:ss ☒ Wait For Temperature

Test Parameters
Strain % 1.0 %
Logarithmic sweep
Angular frequency 100.0 to 0.1 rad/s
Points per decade 5



Display of Frequency Sweep Data



Flow Procedure

1: Flow Sweep

Environmental Control

Temperature °C ☐ Inherit Set Point

Soak Time hh:mm:ss ☒ Wait For Temperature

Test Parameters

Logarithmic sweep

Shear rate to 1/s

Points per decade

☒ Steady state sensing

Max. equilibration time hh:mm:ss


Sample period hh:mm:ss

% tolerance

Consecutive within

☐ Scaled time average

Display of Flow Data



TA Instruments Rubber Testing Products




- Rheometer
- Viscometer
- Density Meter
- Hardness Tester



RPA Elite



RPA elite Rubber & Polymer Analyzer



We obtained the RPA and others for rubber testing when we acquired Scarabaeus recently.


Densitometer

Shore A
Hardness



TA | TAINSTRUMENTS.COM

ASTM D5289: MDR/RPA Isothermal Curing




TA | TAINSTRUMENTS.COM

Need Some Assistance?

- Instrument Manuals
- Help Feature in Trios
- TA Instruments website
- TA Instruments Rheology Helpline
 - rsupport@tainstruments.com




Website: www.tainstruments.com



Website: www.tainstruments.com

Products ▾ About TA ▾ Videos ▾ **Support ▾** Training ▾ News & Events ▾ Sales Promos

Service Support	Application Support	Software Downloads & Support	Support Plans
Service Support Helpline	Applications Support Helpline	Software Downloads	Lifetime Support Plan
Site Preparation Guides	Tech Tips	Instruments sorted by software	Premium Support Plan
The IQ/OQ Product Offering	Applications Notes Library	Software Sorted by Instruments	Plus Support Plan
Calibration with Certified Standards	Training	Report a Bug	Basic Support Plan
Safety Data Sheets		Request a Feature	Performance Maintenance Visit (PMV)
Supported Instruments			Academic Support Plan
Service Shop			ElectroForce Support Plans

 TAINSTRUMENTS.COM

Thank you for attending

Don't forget <http://www.tainstruments.com/>

If you have questions or need help
Contact the TA Instruments Technical Helpline
<http://www.tainstruments.com/support/applications/applications-hotline/>

