

ElectroForce® 5500 Series



Site Preparation Guide

Table of Contents

Table of Contents	2
Ideal Setup.....	3
System Components.....	4
Instrument Measurements.....	5
Utility Requirements	6–7
Power	6
Miscellaneous.....	7
Site Preparation Checklist	8
TA Instrument Offices.....	9



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab



Customer

Ideal Setup



IDEAL PLACEMENT AND BENCH MEASUREMENTS

Select a location with adequate floor and ceiling space and a rigid laboratory bench that is level and is in a vibration-free environment. Bench must be rated to support several hundred pounds.




Table width: 1.5 m (5 ft)


Table depth: 1.2 m (4 ft)

Distance from the wall:
0.15 m (0.5 ft) min.

System Components



MAIN SYSTEM COMPONENTS



- A. Test Instrument
- B. Computer Tower
- C. Computer Monitor
- D. PCI Electronics Box
- E. Power Supply (Axial)
- F. Emergency Stop

Instrument Measurements



5500 SERIES TEST INSTRUMENT



Height: 490.4 mm (19.3 in)

Width: 269.2 mm (10.6 in)

Depth: 203 mm (18 in)

Weight: 16.2 kg (35.6 lbs)



POWER SUPPLY

Height: 230 mm (9 in)

Width: 405 mm (16 in)



Depth: 510 mm (20 in)



Utility Requirements



POWER

Item	Requirement
Instrument Power	<ul style="list-style-type: none">• 104–132V, 50/60 Hz, 15A• 207–250V, 50/60 Hz, 10A• Neutral to Ground (NG) voltage max 0.5 volt• Safety ground per local regulation
PCI Box Power	<ul style="list-style-type: none">• 104–120V, 50/60 Hz Hz, 10 A• 207–230V, 50/60 Hz Hz, 6.3 A
UPS (Optional)	<ul style="list-style-type: none">• 104–120V, 50/60 Hz, 1A• 207–230V, 50/60 Hz, 0.5A
Power cords provided	<ul style="list-style-type: none">• 5-15 plug for 120V systems• 6-20P plug for 230V systems• International: Line power cord provided is based on country <div style="display: flex; justify-content: center; align-items: center; gap: 20px;"><div style="text-align: center;"> 6-20P</div><div style="text-align: center;"> 5-15</div></div>



Use power cords with plugs appropriate for your circuit.



Supply voltages lower than indicated may result in a degradation of performance.



Ensure that the mains assigned do not also supply power to noise generating equipment nearby, such as motors, welders, transformers, etc.



An independent heavy GROUND wire must be provided through the power hookup. Improper grounding may cause severe damage for which the supplier will not accept responsibility. All power strips must be fully grounded and carry the ground through to the sockets into which the computer is plugged.

Utility Requirements



MISCELLANEOUS

Item	Requirement
Operating temperature	18°C (64°F) to 30°C (86°F)
Relative humidity	40–65% (non-condensing)

TA Instruments Offices

For information on our latest products, contact information, and more, see our website at:
<http://www.tainstruments.com>.

To find your local TA Instruments office and contact information, visit
<http://www.tainstruments.com/contact/ta-directory/>

TA Instruments – Waters LLC
Corporate Headquarters
159 Lukens Drive
New Castle, DE 19720
USA

Telephone: 302-427-4000
Fax: 302-427-4001
Email: info@tainstruments.com