




ElectroForce® test instruments featuring patented linear motion technologies and WinTest® controls provide a revolutionary approach to mechanical fatigue and dynamic characterization. The ElectroForce family of test instruments includes a full range of force and performance capabilities for a wide variety of test applications based on the most unique motor design in the industry. The end result is billions of cycles of unmatched reliable performance in a dynamic test instrument with precision, accuracy, and ease of use for a wide range of applications.

ElectroForce® Load Frames	2
ElectroForce TestBenches	4
ElectroForce Instrument Features	6
ElectroForce Linear Motor	8
Software	10
Load Frame Instruments - 3100	12
Load Frame Instruments – 5500	13
Load Frame Instruments - 3200	14
Load Frame Instruments - 3300	16
Load Frame Instruments – 3500	18
TestBench Instruments	20
Planar Biaxial Instruments	22
Applications	24
Accessories	26
Support	28
Specifications	30

ELECTROFORCE® | LOAD FRAMES



ELECTROFORCE® TESTBENCHES






ELECTROFORCE® INSTRUMENT FEATURES

For more than 20 years, ElectroForce test instruments have been successfully utilized in a wide variety of mechanical characterization and fatigue applications.

- Patented linear motor operates without friction, an important feature for ultra-durable and high-precision testing
- Efficient, direct electromagnetic conversion to force results in greater acceleration, high frequencies and high velocities
- Intuitive software design to simplify test setup and a flexible hardware platform for changing test needs
- Powered from a standard electrical outlet, requiring no additional infrastructure, air conditioning or water cooling
- Air-cooled, clean-room compatible and whisper-quiet operation in a compact, space-saving package
- Energy efficient and environmentally friendly by using pollution-free, non-toxic technologies and oil-free design



ELECTROFORCE® LINEAR MOTOR | TECHNOLOGY

Reliability that won't let you down

The flexural suspension is engineered to guide the magnet assembly without contact or lubrication. The magnet, coil and stators are designed to control temperatures to eliminate performance degradation over decades of use and deliver maintenance-free operation that you can count on for your longest running tests.

Unmatched waveform control and fidelity

Without the friction of rolling or sliding bearings, the ElectroForce® motor design provides the control required for the most sensitive of tests. The motor converts even the smallest of increments of current precisely to specimen force, displacement or pressure. This means applied forces can be controlled to gram force ranges and displacements can be controlled to a micron.


Testing in a variety of environmental conditions

With its efficient, quiet and bearing-free design, the ElectroForce motor is engineered for use in a variety of diverse and challenging environments, including clean rooms, humidity chambers and even radioactive hot cells.

Sized to meet your requirements

The architecture of the ElectroForce motor can be scaled to accommodate a wide range of forces and displacements:

- Maximum forces from 22 N up to 15 kN
- Maximum displacements from 5 mm up to 50 mm.
- Extended Stroke Options add 150mm displacement.




Stationary Coils

Design simplicity provides unmatched performance & billions of maintenance-free cycles

Fast or Slow

The Most Dynamic & Controllable in the **Industry**



Bandwidth for any Test:

ElectroForce[®] motors excel at the full spectrum of testing speeds, due to the low mass of the rare earth moving magnet and stationary coil design. Testing speeds can range from static tests to one cycle per day, and up to frequencies of 300 Hz.

SOFTWARE TECHNOLOGY



The Most Flexible **Control System Available**

Single comprehensive software with intuitive user interface for instrument waveform control and data analysis.

- Comprehensive tool enables quick and easy generation of periodic waveforms (fatigue), ramps (monotonic tests), and more complex block-multiple waveform segments.
- Integrated data acquisition algorithms provide for a variety of data collection techniques including timed acquisition, peak/valley capture, level-crossing, and others
- High speed real-time measurements enable dynamic and precise motor control and up to 10kHz data acquisition rates of sensor and calculated channels
- Advanced controls including multi-channel synchronization of phase and amplitude, and cross-channel compensation
- Advanced limit monitoring and conditional logic features enable users to monitor test performance and automatically change test activity in response to real-time test measurements
- Additional options include:
- Dynamic Mechanical Analysis
- External Waveform Input
- Dynamic Link Library Interface for External Applications

- Provides excellent re-creation of desired waveforms, allowing for improved test control and ultimately better test results
- Allows users to quickly and confidently tune controllers for dynamic linear motors and torsion motors

- Tan Delta, δ
- Glass Transition

TRIOS for DMA:

- Time Temperature Superposition (TTS)
- Peak analysis
- Peak integration
- User-defined Variables

WinTest[®] Tune IQ[®]

The Most Accurate Closed Loop Control Algorithms

Tune IQ software uses advanced proprietary algorithms to simplify the tuning process

 Advanced methods that analyze the dynamic response of the system over a widerange of frequencies, characterizing instrument, sensor and sample for optimal control. This patented multi-frequency approach is superior to a single-point measurement that fails to factor in specimen dynamics

WinTest DMA (Dynamic Mechanical Analysis) A flexible application for advanced viscoelastic property measurements, including:

- Storage Modulus, E'
- Loss Modulus, E"

- Secondary Transitions
- Viscoelastic changes during fatigue
- Hysteresis Loop Raw Data


Most Versatile Analysis and Plotting Package

- Onset point analysis

(HADS) High Accuracy Displacement Sensor The Most Accurate and Precise **Displacement Sensor on the Market**

- 1 nanometer resolution and micron-level of accuracy
- Far exceeds Class A. ASTM E2309 calibrated accuracy
- Extremely low noise to eliminate the need to filter data
- High responsiveness extends the dynamic performance of system
- Single displacement channel to provide both absolute and highresolution measurements
- Linear HADS measurement is standard on 3200, 3300 and 3510 Instruments
- Rotation HADS measurement is standard on 3200-AT and 3300-AT instruments

Sub-micron level noise and control precision, without data filters



The 3100 is a highly sensitive load frame that delivers the reliability and performance of the ElectroForce linear motor in a compact and portable package.

Features and Benefits:

- ElectroForce linear motor capable of precisely controlling force, displacement or strain over a wide range of frequencies in tension or compression
- Compact package with integrated power supply making it the smallest and most portable ElectroForce instrument
- Open test space, micro-adjust assembly and convenient fixturing interfaces provide flexibility for many specimen types and geometries
- Interchangeable force sensors and precise power amplifier enable force control from \pm 0.002 N to ± 22 N

Specification	Model 3100	
Force Range	0.002 – 22 N	
Displacement Range	0.010 – 5 mm	
Frequency Range	0.0001 – 100 Hz	
Instrument Dimensions (HxWxD)	518 x 292 x 172 mm (20.4 x 11.5 x 6.8 in)	
Instrument Weight	15 kg (32 lbs)	



LOAD FRAME INSTRUMENTS 5500



The 5500 is a unique single-beam support load frame system offering open architecture for adaptability to the broadest range of applications. It offers a broad force range, durability, and versatility in an affordable design making it an ideal choice for academic laboratories.

Features and Benefits:

- ElectroForce linear motor capable of precisely controlling force, displacement or strain over a moderate range of frequencies in tension or compression
- Compact package and clean operation for placement into small enclosures such as incubators and glove boxes.
- 3-sided test space access, micro-adjust assembly and convenient fixturing interfaces provide flexibility for many specimen types and geometries
- Optional 24-well Compression Fixture enables many specimens to be loaded simultaneously in a common 24 well cell culture plate
- Interchangeable force sensors and a powerful amplifier enable force control from $\pm\,0.002$ N to ±200 N

Specification	Model 5500	
Force Range	0.002 – 200 N	
Displacement Range	0.025 – 13 mm	
Frequency Range	0.0001 – 20 Hz	
Instrument Dimensions (HxWxD)	490 x 203 x 269 mm (19.3 x 8.0 x 10.6 in)	
Instrument Weight	16 kg (36 lbs)	

The 3200 Series III load frames are the highest precision and most versatile instruments in their force range. With the highest technology sensors and most advanced frame features, they are the premier choice for demanding test applications requiring wide ranges of force, displacement and frequency.

Features and Benefits:

- 225 N or 450 N ElectroForce linear motor options capable of precisely controlling force, displacement or strain over a wide range of frequencies in tension or compression
- A wide range of accessories make it one of the most versatile ElectroForce instruments including numerous environmental accessories and motor options
- A tabletop frame that's clean, durable and quiet allow it to be used in virtually any space from office to lab to cleanroom to manufacturing floor.
- The optional torsion motor integrates a high-resolution optical encoder for control and measurements up to 62 turns
- The optional extended stroke motor adds 150 mm of linear motion for high-elongation tests
- Interchangeable force sensors and precise power amplifier enable force control from \pm 0.002 N to ±450 N

Configurations

- 225 N or 450 N Axial Instruments
- Axial Torsion with 5.6Nm torsion motor
- Axial with 150mm Extended Stroke
 motor
- Multi-Specimen Fatigue
- DMA 3200 with integrated oven for the best DMA results

Accessories

- Ovens, fluid baths and bioreactors for simulating a variety of environments
- Lower force sensors to improve data quality for low-force tests
- Numerous grips and fixtures for a variety of specimen geometries
- DMA and other software options



Axial Configuration

Speci

Force

Displacer

Frequen

Instrument Dim

Instrume Axial Col Fluid Bath Accessory



ification	Model 3220	Model 3230	
e Range	0.002 - 225 N	0.002 - 450 N	
ement Range	0.002 – 13 mm	0.002 – 13 mm	
ency Range	0.0001 - 300 Hz	0.0001 - 300 Hz	
mensions (HxWxD)	1051 x 579 x 522 mm (42 x 23 x 21 in)	1051 x 579 x 522 mm (42 x 23 x 21 in)	
nent Weight Configuration	98 kg (215 lbs)	105 kg (230 lbs)	



The 3300 Series III load frames are highly flexible offering multiple loading ranges and frame configurations. The high performance motor and sensors combined with flexible and convenient frame features make the 3300 the most widely used instrument for testing a broad of materials, components and devices.

Features and Benefits:

- 1,000 N or 3,000 N ElectroForce linear motor options precisely control force, displacement or strain over a wide range of frequencies in tension or compression
- Multiple frame, motor and environmental accessory combinations allow the instrument to be tailored to a wide variety of applications and sample types.
- Test setup is quick and easy with convenience features of motor-adjustable test space, crosshead locks, Micro-adjust and a T-slot baseplate.
- The optional torsion motor integrates a high-resolution optical encoder for control and measurements up to 62 turns or the Extended Stroke motor option adds 150 mm of linear motion
- Interchangeable force sensors and amplifiers enable force control from ±0.02 N to ±3,000 N

Configurations

• 1,000 N or 3,000 N Axial Instruments

- Axial Torsion with 14Nm or 25Nm torsion motor
- Axial with 150mm Extended Stroke motor
- Multi-Specimen Fatigue

Accessories

• Ovens and fluid baths for simulating a variety of environments

with Oven Accessory

- Lower force sensors to improve data quality for low-force tests
- Numerous grips and fixtures for a variety of specimen geometries
- DMA and other software options



Speci

Force


Displacer

Frequer

Instrument Dim Tabletop Axio

> Instrume Tabletop Axio

ification	Model 3310	Model 3330
e Range	0.02 - 1000 N	0.2 - 3000 N
ement Range	0.005 – 25 mm	0.005 – 25 mm
ency Range	0.0001 – 100 Hz	0.0001 – 100 Hz
mensions (HxWxD) ial Configuration	1283 x 686 x 555 mm (50.5 x 27 x 22 in)	1283 x 686 x 555 mm (50.5 x 27 x 22 in)
nent Weight ial Configuration	185 kg (408 lbs)	185 kg (408 lbs)



Extended Stroke Configuration



The 3500 series includes the highest force electrodynamic load frames available. With multiple force and frame configurations they deliver up to 15,000N of force for strength and dynamic characterization test of a wide variety of materials, components and devices.

Features and Benefits:

- 7,500 N or 15,000 N ElectroForce linear motor options for reliably controlling force, displacement or strain over a wide range of frequencies in tension or compression
- Multiple frame, motor and environmental accessory combinations allow the instrument to be tailored to a wide variety of applications and sample types.
- HADS displacement sensor is standard on the 3510 Series II for even higher resolution and accuracy of displacement measurements.
- Torsion motors offer multi-turn capability with measurements up to 40 turns.
- Interchangeable force sensors and amplifiers enable force control from \pm 1 N to $\pm 15,000$ N

Configurations

- Upper motor 7,500 N (Model 3510)
- Lower motor 7,500 N (Model 3520)
- Lower motor 15,000 N (Model 3550)
- Axial Torsion with 49Nm or 70Nm torsion motor is available for each model

Accessories

- Ovens and fluid baths for simulating a variety of environments
- Lower force sensors to improve data quality for low-force tests
- Numerous grips and fixtures for a variety of specimen geometries
- DMA and other software options

Specification	Model 3510	Model 3520	Model 3550
Force Range	1 – 7,500 N	1 – 7,500 N	1 – 15,000 N
Displacement Range	0.010 – 50 mm	0.025 – 50 mm	0.025 – 50 mm
Frequency Range	0.0001 – 100 Hz	0.0001 – 50 Hz	0.0001 – 50 Hz
Instrument Dimensions (HxWxD) Axial Configuration	2510 x 979 x 813 mm (99 x 39 x 32 in)	2498 x 861 x 756 mm (98 x 34 x 30 in)	2498 x 861 x 756 mm (98 x 34 x 30 in)
Instrument Weight Axial Configuration	1000 kg (2200 lbs)	644 kg (1420 lbs)	816 kg (1800 lbs)



3510 Axial Configuration



3550 Axial Torsion with Oven Accessory

TESTBENCH INSTRUMENTS

The TestBench series is a highly flexible mechanical test platform built upon a breadboard that enables user to relocate and orient motors and force sensors in multiple configurations. This offers a uniquely open platform for demanding and creative test solutions across a variety of applications.

Features and Benefits:

- 200 N, 400 N, 3000 N ElectroForce linear motors capable of precisely controlling force, displacement or strain over a wide range of frequencies in tension or compression
- Multiple motor combinations offering dynamic axial, long-stroke axial and rotational motions that can be oriented and synchronized in virtually any combination.
- Interchangeable force sensors and powerful amplifiers enable force control from $\pm\,0.002$ N to $\pm3,000$ N

Configurations

- 200 N TestBench with Optional 5.6 Nm Torsion or 150 mm Stroke Extension motor
- 400 N TestBench with Optional 5.6 Nm Torsion or 150 mm Stroke Extension motor
- 3000 N TestBench
- Multi-Station configurations

Accessories

- Numerous grips and fixtures for a variety of specimen geometries
- Fluid bath and vertical mount legs available for 200 N and 400 N models
- Lower force sensors to improve data quality for low-force tests
- DMA and other software options





PLANAR BIAXIAL INSTRUMENTS


The Planar Biaxial instruments are special configurations of the TestBench series designed for precise characterization of planar samples such as membranes, textile, skin and pericardium. Multiple motors can be synchronized in various control modes and phase for almost unlimited combinations of planar loading.

Features and Benefits:


- Two or Four 200 N ElectroForce[®] linear motors arranged orthogonally for precise quasi-static or dynamic loading of planar samples
- Direct control of each axes' displacement or force plus phase enables a wide variety of test methods
- The 4-motor configuration offers the additional capability of stationary specimen center with each pair of motors acting equal-and-opposite.
- Interchangeable force sensors and precise power amplifiers enable force control from \pm 0.002 N to ± 200 N

Accessories

- Hook grips for isolating loading axes
- Fluid bath and bioreactor chamber
- Lower force sensors to improve data quality for low-force tests
- Digital Video Extensometer (DVE) for real-time strain measurements
- Interface kits for synchronizing with 3rd party DIC solutions



BioDynamic[®] Bioreactor Chamber Accessory



Specification	2-Motor Planar Biaxial	4-Motor Planar Biaxial	
Force Range	0.002 – 200 N	0.002 – 200 N	
Displacement Range	0.025 – 13 mm	0.050 – 26 mm	
Frequency Range	0.0001 – 100 Hz	0.0001 – 100 Hz	
Instrument Dimensions (HxWxD)	267 x 1000 x 1000 mm (10.5 x 39 x 39 in)	286 x 1086 x 1086 mm (10.5 x 43 x 43 in)	
Instrument Weight	102 kg (225 lbs)	110 kg (243 lbs)	

APPLICATIONS FATIGUE, DURABILITY & MATERIAL CHARACTERIZATION

Medical Devices Electronics Elasiomers Aerospace Biomaterials Composities Polymers Automotive lissue Engineering



WWW






The design of new materials and products requires a thorough assessment of material properties and complete performance evaluation within the customer's intended end-use. A variety of basic and advanced testing techniques are available to meet this need.

- Bending
- Torsion
- Shear
- Pulsatile
- Multi-axial



- Tension/Compression

- Failure Testing
- Fatigue
- Dynamic Characterization
- Creep
- Stress Relaxation
- Accelerated Life Testing








Your **Success** our **Mission**™





ACCESSORIES

ElectroForce® test instruments can be integrated with a variety of specimen fixtures, sensors, environmental chambers, fluid baths, and optional software.

Grips/Platens

Tension/Torsion Grips Wedge Grips DMA Grips Tissue Grips: • Thermal-Electrically Cooled • BioDynamic® Tensile Grips Compression Platens BioDynamic Compression Platens 3- and 4-Point Bend

Sensors


Force/Torque Displacement/Rotation Strain Pressure Acceleration Submersible Force Sensors

Fixtures and Chambers

Multi-specimen Fixture Saline Baths BioDynamic Chambers Hot/Cold Chambers Furnaces 24-well Plate Fixture T-Slot

Upgrade Options

Axial Torsion Extended Stroke Pulsatile Verical Mount System Status Indicator (SSI) Lights Battery Backup





Industry-Leading Sales & Support

RAINING

TA Instruments' leadership position results from the fact that we offer the best overall product in terms of technology, performance, quality, and customer support. While each is important, our demonstrated commitment to after-sales support is a primary reason for the continued loyalty of our customers. To provide this level of support, TA Instruments has assembled the largest worldwide team of field technical and service professionals in the industry. Others promise good service. Talk to our customers and learn how TA Instruments consistently delivers on our promise to provide exceptional service.

With direct support staff in **24 countries** and **5 continents**, TA Instruments can extend its exceptional support to you, wherever you are.

UPGRADES



SPECIFICATIONS

	3100	5500	3200	3300
Linear Motor				
Standard				
Peak/Max Sine	± 22 N	± 200 N	± 225 N	± 1000 N
Static or RMS (continuous)	± 22 N	± 140 N	± 160 N	± 700 N
High Force Option				
Peak/Max Sine	-	—	± 450 N	± 3000 N
Static or RMS (continuous)	-	-	± 320 N	± 2100 N
Displacement	5 mm	13 mm	13 mm	25 mm
Extended Stroke Option	_	_	150 mm	150 mm
Linear Velocity	0.0025 µm/s – 1.0 m/s	0.0065 µm/s – 0.80m/s	0.0065 µm/s – 3.2 m/s	0.013 µm/s – 1.5 m/s [1] 0.013 µm/s – 2.0 m/s [2]
Frequency	0.00001 Hz – 100 Hz	0.00001 Hz – 20 Hz	0.00001 Hz – 300 Hz	0.00001 Hz – 100 Hzz
Torsional Motor Option				
Standard				
Peak/Max	_	_	± 5.6 N-m	± 14 N-m ^[3] / ± 25 N-m ^[4]
Static or RMS (continuous)	—	—	± 5.6 N-m	± 14 N-m ^[3] / ± 25 N-m ^[4]
High Torque Option				
Peak/Max	—	_	—	_
Static or RMS (continuous)	-	-	_	_
Rotation	_	-	Multi-turn	Multi-turn
			62 revolutions	62 revolutions
Thermal Chamber Option	-	_	-150 to 315 °C -150 to 600 °C ^[5]	-150 to 350°C
Fluid/Saline Bath Option	Ambient to 40 °C	Ambient to 40 °C	Ambient to 40 °C	Ambient to 40 °C

Linear Motor Standard Peak/Max Sine Static or RMS (High Force O Peak/Max Sine Static or RMS (Displacemen Extended Stro

Linear Veloci

Frequency

Torsional Mo Standard Peak/Max Static or RM High Torque Peak/Max Static or RM

Rotation

Thermal Cha Fluid/Saline

— Not Available

Notes: ^[1] Linear Velocity on ElectroForce® 3310 ^[2] Linear Velocity on ElectroForce 3330 ^[3] Torque capacity on ElectroForce 3330 ^[4] Torque capacity on ElectroForce 3330 ^{5]} Temp Range for DMA 3200 with FCO Oven

Specifications are subject to change

	-			
	3510	3520/3550	200N/400N TestBench	3000N TestBench
or				
ine	± 7500 N	± 7500 N	± 200 N	± 3000 N
S (continuous)	± 5300 N	± 5300 N	± 140 N	± 2100 N
Option				
ine	_	± 15000 N	± 400 N	—
S (continuous)		± 10600 N	± 280 N	
ent	50 mm	50 mm	13 mm	25 mm
itroke Option	_		150 mm	
ocity	0.025 µm/s – 1.5 m/s	0.025 µm/s – 1.5 m/s	0.0065 µm/s – 3.2 m/s	0.013 µm/s – 2.0 m/s
	0.00001 Hz – 100 Hz	0.00001 Hz – 50 Hz	0.00001 Hz – 100 Hz	0.00001 Hz – 100 Hz
lotor Option				
x	± 49 N-m	± 49 N-m	± 5.6 N-m	_
RMS (continuous)	± 42 N-m	± 42 N-m	± 5.6 N-m	
e Option				
x	± 70 N-m	± 70 N-m	—	_
RMS (continuous)	± 50 N-m	± 50 N-m	_	_
	Multi-turn	Multi-turn	Multi-turn	
	40 revolutions	40 revolutions	40 revolutions	
hamber Option	-150 to 350 °C	-150 to 350 °C	_	_
e Bath Option	Ambient to 40 °C	_	Ambient to 40 °C	_



AMERICAS

New Castle, DE USA Lindon, UT USA Wakefield, MA USA Eden Prairie, MN USA Chicago, IL USA Costa Mesa, CA USA Montreal, Canada Toronto, Canada Mexico City, Mexico São Paulo, Brazil

EUROPE

Bochum, Germany Eschborn, Germany Wetzlar, Germany Elstree, United Kingdom Brussels, Belgium Etten-Leur, Netherlands Paris, France Barcelona, Spain Milano, Italy Warsaw, Poland Prague, Czech Republic Sollentuna, Sweden

Copenhagen, Denmark

Hüllhorst, Germany

Shanghai, China
Beijing, China
Tokyo, Japan
Seoul, South Korea
Taipei, Taiwan
Guangzhou, China
Petaling Jaya, Malaysia
Singapore
Bangalore, India
Sydney, Australia

