

# Discovery Hybrid Rheometer



## Site Preparation Guide for HR 10/20/30 Series

# Table of Contents

Table of Contents .....	2
Ideal Setup .....	3
System Components.....	4
Instrument Measurements .....	5
Utility Requirements.....	6
Power .....	6
Gas .....	7
Water .....	7
Computer Requirements .....	8–9
Hardware .....	8
Software.....	9
Temperature Systems.....	10–12
EHP, ETC, Peltier Plate/Peltier Concentric Cylinder .....	10
UHP .....	11-12
Site Preparation Checklist .....	13
TA Instrument Offices.....	14



Circulator



Power



Cooling



Gas



LN<sub>2</sub>



Fluid



Light



Hardware



Software



Temp



Lab




Customer

# Ideal Setup



## IDEAL PLACEMENT AND BENCH MEASUREMENTS

Select a location with adequate floor space and a rigid laboratory bench that is level and is in a vibration-free environment. A marble table is recommended.



Bench width: 127 cm (50 in)

Marble table width: 60 cm (24 in)

Bench depth: 76 cm (30 in)

Marble table depth: 76 cm (30 in)

Distance from the wall: 30.5 cm (12 in) min.

# System Components



## MAIN SYSTEM COMPONENTS



- A. Instrument
- B. Computer
- C. Electronics Control Module

# Instrument Measurements



## MAIN INSTRUMENT



Height: 91 cm (35.8 in)

Width: 34 cm (13.4 in)

Depth: 48 cm (18.5 in)

Weight: 42.4 kg (93.5 lbs)



## ELECTRONICS CONTROL MODULE

Height: 36 cm (14.2 in)

Width: 23 cm (9 in)

Depth: 40 cm (15.7 in)



Weight: 10.7 kg (23.6 lbs)



# Utility Requirements



## POWER

Item	Requirement
Power	<ul style="list-style-type: none"><li>• 110–240 VAC, 47–63 Hz, 1.4 kW</li><li>• Neutral to Ground (NG) voltage max 0.5 volt</li><li>• Safety ground per local regulation</li></ul>
Power cords provided	<ul style="list-style-type: none"><li>• NEMA 5-15 plug</li><li>• Type F plug</li></ul>  Type F  NEMA 5-15



**Use power cords with plugs appropriate for your circuit.**



Supply voltages lower than indicated may result in a degradation of performance.



Ensure that the mains assigned do not also supply power to noise generating equipment nearby, such as motors, welders, transformers, etc.



An independent heavy GROUND wire must be provided through the power hookup. Improper grounding may cause severe damage for which the supplier will not accept responsibility. All power strips must be fully grounded and carry the ground through to the sockets into which the computer is plugged.

# Utility Requirements



## GAS

### Air Bearing (Air or Nitrogen):

Gas Pressure	Compressed at 345–1034 kPa gauge (70–150 psig)
Flow Rate	Air flow volume capacity requirements are specific to the Temp Control Options purchased with the DHR instrument. Refer to the DHR Accessories Site Prep Guide for additional information.
Dew point	-20°C or better
Conditions	<ul style="list-style-type: none"><li>• Must be dry</li><li>• Must be free from oil and dirt<sup>1</sup></li></ul>
Other	<ul style="list-style-type: none"><li>• ¼ NPT female connection required for DHR main air supply (not provided)</li></ul>

### <sup>1</sup>Compressed Air Quality Requirements

Dew point	Ideal: -40°C Minimum: -20°C
Dirt particle	5µm
Oil including vapor	0.01 mg/m <sup>3</sup>



## WATER

Item	Requirement
Water	Fluid circulator with cooling ability for Peltier and UHP temperature systems

# Computer Requirements



## HARDWARE REQUIREMENTS

Description	Requirement
Processor	<ul style="list-style-type: none"><li>• Intel® Core™ i5 8400 or better</li><li>• 2.8 GHz with 9 MB L2 cache</li></ul>
Memory	≥ 16 GB RAM DDR4 2666 SDRAM
Hard drive	≥ 80 GB free space <ul style="list-style-type: none"><li>• 1.5 GB required for Full version of TRIOS</li><li>• 675 MB required for Lite version of TRIOS (without Online help)</li></ul>
DVD (optional)	≥ 48x CD-ROM or DVD (optional for installing TRIOS)
Screen resolution	Required: 1280 x 1024 with 24-bit colors Recommended: 1920 x 1080 with 24-bit colors
Graphic memory	128 MB
Screen (LCD) size	Required: 19" or greater Recommended: 24" wide screen
USB II port	Required with ETC and Peltier Camera Viewer options, SALS accessory, and Automatic Asphalt calibration kit
Network card	Ethernet 10Base T/100 Base TX
Additional Ethernet card(s)	Necessary if connecting the instrument directly and access is needed to the Corporate LAN. Also required for Modular Microscope Accessory.
Ethernet Cabling	10/100BaseTX Ethernet hub/switch. Must be EIA-568B Category 5+ UTP
Client-Server Protocol	DHCP
Image Capture (Camera Option)	DirectX 9.0 or higher
Second Monitor	Recommended for SALS Accessory image viewing and Modular Microscope Accessory
TCP/IP ports used	<ul style="list-style-type: none"><li>• TCP: 20010, 20011</li><li>• UDP: 5050, 5056</li></ul>



# Computer Requirements







## SOFTWARE REQUIREMENTS

Item	TRIOS
Operating System	<ul style="list-style-type: none"><li>• Windows 7, 8, 10 Ultimate &amp; Professional</li><li>• Home version not supported</li><li>• ≥ 64-bit version</li></ul>
Internet	<b>Internet connection is strongly recommended for ongoing support after installation</b>
Service Pack	Microsoft Operating System Service Pack
Updates	Windows Operating System and associated Microsoft updates must be up to date
Network	<i>A second network card for corporate connection is recommended. TA Instruments is not responsible for resolving issues associated with connections to your corporate network.</i>
Conflicts	<i>TA Instruments is not responsible for resolving hardware/software conflicts created by the addition of third-party hardware or software to the computer.</i>

# Temperature Systems

The cooling rate and minimum temperature will depend on the source of cooling.

For accessory requirements not listed, refer to the DHR Accessory Requirements guide.

Accessory	Smart Swap Requirements	
Electrically Heated Plates (EHP)		<ul style="list-style-type: none"> <li>Purge flow of 5 L/min (305 in<sup>3</sup>/min) inert gas</li> <li>Motor cooling gas flow of 10 L/min for temperatures above 250°C. Air pressure of 50–100 psig.</li> </ul>
		<ul style="list-style-type: none"> <li>Optional controlled cooling with <b>GCA</b>. Refer to the DHR Accessories Requirement guide for GCA requirements.</li> <li>Crash cooling pressure of 50–100 psig and a flow of ~2.5 scfm (70 L/min)</li> </ul>
Environmental Test Chamber (ETC)		Purge gas flow rate should be 10 L/min (610 in <sup>3</sup> /min) at 206–690 kPa (30–100 psig)
Peltier Plate/ Peltier Concentric Cylinder		<ul style="list-style-type: none"> <li>Recirculating water bath (not supplied) at 0.5 L/min (30.5 in<sup>3</sup>/min)</li> </ul>



Circulator



Power



Cooling



Gas



LN<sub>2</sub>



Fluid



Light



Hardware



Software



Temp









Lab







Customer


# Temperature Systems




Accessory	Smart Swap Requirements	
Upper Heated Plate (UHP)		<p><b><u>Option 1: Standard Cooling (temps above 10°C)</u></b></p> <ul style="list-style-type: none"> <li>• Circulation fluid: Koolance (2 bottles)</li> <li>• Fluid cooling: Supply should be 5°C below the minimum required temperature at minimum required flow rate through the system of 0.5 L/min</li> </ul>
		<p><b><u>Option 2: Standard Cooling Accessory</u></b></p> <ul style="list-style-type: none"> <li>• TA-supplied <b><u>Air-Cooled Circulator</u></b>  </li> <li>• Recommended fluid: Koolance (2 bottles)</li> </ul>
		<p><b><u>Option 3: Low Temperature Cooling Accessory–TCube Edge</u></b></p> <ul style="list-style-type: none"> <li>• TA-supplied <b><u>TCube Edge Model 5A</u></b> (PN 404500.901)</li> <li>• Fluid: Koolance (supplied with the TCube)</li> </ul> 



-   
Circulator
-   
Power
-   
Cooling
-   
Gas
-   
LN<sub>2</sub>
-   
Fluid
-   
Light
-   
Hardware
-   
Software
-   
Temp
-   
Lab
-   
Customer





# Temperature Systems

Accessory	Smart Swap Requirements	
Upper Heated Plate (UHP)	  	<p><b><u>Option 4: Low Temperature Cooling Accessory–Customer-supplied</u></b></p> <ul style="list-style-type: none"> <li>Customer-supplied refrigerated and heating circulator and appropriate fluid (ie. silicone fluid)</li> </ul> <p> <b>DO NOT USE WATER AS CIRCULATION FLUID</b></p> <ul style="list-style-type: none"> <li>Supply: 5°C below the minimum required temperature at a minimum flow rate through the system of 0.5 L/min (12.2 in<sup>3</sup>/min)</li> </ul>
		<p><b><u>Option 5: Low Temperature Cooling Accessory–Vortex</u></b></p> <ul style="list-style-type: none"> <li>TA-supplied <b><u>Vortex Cooler</u></b> (PN 545809.901)</li> <li>Air: Clean, dry, oil-free, compressed air 200L/min at 552–690 kPa gauge (80–100 psig)</li> <li>Dew point: -30°C or better</li> </ul>



Circulator   Power   Cooling   Gas   LN<sub>2</sub>   Fluid   Light   Hardware   Software   Temp   Lab   Customer



# TA Instruments Offices

For information on our latest products, contact information, and more, see our website at:  
<http://www.tainstruments.com>.

To find your local TA Instruments office and contact information, visit  
<http://www.tainstruments.com/contact/ta-directory/>

TA Instruments – Waters LLC  
Corporate Headquarters  
159 Lukens Drive  
New Castle, DE 19720  
USA

Telephone: 302-427-4000  
Fax: 302-427-4001  
Email: [info@tainstruments.com](mailto:info@tainstruments.com)