

Thermal Conductivity Instruments

DTC-300



Site Preparation Guide

Table of Contents

Table of Contents	1
Ideal Setup	2-3
System Components.....	3
Instrument Measurements	4
Utility Requirements.....	5
Power	5
Gas	6
Coolant.....	7
Water	8
Computer Requirements	9
Site Preparation Checklist	10
TA Instrument Offices.....	11

Ideal Setup



IDEAL PLACEMENT AND BENCH MEASUREMENTS

Select a location with adequate floor space and a rigid laboratory bench that is level and is in a vibration-free environment.



Bench length: 1.8 m (6 ft)

Bench depth: 76 cm (30 in)

Ideal Setup and Components



COMPONENTS



- A. Chiller
- B. Instrument
- C. Computer (Controller)



Place the chiller on the **same level** as the instrument.

Instrument Measurements



INSTRUMENT MEASUREMENTS



Height: 94 cm (37 in)

Width: 64 cm (25 in)

Depth: 71 cm (28 in)

Weight: 44 kg (98 lbs)

Utility Requirements



POWER



Instrument

- 100–120 VAC, 10 A, 60 Hz (US)
- 200–240 VAC, 10 A, 50/60 Hz (Non-US)
 - Installation of a 16A Class B or C main fuse is recommended.

Power cords provided



Use power cords with plugs appropriate for your circuit.

- NEMA 5-15 plug for 100–120 VAC
- OR
- NEMA 6-15 plug for 200–240 VAC



NEMA 6-15 plug



NEMA 5-15 plug



Supply voltages lower than indicated may result in a degradation of performance.



Ensure that the mains assigned do not also supply power to noise generating equipment nearby, such as motors, welders, transformers, etc.



An independent heavy GROUND wire must be provided through the power hookup. Improper grounding may cause severe damage for which the supplier will not accept responsibility. All power strips must be fully grounded and carry the ground through to the sockets into which the computer is plugged.

Utility Requirements



GAS

Pneumatic Ram

Gas Requirements for Pneumatic Ram	
Type	Must be nitrogen or air
Source	Must be from a gas cylinder, Grade 5 purity
Inlet Pressure	Minimum: 40 psig (2.75 bar) Maximum: 50 psig (3.45 bar)
Tubing	<ul style="list-style-type: none"> • ¼" tubing is supplied to connect from instrument to gas regulator • Tubing is rated to 100 psig (7 bar)



Purge Gas

Gas Requirements for Purge	
Type	Must be nitrogen ; at sub-ambient temperatures to prevent frost and moisture buildup on internal metal components
Inlet Pressure	Less than 5 psig needed—it is simply a low, steady flow to flush the furnace cavity
Tubing	<ul style="list-style-type: none"> • Purge gas port requires 1/8" tubing (provided) • A 1/8" stem to ¼" tubing adapter is also provided




The customer is required to supply either compressed air, dried and filtered, to 10 microns, or laboratory grade inert gas (such as nitrogen from a high pressure cylinder).

Utility Requirements



COOLANT

	Requirements
Pressure	4.3 psig (0.30 bar) maximum
Flow Rate	20 L/min (5 gal/min) maximum  Operating close to or below min. pressure will result in erratic operation.
Cooling Capacity	<ul style="list-style-type: none"> At -20°C = 265W At 0°C = 650W At +20°C = 1000W
Coolant Liquid	<ul style="list-style-type: none"> 50/50 mix of DI water and laboratory grade ethylene glycol for testing up to 300°C Denaturalized alcohol for sub-zero testing
Connections	Supplied: Two hoses (6 feet / 1.8 meters each) to connect to the Coolant Inlet and Outlet ports on the back of the furnace. The other end of each hose is blank; a male 3/8" barb with 1/4" threads (provided) must be installed to connect to the coolant source.



Coolant Inlet and Outlet ports and hoses

P/N 202612.001 (120V for 60Hz power)
 or
 202442.001 (240V for 50Hz power)



Place the chiller on the **same level** as the instrument.

Utility Requirements



WATER



Excessively cold water will result in “sweating” and corrosion of cooled metal surfaces. The purge of nitrogen gas through the furnace will prevent this.



Operating below the minimum chiller coolant level will result in erratic operation. A chiller/circulator is recommended for this instrument.



If a chiller/circulator is to be used, it **must** be placed at the same level as the instrument.

Computer Requirements



HARDWARE REQUIREMENTS

- Unused RS-232 (serial) port
- Unused USB port



Instrument drivers and software are provided on a CD.



Computer should not be attached to other analytical instruments or LAN.







SOFTWARE REQUIREMENTS

Item	Requirement
Operating System	<ul style="list-style-type: none">• Windows 7 or 10, 32- or 64-bit, Ultimate, Enterprise & Professional• Home version not supported
Network	<ul style="list-style-type: none">• <i>TA Instruments is not responsible for resolving issues associated with connections to your corporate network.</i>• <i>Network cards and/or certain network operation frequently interfere with the operation of the instrument and software.</i>
Conflicts	<i>TA Instruments is not responsible for resolving hardware/software conflicts created by the addition of third party hardware or software to the computer.</i>

Site Preparation Checklist



Thermal Conductivity Instruments: DTC-300

	<p>Sufficient bench space for instrument, computer, and chiller (if needed):</p> <ul style="list-style-type: none"> <input type="checkbox"/> Length: 1.8 m (6 ft) <input type="checkbox"/> Depth: 76 cm (30 in)
	<p>Instrument power is</p> <ul style="list-style-type: none"> <input type="checkbox"/> 100–120 VAC, 10 A max, 60 Hz (US) <input type="checkbox"/> 220–240 VAC, 10 A max, 50/60 Hz (Non-US) <p>Computer, monitor, and chiller power is</p> <ul style="list-style-type: none"> <input type="checkbox"/> 120 V (US) <input type="checkbox"/> 220–240 V, 6.4 A max, 50/60 Hz (Non-US)
	<p>Pneumatic Ram Gas – Air or Nitrogen</p> <ul style="list-style-type: none"> <input type="checkbox"/> Regulator to allow 40–50 psig (2.75–5.50 bar) <p>Purge Gas – Dry Nitrogen</p> <ul style="list-style-type: none"> <input type="checkbox"/> Under 5 psig (2.75–5.50 bar)
	<p>Water Circulation</p> <ul style="list-style-type: none"> <input type="checkbox"/> Nominal flow rate of 20 L/minute (5 gal/min) <input type="checkbox"/> Maximum pressure of 4.3 psig (0.30 bar) <input type="checkbox"/> Sufficient cooling capacity <input type="checkbox"/> 50/50 mix of DI water and laboratory grade ethylene glycol (for temperatures up to 300°C) <input type="checkbox"/> Denaturalized alcohol (for sub-zero measurements)

I hereby acknowledge that all utility requirements have been met per the checklist above and that they will be ready at the agreed time of installation.

If all utility requirements are not met at the agreed time of installation, additional charges may be incurred for a return Service trip.

_____ / /
Customer *DD* *MM* *YYYY*

_____ / / /
Company *City* *State* *Country*

Please send a signed copy of the completed checklist to your local Service representative.

TA Instruments Offices

For information on our latest products, contact information, and more, see our website at:
<http://www.tainstruments.com>.

To find your local TA Instruments office and contact information, visit
<http://www.tainstruments.com/contact/ta-directory/>

TA Instruments – Waters LLC
Corporate Headquarters
159 Lukens Drive
New Castle, DE 19720
USA

Telephone: 302-427-4000

Fax: 302-427-4001

Email: info@tainstruments.com