
HEATING MICROSCOPE - HSM

INSTANTANEOUS HEATING

- MELTING
- SOFTENING
- SPHERE POINT
- HALF SPHERE POINT
- FUSIBILITY
- BLOATING
- CONTACT ANGLE
- SURFACE TENSION
- THEORETICAL VISCOSITY

Misura® HSM

A SINGLE INVESTMENT FOR A COMPLETE CHARACTERIZATION OF MATERIALS

DOUBLE, TRIPLE OR QUADRUPLE THE NUMBER OF TESTS PER DAY

ANALYSIS WITH INDUSTRIAL OR INSTANTANEOUS FIRING CYCLES

TIME SAVING IN THE SAMPLE PREPARATION

Heating Microscopes Misura® HSM are the result of more than twenty years of research and development: improved and developed to satisfy customers' needs, they have characteristics that make possible their use out of classical boundaries of the heating microscopy. Misura® Heating Microscopes can simulate the fast or instantaneous firing cycles, in order to get a more complete characterization of the evaluated materials. The analyses can be carried out on one, two, up to four samples at the same time, to identify automatically:

- characteristic temperatures: Sintering Beginning, Softening, Sphere, Half Sphere, Melting;
- fusibility according to International Standards;
- flattening curve;
- sintering rate of the materials;
- contact angle curve;
- sample area variation curve;
- ratio curve between base and height;
- bloating effects;
- combustion;
- theoretical viscosity;
- surface tension at high temperatures.

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

IMAGES ACQUISITION, TEMPERATURE OR TIME BASED

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TECHNICAL CHARACTERISTICS

HSM


Optical measuring system	single beam
Expansion/ Contraction range (according to sample dimensions)	+ 50% /- 100%
Sample displacement	vertical
Sample number	1 or 2
Standard sample dimensions:	
- cylinder (height/diam.)	3 x 2 mm
International Standards	ASTM D1857-68, BS 1016:Part 15:1960, CEN/TR 15404:2010, CEN/TS 15370-1:2006, DIN 51730,DM 05-02-1998, IS 12891:1990, ISO 540:1995, NF M03-048
Temperature range on specimen	from RT up to 1650 °C
Max heating rate per minute	0,1- 80 °C/min
Sample thermocouple	S or B type
Kiln thermocouple (models from 1600 °C)	B type
PC interface	USB and RS-232
Software	Misura® thermal analysis

CONFIGURATION OPTIONS

Optional sample number	up to 4 samples simultaneously
Optional sample dimensions:	
- cylinder (height/diam.)	up to 20 x 20 mm
- cube	up to 20 x 20 x 20 mm
- pyramid (height/side - angle)	up to 19 x 6,4 mm - 60°
Granule (microns)	100-800
Instantaneous heating (Flash)	350 °C/min
DTA:	optional
- sample type	pressed powders (no crucible), bulk (no crucible), powders (crucible)
Atmosphere:	
- oxidant	static or dynamic
- inert gas	dynamic
Cooling	optional


The advertising material is not valid for contractual purposes and technical specifications could be subjected to changes

COMPARISON BETWEEN FUSIBILITY ANALYSES OF ASHES CARRIED OUT WITH DIFFERENT HEATING RATES

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FLATTENING CURVE AND CONTACT ANGLE

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