# Materials Characterization by Thermal Analysis (DSC & TGA), Rheology, and Dynamic Mechanical Analysis (Part 2)

Sarah Cotts, Charles Potter & Don DiPietro Applications Scientists & Sales Representative





# Afternoon Agenda: Rheology and DMA

#### Introduction to Rheology

- Fundamentals of Stress, Strain and Shear Rate
- Instrumentation
- Geometries
- Flow Testing
  - Shear Thinning and Shear Thickening
  - Thixotropy
  - Yield Stress
- Oscillation Testing
  - Viscoelasticity
  - "Hands On" demonstration
  - Characterizing Thermoplastics

#### Thermosets and Gels

- Viscoelastic properties of Gels
- Characterizing Gelation
- Experimental considerations
- Specialized Accessories
- Dynamic Mechanical Analysis
  - Characterizing Glass Transition
  - Comparison of Deformation Modes
- •Case Study:
  - Rubber Characterization using DSC, TGA and DMA.



What is Rheology?

# Rheology is the study of **flow** and **deformation**.









© TA Instruments

# Rheology: the study of flow and deformation

#### Viscosity

- Non-Newtonian Viscosity
  - Shear thinning
  - Shear thickening
  - Thixotropy
  - Yield Stress
- Viscosity under processing conditions

#### Modulus

- Measure viscoelastic properties
  - Storage Modulus
  - Loss Modulus
  - Tan Delta
- Changes with time, temperature







#### **Discovery Hybrid Rheometer Technology**



## **ARES G2: Separate Motor and Transducer**



Torque Transducer maintains the null position as the sample is deformed.

Sample torque to be measured directly, without contributions from motor friction or inertia.

Normal Force Transducer provides highly accurate normal force measurements

High stiffness for precise gap control.

Direct Drive motor applies accurate and precise rotational deformation, without contributing to measured torque.



#### **Torsion Flow in Parallel Plates**



 $\begin{aligned} r &= \text{plate radius} \\ h &= \text{distance between plates} \\ M &= \text{torque } (\mu \text{N.m}) \\ \theta &= \text{Angular motor deflection (radians)} \\ \Omega &= \text{Motor angular velocity (rad/s)} \end{aligned}$ 

Stress ( $\sigma$ ) $\sigma = \frac{2}{\pi r^3} \times M$ Strain ( $\gamma$ ) $\gamma = \frac{r}{h} \times \theta$ Strain rate ( $\dot{\gamma}$ )  $\dot{\gamma} = \frac{r}{h} \times \Omega$ 



# **Geometry Options**





# Choosing a Geometry Size



- Assess the 'viscosity' of your sample
- Select diameter appropriate for viscosity of sample
  - Low viscosity (milk) 60mm geometry
  - Medium viscosity (honey) 40mm geometry
  - High viscosity (caramel) 20 or 25mm geometry
- Examine data in terms of torque/displacement/speed and modify geometry choice to move into optimum working range



# **Correct Sample Loading**





### Shear Rate varies across a Parallel Plate

 For a given angle of deformation, there is a greater arc of deformation at the edge of the plate than at the center



 $\gamma = \frac{dx}{h}$  dx increases further from the center, h stays constant

Single-point correction for the parallel plate geometry (0.76 radius) [M.S. Carvalho, M. Padmanabhan and C.W. Macosko, *J. Rheol.* 38 (1994) 1925-1936]



© TA Instruments

# Shear Rate is Uniform across a Cone

The cone shape produces a smaller gap height closer to the center, so the shear on the sample is constant



$$\gamma = \frac{dx}{h}$$

h increases proportionally to dx,  $\gamma$  is uniform



# Limitations of Cone and Plate



Gap must be > or = 10 [particle size]!!



#### What Geometry should we use?





# Should have used a parallel plate!





#### When to Use Concentric Cylinders



**Peltier Concentric Cylinder** 

- Low to Medium Viscosity Liquids
- Unstable Dispersions and Slurries
- Minimize Effects of Evaporation
- Easy Sample Loading
- Weakly Structured Samples (Vane)
- Low Shear Rates



# **Geometry Summary**

Parallel Plates	Cone and Plate	Concentric Cylinders
1-2 mL	50-500 μL	7-25 mL
Liquids, gels, soft solids, dispersions, etc.	Unfilled liquids, isothermal tests	"Pourable" liquids, low viscosities, dispersions
Used for all samples. Roughened surfaces available to prevent slip.	Most accurate measurement of non- Newtonian Viscosity, small sample volume.	Least effected by sample loading technique or evaporation.



## **Flow Testing**



## Newtonian and Non-Newtonian Fluids

 Newtonian Fluids - Viscosity does not change with changes in shear rate or time.

(examples: water, oil, honey)

 Non-Newtonian Fluids - Viscosity is time or shear rate dependent

(examples: mayonnaise, paint, polymer, asphalt)

- Shear Thinning: viscosity decreases as shear rate increases
- Shear Thickening: viscosity increases as shear rate increases.



#### Viscosity of Honey and Mayonnaise



#### Viscosity is curve, not a single value!



## **Viscosity Flow Curve**



shear rate (1/s)



## **General Viscosity Curve for Suspensions**



 $\log \dot{\gamma}$ 

Reference:Barnes, H.A., Hutton, J.F., and Walters, K., <u>An Introduction to Rheology</u>, Elsevier Science B.V., 1989. ISBN 0-444-87469-0



## Flow Curve for Suspensions



TAINSTRUMENTS.COM

#### Viscosity vs. Shear Rate- Shampoo



#### Viscosity vs. Shear Rate- Shampoo



# Viscosity of Water: effect of geometry



# Wall Slip

#### No slip condition

Ideal, Assumed Velocity Profile



- Can manifest in many ways:
  - Apparent Double Yielding
  - Low Yield Stress
  - Low Viscosity
  - High tan δ
- Test by running the same experiment at different gaps
- For samples that don't slip, the results will be independent of the gap
- Can happen with any geometry: plate, cone, concentric cylinder, etc...

#### Wall Slip Incorrect Velocity Profile



Yield Stress Measurements on Toothpaste





# **Edge Failure**



- Results in apparent drop in viscosity
- Remedy: decrease gap to increase stabilizing influence of surface tension



#### Shear Thinning or Sample Instability?



# Thixotropy

The thixotropy characterizes the time dependence of reversible structure changes in complex fluids. The control of thixotropy is important to control:

- process conditions; for example, to avoid structure build up in pipes during rest periods
- sagging and leveling; gloss of paints and coatings





## **Thixotropy of Paint**



© TA Instruments

## Thixotropic Loop-Mayonnaise



#### **Thixotropic Material**

- Up and down rate ramps do not superimpose
- Area under the curve is a measure of thixotropy



# What is Yield Stress?

- Some Structured Fluids behave like "solids" at rest.
- A critical stress must be applied for these materials to flow.
- •What does Yield Stress do?
  - Stabilize against sedimentation of separation
  - Improve ease of use
  - Prevent dispensing of product





#### **Yield Stress- Ketchup**



### **Paint Formulations- Yield Stress**




#### **Yield Stress: Orange Juice**



#### **Yield Stresses of Protein Solutions**



#### **Protein Solutions: Low Torque Sensitivity**



# **Understanding Torque**

- •What is a **N.m**?
  - An apple (about 150 g) on the end of a meter stick
  - 1.5 N.m
- •What is a **µN.m**?
  - A grain of salt (about 0.5 mg) on the end of a meter stick)
  - **5 μN.m**

#### •And a **nanoN.m**?

- A speck of dust (1 µg) on the end of a meter stick
- 10 nanoN.m





# Oscillation



# What is Oscillation?



Dynamic stress applied sinusoidally User-defined Stress or Strain amplitude and frequency



# **Oscillation Testing**

- An oscillatory (sinusoidal) deformation is applied to a sample.
- The material response (stress) is measured.
- The phase angle  $\delta$ , or phase shift, between the deformation and response is measured.





#### Oscillation Testing: Response for solids and liquids





#### **Oscillation Testing: Viscoelastic Material**





# **Viscoelastic Parameters**

<u>The Modulus:</u> Measure of materials overall resistance to deformation.

<u>The Elastic (Storage) Modulus:</u> Measure of elasticity of material. The ability of the material to store energy.

<u>The Viscous (loss) Modulus:</u> The ability of the material to dissipate energy. Energy lost as heat.

#### <u>Tan Delta:</u>

Measure of material damping such as vibration or sound damping.

$$G^* = \left(\frac{\text{Stress}^*}{\text{Strain}}\right)$$

$$G' = \left(\frac{\text{Stress}^*}{\text{Strain}}\right)\cos\delta$$

$$G'' = \left(\frac{Stress^*}{Strain}\right) \sin \delta$$

$$\tan \delta = \left(\frac{G'}{G'}\right)$$



#### Storage and Loss of a Viscoelastic Material





# **Viscoelasticity Defined**

Range of Material Behavior Liquid Like----- Solid Like Ideal Fluid ----- Most Materials -----Ideal Solid Purely Viscous ----- Viscoelastic ----- Purely Elastic

# *Viscoelasticity*: Having both viscous and elastic properties

 Materials behave in the linear manner, as described by Hooke and Newton, only on a small scale in stress or deformation.



# **Oscillation Testing- Linear Viscoelastic Region**



Strain (amplitude)



#### **Time-Dependent Viscoelastic Behavior**



- •Short deformation time: pitch behaves like a solid
- •Long deformation time: pitch behaves like a highly viscous liquid
  - 9<sup>th</sup> drop fell July 2013



Started in 1927 by Thomas Parnell in Queensland, Australia

http://www.theatlantic.com/technology/archive/2013/07/the-3-most-exciting-words-in-science-right-now-the-pitch-dropped/277919/



#### **Time-Dependent Viscoelastic Behavior**





# **Viscoelastic Putties**

•Compare with your neighbors!

•How do these differ in their rheology?

- Viscosity
- Modulus
- Phase angle
- Relaxation Time
- •How can you test it?
  - "Manual Rheology"
  - With a Rheometer



## Influence of Molecular Weight on G' and G"

The intersection of G' and G" shifts to lower frequency as MW increases.



INSTRUMENTS.COM

# Influence of MWD on G' and G"



FAINSTRUMENTS.COM

# Influence of MW on Viscosity

The zero shear viscosity increases with increasing molecular weight. TTS is applied to obtain the extended frequency range.



AINSTRUMENTS.COM



## Melt Rheology: MW Effect on Zero Shear Viscosity

- Sensitive to Molecular Weight, MW
- For Low MW (no Entanglements)  $\eta_0$  is proportional to MW
- For MW > Critical MW<sub>c</sub>,  $\eta_0$  is proportional to MW<sup>3.4</sup>



Ref. Graessley, Physical Properties of Polymers, ACS, c 1984.



#### Idealized Flow Curve – Polymer Melts





## Viscosity Measurements of LDPE at 190°C



# **High MW Contributions**



Macosko, TA Instruments Users' Meeting, 2015



## **Example: Surface Defects during Pipe Extrusion**



# **Tack and Peel of Adhesives**



- Bond strength is obtained from peel (fast) and tack (slow) tests
- Tack and Peel are a function of viscoelastic properties at different frequencies

Tack and Peel performance of a PSA

## **Thermosets and Gels**



# What is a Gel?

 A soft solid that contains a polymeric network and a substantial fraction of solvent

- Latin: gelatus (frozen; immobile)
- •"A substantially dilute crosslinked system that exhibits no flow in the steady state."
  - J.D. Ferry, Viscoelastic Properties of Polymers, 1980.
- •Chemical Gel: covalent network
- Physical Gel: non-covalent network
- Hydrogel: network with significant water content



# Gelation





# Physical, Chemical and Biological Hydrogels



E.S. Place, J.H. George, C.K. Williams, M.M. Stevens, Chem. Soc. Rev. 2009, 38, 1139–1151



# Viscoelastic Properties: "Is it a Gel?"





Y. He, P. G. Boswell, P. Bühlmann, T. P. Lodge, J. Phys. Chem. B, <u>111</u>, 4645, (2007) CTA

## **Sample Preparation**



## Isothermal Gelation: Hyaluronic Acid Gels

- Hyaluronic acid gels are used as lubricating agent during abdominal surgeries to prevent adhesion and also for join lubrication, wound healing etc.
- Rheology can monitor gelation and evaluate gel





Complex viscosity n\*

(Pa.s)

10

101

100

10

10-2

10-3

0

(Pa) 🔵

Loss modulus G"

Storage modulus G' (Pa) 🔘

# Rapid Gelation of 2-component system

- Experimental Challenge:
  - 2 components form a gel upon contact
  - Used to create small encapsulating beads and drug delivery devices.
- Upon mixing, gelation occurs in seconds and cannot be easily loaded.
- Decreased temperature to delay gelation.







#### **Characterizing short Gel Time**



© TA Instruments

#### **Thermo-Reversible Gelation: Gelatin**



# **UV-Cured Hydrogels**

- Widely used systems
  - 3D curing
  - Photolithography
  - On-demand gelation
- UV accessory attachment for Rheometer
  - Allows for in-situ gelation and optimum measurement of modulus.
  - Characterize gelation kinetics with controlled:
    - irradiance
    - wavelength
    - exposure time
    - temperature








## UV Gelation, 365 nm LED



## **Powder Coating: Cure Test**



## **Cure Testing- Dimensional Change**



# **Change in Mechanical Properties During Drying**

- Relative Humidity and Temperature Controlled Chamber
- Quantitative measurement of modulus during drying of the bulk material.
- Characterize time of drying
  - Time needed to "set" (Crossover point)
- Determine conditions needed to achieve drying
- Test Method: constant temperature and humidity







## **Example: Glue Drying under Different Humidity**





#### **Moisture-Cured System- Humidity Control**



# **Dynamic Mechanical Analyzers**



# Is DMA Thermal Analysis or Rheology?

#### Thermal Analysis

measurement as a *function of temperature or time*.

#### Rheology

- the science of *stress* and *deformation* of matter.
- DMA mechanically deforms a sample and measures the sample response. The response to the deformation can be monitored as a function of temperature or time.



### **DMAs from TA Instruments**







**Q800** 



## **DMAs from TA Instruments**

RSA G2 Separate Motor & Transducer

#### DMA850 and Q800

#### **Combined Motor & Transducer**





## DMA Mode on DHR and ARES-G2



## DMA Q800: Schematic





## **RSA G2 Schematic: Dual Head Design**





## **DMA Viscoelastic Parameters**

<u>The Modulus:</u> Measure of materials overall resistance to deformation.

<u>The Elastic (Storage) Modulus:</u> Measure of elasticity of material. The ability of the material to store energy.

<u>The Viscous (loss) Modulus:</u> The ability of the material to dissipate energy. Energy lost as heat.

#### Tan Delta:

Measure of material damping - such as vibration or sound damping.

$$E^* = \left(\frac{Stress^*}{Strain}\right)$$

$$\mathbf{E'} = \left(\frac{\mathbf{Stress}^*}{\mathbf{Strain}}\right)\cos\delta$$

$$E'' = \left(\frac{Stress^*}{Strain}\right) \sin \delta$$

$$\tan \delta = \left(\frac{E'}{E'}\right)$$



# **Dynamic Temperature Ramp**



#### Temperature



# **Dynamic Mechanical Analysis of Glass Transition**

- •DMA more sensitive to Tg than DSC, directly measures changes in mechanical and viscoelastic properties as a function of temperature.
- Materials whose glass transitions cannot be resolved by DSC can often be measured easily in DMA
  - Semi-crystalline materials with low amorphous content
  - Composites in which the polymer weight fraction is small
  - Glass Transitions that occur over a wide range, or overlap with other thermal events
- Glass Transition measurement by DMA particularly relevant to characterizing materials for their end-use properties



# Glass Transition E' Onset, E'' Peak, and Tan $\delta$ Peak

#### Storage Modulus E' Onset:

Occurs at lowest temperature, relates to mechanical failure

#### Loss Modulus E" Peak:

- Occurs at middle temperature
- Related to the physical property changes
- Reflects molecular processes the temperature at the onset of segmental motion

#### Tan Delta Peak:

- Occurs at highest temperature; Used historically in literature
- Measure of the "leatherlike" midpoint between the glassy and rubbery states
- Height and shape change systematically with amorphous content.

Turi, Edith, A, Thermal Characterization of Polymeric Materials, Second Edition, Volume I., Academic Press, Brooklyn, New York, P. 980.



#### **Temperature Ramp- Glass Transition**



## PET Film: Effect of Frequency on Tg

• PET film tested at 0.1 Hz, 1Hz and 10 Hz



#### Glass Transition is a Range, not a Temperature



## Crystallinity, Molecular Weight, and Crosslinking





# **DMA Deformation Modes**



# **Tension DMA**

#### • Young's Modulus (E)

 Easy to adjust clamps to accommodate different samples. Allows for thermal expansion or shrinkage.

#### •Dimensions:

- Length can be adjusted directly, measured by instrument.
- Thickness up to 2 mm

#### •Materials:

- Polymer films (Mylar, Kapton)
- Elastomers (o-rings, seals)
- Free films of coatings (dried paint)
- Fibers, bundled or single.





## **Storage Modulus of PET Fiber- Draw Ratios**



Murayama, Takayuki. "Dynamic Mechanical Analysis of Polymeric Material." Elsevier Scientific, 1978. pp. 80.



© TA Instruments

### **Glass Transition of EPDM- DHR DMA Tension**



# **Effect of Solvent**



#### Nylon 6: Effect of Humidity on Glass Transition



#### **Error: Film Sample not Loaded Flat**





Instrument: Q800 Clamp: tension Temperature: 0 ℃ to 180 ℃ Heating rate: 3 ℃/min Frequency: 1Hz Amplitude: 10 µm



### **PET Film Measured at X and Z Direction**





#### Iso-force Temp Ramp: measure shrinkage





## Iso-strain Temp Ramp: measure shrinking force





# **Compression DMA**

- Compression Modulus (E)
- Soft materials with high elasticity.
- Must be compressible, without yielding under deformation.
- •Dimensions:
  - Ideally cut to diameter of the plates. Can also accommodate smaller disks or rectangles.
  - 1-10 mm thick
- •Materials:
  - Foams (mattress, packaging, anti-vibration)
  - Soft Elastomers (above Tg only!)
  - Stiff hydrogels, biological tissue





#### Foam Compression DMA: Temperature Ramp Rate



#### **Effects of Humidity on Glass Transition of Foam**





# **Bending DMA**

#### •Flexural Modulus (E)

•3 Point Bend (unclamped) and Cantilever (clamped)

#### Dimensions

- Fixed lengths: (i.e. 40, 25 and 10 mm 3PB)
- Width up to 12.5 mm
- Thickness ideally less than 1/10 length.

#### Materials

- Unfilled thermoplastics (Cantilever only > Tg)
- Elastomers (Cantilever)
- Thermosets (3PB)
- Composites (3PB)
- Metals (3PB)





## Primary and Secondary Transitions in PC




## What Causes E' Increase after Tg?

- Sample sagging after Tg
- Solution: use cantilever clamp instead of 3-p bending





Instrument: RSA G2 Clamp: 3-p bending Temperature: 50 °C to 180 °C Heating rate: 3 °C/min Frequency: 1Hz Amplitude: 10 µm



## Fiber Reinforced Polymer- 3 Point Bending



## **Torsion DMA**

- •Shear Modulus (G) "Modulus of Rigidity"
- Ideal for very high modulus materials; accommodates wide range of dimensions.

#### •Dimensions:

- Small: 7 mm long, 3 mm wide, 0.5 mm thick
- Large: 40 mm long, 12.5 mm wide, 4 mm thick
- Cylinder: 1.5, 3 or 4.5 mm diameter

#### Materials:

- Thermoplastics and Thermosets
- Elastomers
- Composites
- Metals





### **Fiber Reinforced Polymer- Torsion**



### **Torsion vs. 3 Point Bending**



# Parallel Plate DMA

#### •Shear Modulus (G)

• Full range of viscoelastic behavior (glassy, rubbery and terminal region).

#### •Dimensions:

- 25, 8 or 4 mm parallel plates
- •0.5 3 mm gap thickness

#### Materials:

- Thermoplastics: load above softening point, ramp temperature down.
- Thermosets: cure in place on disposable plates.
- Elastomers: cut disk and glue to plates.
- Adhesives: too soft to test with linear DMA





#### Hot Melt Adhesive: Parallel Plates



### **Rubber DMA: Parallel Plates**



## Thermal Analysis Techniques for Rubber Characterization

"Happy and Sad Balls" TGA, DSC and DMA









## What are Happy and Sad Balls?

- Set of 2 black rubber balls used as demonstration of viscoelastic behavior.
  - When dropped, the Happy Ball bounces and the Sad Ball does not.
- Quality varies between vendors!
- Thermal Analysis can help understand why these materials perform differently.





## Amazon Reviews, "Vendor 1"





## Frequency Sweep: Tan Delta (RSA G2)



## **Tensile Testing (RSA G2)**



© TA Instruments

## **Physical Properties**

	Vend	dor 1	Vendor 2		
	Happy Ball	Sad Ball	Happy Ball	Sad Ball	
Shore A Hardness	60	63	44	42	
Density (g/cm <sup>3</sup> )	1.53	1.39	1.09	1.54	
Tensile Strength (MPa)	2.2	2.7	2.2	3.5	
Elongation at Break	320%	250%	190%	1,290%	
Tan δ (damping) at 50 Hz	0.127	0.232	0.046	1.158	

TAINSTRUMENTS COM

What is the source of the difference in physical properties?

- Chemistry of elastomer?
- Degree of crosslinking?
- Amount of additives or filler?

### FT-IR: Useful, but not conclusive

Vendor 1

#### Vendor 2





## **Thermal Analysis Characterization Techniques**

#### **1.** Thermo-Gravimetric Analysis (TGA)

- Measure weight change while increasing temperature.
- Quantify the content of volatiles, elastomer, carbon black and inorganic filler (ash).

#### 2. Differential Scanning Calorimetry (DSC)

- Measure Heat Flow as a function of temperature.
- Identify glass transition, melting/ crystallization, and any polymerization or curing.

#### 3. Dynamic Mechanical Analysis (DMA)

- Measure viscoelastic properties as a function of temperature.
- More sensitive characterization of glass transition.









### **TGA: Vendor 1, filler**



### **TGA: Vendor 2, filler**



## DSC: Vendor 1



## DSC: Vendor 2



### **DMA: Vendor 1**



#### **DMA: Vendor 2**



### Comparison of Glass Transitions, Vendors 1 & 2



© TA Instruments

# Summary of Results

•TGA

- Quantified the basic components of the compound.
- Saw higher filler loading in rubber from Vendor 1.
- Other applications:
  - Use evolved gas analysis (mass spec, FT-IR) to identify components separately

•DSC

- •Ruled out extent of cure as a root cause for performance problems.
- Saw differences in glass transition, melting in rubber from Vendor 1.

• DMA

- Measured Tan Delta, quantifying the performance difference.
- Observed glass transition in detail, learned the Vendor 1 Sad Ball is a blend of elastomers.
- Other applications:
  - Creep-Recovery and Stress Relaxation, to mimic response in end use.



# **Rheology Applications Resources**

How do I get started using my new rheometer?	•	Quick Start Guides (e-Training) video "Getting Started Guide" (desktop shortcut)
How do I load my sample?	•	TA Tech Tips
What kind of test should I use on my sample?	•	Applications Notes "Practical Approach to Rheology"
What parameters should I use in my procedure?	•	TRIOS Help files TA Tech Tips rheologysupport@tainstruments.com
What do these measurements tell me about my sample?	•	Webinar: "Essential Tools for the New Rheologist," Neil Cunningham
How can I be confident in the quality of my data?	•	<u>"Strategies for Better Data"</u> <u>Webinar:</u> "Experimental Challenges of Shear Rheology; How to Avoid Bad Data" Randy Ewoldt
How do I analyze my data in TRIOS?	•	Quick Start Guide: TRIOS Analysis TRIOS Help files



## **Trios Online Help Manual**





## **Instructional Videos**

From <u>www.tainstruments.com</u> click on Videos, Support or Training



Select Videos for TA Tech Tips, Webinars and Quick Start Courses



See also: <u>https://www.youtube.com/user/TATechTips</u>



## **Instructional Video Resources**

#### **Quickstart e-Training Courses**

		Web based e-Training Courses		
April - State Andread - State Andread - State -	DMA Q800 Quickstart Course – Instrument and Experimental Setup	TA Instruments offers a variety of training opportunities via the Internet. e-Training opportunities include the following: QUICKSTART e-TRAINING COURSES		
	DMA Q800 – Analysis Quickstart	QuickStart e-Training courses are designed to teach a new user how to set up and run samples on their analyzers. These 60-90 minute courses are available whenever you are. These pre-recorded courses are available to anyone at no charge. Typically these courses should be attended shortly after installation. Contact Us for Web based e-Training Courses		
		Search Results for "TRIOS"		
Universal Analysis QuickStart Course	Discovery DSC – TRIOS Data Analysis			
	U.Y	TRIOS Guardian – a tool to aid in 21 CFR 11 compliance		
Universal Analysis Advanced E-Training				





Universal Analysis Custom Report



TRIOS - Analysis Reports

TRIOS – Analysis in Overlay



- San alama dagi sebas Pananana anda Sanaharan dari sahaga

TRIOS – Cox Merz Transformation



## Help with TA Instruments

- Check the manuals, TRIOS help.
- Contact the Applications Helpline
  - <u>rheologysupport@tainstruments.com</u>
  - thermalsupport@tainstruments.com
- Contact the Service Helpline
  - servicehelpline@tainstruments.com
- Call your local Technical or Service Representative
- Visit our Website <u>www.tainstruments.com</u> for training videos, TA Tech tips, application notes and much more!



## **Thank You**

The World Leader in Thermal Analysis, Rheology, and Microcalorimetry

