

## Thermal Analysis and Stability of Biomaterials

Louis Waguespack  
Applications Laboratory Manager  
TA Instruments – Waters LLC  
[lwaquespack@tainstruments.com](mailto:lwaquespack@tainstruments.com)



©2017 TA Instruments

## Agenda

- **9 - 10:15am: Thermal analysis and stability of biomaterials**
- 10:15-10:30am: Break
- 10:30-noon: Using rheology to characterize flow and viscoelastic properties of hydrogels, adhesives and biopolymers
- 12-1pm: Lunch
- 1-2:15pm: Mechanical testing of medical devices
- 2:15-2:30 - break
- 2:30-3:15pm: Mechanical testing of engineered tissues and biomaterials
- 3:15-4pm – Q&A with TA Instruments Applications Engineers

## What is a Biomaterial?

- Any synthetic material used to replace or restore function to a body tissue
  - Biomaterials are continuously in contact with body fluids
  - Exposure to body fluids implies that the biomaterial is placed within the interior of the body
    - This places restrictions on the types of biomaterials used




3

 TAINSTRUMENTS.COM

## Ideal Biomaterial Properties

- Ideal material combination for biomaterials exhibits these qualities
  - A biocompatible chemical composition to avoid negative tissue reactions
  - Resistance to degradation
    - Corrosion resistance for metals
    - Resistance to biological degradation in polymers
  - Acceptable strength to sustain cyclic loading by the joint
  - A low modulus to minimize bone resorption
  - High wear resistance to minimize wear-debris generation



4

 TAINSTRUMENTS.COM

## Types of Biomaterials

### • Metals

- Select metals are biocompatible and capable of long-term success as bodily implants
  - Nontoxic elements are used as alloying elements in biomedical alloys
  - Principal metals used: Stainless Steels, Pure titanium and titanium alloys, and Cobalt-base alloys
    - Ability to bear significant loads, withstand fatigue loading, and undergo plastic deformation prior to failure
- Orthopedic surgeries commonly involve metallic implantation
  - Metallic implants are used in maxillofacial surgery, cardiovascular surgery and dental materials




5

 TAINSTRUMENTS.COM

## Types of Biomaterials

### • Polymers

- Used for: facial prostheses, tracheal tubes, heart components, kidney components, liver components, dentures, hip and knees joints
- Medical adhesives
- Coatings and sealants



Application	Polymer Used
Knee, Hip, & Shoulder Joints	Ultrahigh molecular weight polyethylene (UHMWPE)
Finger Joints	Silicone
Sutures	Polylactic and Polyglycolic acid, Nylon
Tracheal Tubes	Silicone, Acrylic, Nylon
Heart Pacemaker	Acetal, Polyethylene, Polyurethane
Blood Vessels	Polyester, Polytetrafluoroethylene, PVC
Gastrointestinal Segments	Nylon, PVC, Silicones
Facial Prostheses	Polydimethyl siloxane, Polyurethane, PVC
Bone Cement	Polymethyl methacrylate

6

 TAINSTRUMENTS.COM

## Biopolymers

Polymer	Typical Use
Polyethylene	Tubing, Connectors and Bottles, Plastic Surgery Implants
Polypropylene	Disposable Syringes, Connectors, Finger-joint Prostheses, Nonabsorbable Sutures
Polytetrafluoroethylene, Polydifluoromethylene	Vascular Graft Prostheses, Heart Patches, Retinal Detachment Treatment
Polystyrene	Disposable Laboratory Items
Polymethyl methacrylate	Bone Cement, Artificial and Implanted Teeth, Denture Materials and Fillings, Intraocular Lens
Polyvinyl chloride	Disposable Medical Items, Blood Tubing Line, Cardiac Catheters, Artificial Limb Materials
Polydimethylsiloxanes (Silicone rubber)	Plastic Surgery Implants, Artificial heart and Heart-assist Pump Materials, Atrioventricular Shunts, Finger-joint Repair
Polyacrylonitrile	Membrane for dialysis
Polyurethane	Balloon for Intra-aortic Pump, Heart Valve Prostheses, Tubing
Polyvinyl alcohol	Drug-delivery System
Polyamides (nylons)	Nonabsorbable sutures, Tendon prosthesis, Drug-delivery System, Tracheal Tubes
Polyethylene terephthalate	Tendon and Ligament Reconstruction, Tracheal Replacement, Surgical Mesh Fabric
Hydrogels	Contact Lenses, Wound Dressings, Ophthalmic Implants, Drug-delivery System

7



## Types of Biomaterials

- Ceramics

- Ideal as biological implants; bone bonds well to them and they exhibit inertness within the body, high stiffness, and low friction and wear
  - Main drawback is their brittle nature, low impact resistance
- Used for restorative materials in dentistry



- Composites

- Combination of low density/weight and high strength make them ideal for prosthetic limbs

8



## Adhesives

- Adhesive biomaterials typically fall into three categories
  - Medical device assembly
    - Manufacturing of life-support equipment, sterile disposable items, sterile reusable items, and devices used for sensing, monitoring, and reporting
  - Hard tissue attachments
    - Orthopedics and dentistry
  - Soft tissue attachments
    - Wound closure




9

 TAINSTRUMENTS.COM

## Hydrogels

- Hydrogels are flexible polymers used in a wide variety of applications including tissue engineering, drug delivery, contact lenses and superabsorbent materials.
- When a hydrogel is heated, the structure dissociates and the gel “melts” or collapses. Understanding and quantifying the chemistry of gelation on cooling and collapse on heating is important for consistent hydrogel formation.
- Hydrogels sense changes of pH, temperature, and concentrations of metabolite and these properties are used in their applications.



10

 TAINSTRUMENTS.COM

## TA Instruments – Thermal Analysis, Calorimetry, Rheology, Mechanical Testing

- Heat of Reactions
- Phase Changes
- Transition Temperatures
- Reaction Kinetics
- Thermal, Oxidative Stability

Calorimetry



- Residual Solvent Content
- Sorption/Desorption Profiles
- Thermal, Oxidative Stability
- Decomposition Kinetics, Lifetime Plots

Thermogravimetric Analysis



- Dimensional Changes
- Coefficient of Thermal Expansion
- Softening Point
- Annealing Characteristics

Dilatometry/Thermo-Mechanical Analysis (TMA)



- Visco-Elastic Properties
- Structure-Property Relationship
- Process Conditions
- End Product Performance

Rheology/Rubber Testing Products



- Fatigue Life
- Failure Analysis
- Viscoelastic Properties
- Physiologic Simulation

Mechanical Testing



- Thermal Diffusivity
- Thermal Conductivity
- Specific Heat Capacity
- Thermal Resistivity

Thermal Conductivity




11

TA | TAINSTRUMENTS.COM

## What is Thermogravimetric Analysis (TGA)?

- Thermogravimetric Analysis (TGA) measures weight/mass change (loss or gain) and the rate of weight change as a function of temperature, time and atmosphere.
- Measurements are used primarily to determine the composition of materials and to predict their thermal stability. The technique can characterize materials that exhibit weight loss or gain due to sorption/desorption of volatiles, decomposition, oxidation and reduction.



12

TA | TAINSTRUMENTS.COM

## Thermogravimetric Analyzers

- A TGA must accurately:
  - control heating rate (furnace)
  - measure the change in temperature (thermocouple)
  - measure the mass of a sample and the change in mass as it is heated or held at an isothermal temperature (balance)



13

 TAINSTRUMENTS.COM


## TGA Applications

- Thermal stability of materials
- Oxidative stability of materials
- Composition of multi-component systems
- Decomposition mechanism when coupled with evolve gas analysis techniques (FTIR, MS)
- The effect of reactive or corrosive atmospheres on materials
- Moisture and volatiles content of materials
- Estimated lifetime of a product
- Decomposition kinetics of materials

14

 TAINSTRUMENTS.COM

## Thermal Stability of Polymers



## Thermal Decomposition and Degradation


- Solid polymeric materials undergo both physical and chemical changes when heat is applied. This usually results in undesirable changes to the properties of the material.
- ASTM provides a distinction to the terms thermal decomposition and thermal degradation:
  - **Thermal decomposition:** the process of extensive chemical species change caused by heat
  - **Thermal degradation:** a process whereby the action of heat or elevated temperature on a material, product or assembly causes a loss of physical, mechanical or electrical properties.



## Thermal Stability of Polymers

Upon degradation, typical property changes include:

- Reduced ductility and embrittlement
- Chalking
- Color changes
- Cracking
- General reduction in most other desirable physical properties




Thermal degradation generally involves changes to the molecular weight (and molecular weight distribution) of the polymer.

17




## TGA of Polyvinylpyrrolidone (PVP) - Residual Content (Aqueous Polymer Solution)



18




## 100-Fold Magnification



19



## 1,000-Fold Magnification Accurate analysis in the microgram range




20



## Kinetics by TGA

$$\text{Conversion: } \alpha = \frac{w_{in} - w_T}{w_{in} - w_f}$$




21

TAINSTRUMENTS.COM

## Polyactic Acid (PLA) Decomposition ASTM E1641 (Nitrogen Atmosphere)


TGA Analysis of PLA in Nitrogen  
2, 4, 6, 8 and 10 °C/min



22

TAINSTRUMENTS.COM

## PLA Decomposition by ASTM E1641 Lifetime Plot Using a 2% Conversion




If processing is performed at 200 °C (above the melting temperature, under a nitrogen blanket), PLA is predicted to lose 2% weight due to decomposition if held for 66 hours at that elevated temperature.

23



## Simple Heat Flux DSC Cell Schematic



24



## DSC Heat Flow

$$\frac{dH}{dt} = C_p \frac{dT}{dt} + f(T, t)$$

Heat Capacity


Glass Transition  
Specific Heat Capacity

Kinetic

Crystallization  
Cure reactions  
Volatilization  
Decomposition  
Denaturation

25



## Endothermic Heat Flow Heat Absorbed by Substance



### Endothermic Events

- Glass transition
- Melting
- Evaporation/volatilization
- Enthalpic recovery
- Polymorphic transitions
- Some decompositions

26



## Determining the Glass Transition Temperature, $T_g$


- ASTM Standards
  - E3418 Standard test method for transition temperatures and enthalpies of fusion and crystallization of polymers by differential scanning calorimetry
  - E1356 Standard test method for assignment of the glass transition temperatures by differential scanning calorimetry
- ISO Standard
  - ISO11357-2 Determination of glass transition temperature and glass transition step height

29



## The Glass Transition ( $T_g$ )


- The glass transition is a change in the free volume and molecular mobility in the amorphous phase of a material that results in a step change in heat capacity.



30



## Reporting the Glass Transition Temperature




## Miscibility of Polymer Blends – Impact on Tg

- **Miscible polymer blends** (homogeneous polymer blend): Polymer blend that is a single-phase structure. In this case, **one glass transition temperature** will be observed (Fox equation).

$$\frac{1}{T_g} = \frac{W_a}{T_{g,a}} + \frac{W_b}{T_{g,b}}$$

Fox Equation

- $T_{g,a}$ : glass transition of component a
- $T_{g,b}$ : glass transition of component b
- $w_a$ : weight fraction of component a
- $w_b$ : weight fraction of component b



Compound *	Tg (°C) - Measured	Tg (°C) – Fox Calculated
Trehalose	120	
70% Trehalose/30% Dextran	140	140
30% Trehalose/70% Dextran	182	178
Dextran	225	

32

\* Newman et al. J Pharm Sci. 2008, 97, 4840-4856



## What Affects the Glass Transition?

- Heating Rate
- Heating & Cooling
- Aging
- Molecular Weight
- Plasticizer (compatibility)
- Filler
- Crystalline Content
- Copolymers
- Side Chains
- Polymer Backbone
- Hydrogen Bonding

Anything that effects the mobility of the molecules, affects the Heat Capacity and, in turn, the Glass Transition

33



## Study of melting/crystallization using a DSC

- Melting is the process of converting solid, crystalline structure (lower energy) to a liquid amorphous structure (higher energy).
- Crystallization – The process of converting either solid amorphous structure (cold crystallization on heating) or liquid amorphous structure (cooling) to a more organized solid crystalline structure
- Melting:
  - low energy state → high energy state; requires input of energy; Endothermic peak
- Crystallization:
  - high energy state → low energy state; releases energy; Exothermic peak
- We integrate these peaks, on a time basis to determine the Heat of Fusion (melting) and Heat of crystallization

34



## Experimental Method

### A Heat Cool Ramp Method

- 1) Ramp 10 °C/min to -90 °C or Equilibrate to -90 °C\*
  - 2) Ramp 10 °C/min to 200 °C
  - 3) Ramp 10 °C/min to -90 °C or Equilibrate to -90 °C\*
  - 4) Ramp 10 °C/min to 200 °C
- Start test at least 30 °C below the expected T<sub>g</sub>
  - End test at least 50 °C above the expected T<sub>g</sub> for amorphous solids; stay below the decomposition temperature
  - Increase heating rate and/or mass if T<sub>g</sub> is barely detectable; this increases sensitivity

35



## Basis of Multiple Heat Cool Heat Cycles


- 1<sup>st</sup> heating scan: Called the initial run. It reveals information about the current condition of the specimen; i.e., the thermal and mechanical history as influenced by processing conditions, crystallinity and curing, service temperatures etc.
- Cooling scan: Subsequent controlled cooling rates create a new, known specimen history
- 2<sup>nd</sup> and 3<sup>rd</sup> heating scans: Used for determining the characteristic properties of the material. In the case of reactive resins, a third heat cycle may be performed to validate the completion of the reaction.

\* Reference: *Thermal Analysis of Plastics: Theory and Practice* by G.W. Ehrenstein, G. Riedel, P. Trawiel 2004


36





## Poly(lactic Acid) (PLA): Crystallinity After a Controlled Cooling Rate From the Melt



## Poly(lactic Acid) (PLA): Crystallinity After an Uncontrolled Cooling Rate From the Melt






## Crosslinking reactions – Cure of elastomers, adhesives and epoxies

- A “thermoset” is a cross-linked polymer formed by an irreversible exothermic chemical reaction.
- Crosslinking reactions are generally exothermic. As the chemical reaction takes place, it is almost always accompanied by a release of heat.
- The reactions can be easily monitored using a DSC
  - Heat of reaction
  - Residual cure
  - Glass transition
  - Heat capacity




## Cure Reaction of Liquid Silicone Elastomer Running Integral Analysis



42


## Cure Reaction of an Implantable Elastomer Kinetics Analysis



43



## Cure Reaction of an Implantable Elastomer Kinetics Analysis – Half Life Prediction




Temp °C	Temp K	1/Temp 1000/K	Half-Life min
95.00	368.15	2.72	450.0
101.00	374.15	2.67	56.8
107.00	380.15	2.63	7.65
113.00	386.15	2.59	1.10
119.00	392.15	2.55	0.167
125.00	398.15	2.51	0.0269
131.00	404.15	2.47	0.00458
137.00	410.15	2.44	0.000820

44




## Cure Reaction on a Rubber Process Analyzer (RPA)

- Measures **material response** to shear deformation or force as a function of **time, temperature, frequency, or deformation**
  - Typically reports viscoelastic properties of storage modulus ( $G'$ ), loss modulus ( $G''$ ), and tan delta
- Common Uses:
  - Complete pre and post cure viscoelastic characterization
  - Measure cure profile of materials during isothermal or anisothermal tests
  - Effects of filler/vulcanization network
    - Payne Effect
  - Identifying differences in viscoelastic material properties
    - Frequency dependence
    - Strain dependence
    - Stress relaxation




TA | TAINSTRUMENTS.COM

## Cure Reaction of an Implantable Silicone Elastomer Physical Property Measurement – Anisothermal



## Cure Reaction of an Implantable Silicone Elastomer Physical Property Measurement – Isothermal



## Thermal Analysis Review

- Thermogravimetric Analysis
  - Stability
  - Volatile content
  - Decomposition kinetics
- Differential Scanning Calorimetry
  - Transition temperatures
  - Crystallinity
  - Heat of reactions
  - Reaction kinetics

48



## Agenda

- 9 - 10:15am: Thermal analysis and stability of biomaterials
- **10:15-10:30am: Break**
- 10:30-noon: Using rheology to characterize flow and viscoelastic properties of hydrogels, adhesives and biopolymers
- 12-1pm: Lunch
- 1-2:15pm: Mechanical testing of medical devices
- 2:15-2:30 - break
- 2:30-3:15pm: Mechanical testing of engineered tissues and biomaterials
- 3:15-4pm – Q&A with TA Instruments Applications Engineers

49



## Thank You

The World Leader in Thermal Analysis,  
Rheology, and Microcalorimetry



50