





Dynamic Mechanical Analysis Basic Theory & Applications Training Day 2



Course outline

Part 1

- Basic Theories of Dynamic Mechanical Analysis
- DMA Instrumentation and Clamps
- Introduction to DMA Experiments
 - o Dynamic tests
 - o Transient tests
- Appendix: Screenshots from the instrument control software

Part 2

- Recap of Part 1
- DMA Applications and data interpretation
- Troubleshooting experimental Issues
- Time-temperature superposition (TTS)



Recap of Day 1



Recap: Viscoelasticity and mechanical tests



Recap: DMA instrumentation

RSA G2

Discovery DMA850

Electroforce series (high load frame, fatigue)









850/800 clamps

RSA-G2 clamps

Film/Fiber Tension





Submersible Tension



Submersible Bending





Submersible Compression



3-Point Bending

Shear Sandwich





Recap: DMA dynamic tests

Strain Sweep



Time sweep

Stress or strain



Frequency Sweep



Temperature ramp/sweep



ΤA

Recap: DMA transient tests



DMA Applications and Data Interpretation



The DMA results can correlate to...





Most common DMA test – temperature ramp

- The most Common DMA measurement is a dynamic temperature ramp
- The test results report modulus (E*, E', E''), damping factor (Tan δ) and transition temperatures (T_g)
- Provide information to polymer's structure-property relationship





What is Glass Transition (T_g)?

- A transition over a <u>range of temperature</u> from a glassy state to a rubber state in an amorphous material
- Mechanical:
 - Below the Glass Transition, the material is in a brittle, glassy state, with a modulus of 10⁹ Pa
 - Above the Glass Transition, the material becomes soft and flexible, and the modulus decreases two to three decades
- Molecular:
 - Below the Glass Transition, polymer chains are locked in place, without sufficient energy to overcome the barrier for rotational or translational motion.
 - At temperatures above the Glass Transition, there is molecular mobility, and chains can slide past each other



How to Define T_g in DMA?



Secondary Transitions

<u>Glass Transition</u> (T_g)

14

 Cooperative motion among a large number of chain segments, including those from neighboring polymer chains

<u>Secondary Transitions</u> (T_{β}, T_{γ})

- Local Main-Chain Motion intra-molecular rotational motion of main chain segments four to six atoms in length
- Side group motion with some cooperative motion from the main chain
- Internal motion within a side group without interference from side group.
- Motion of, or within, a small molecule or diluent dissolved in the polymer (eg. plasticizer.)



Primary and Secondary Transitions in PC



What Will Affect T_g and Modulus?

- Heating rate
 - Thermal lag
- Test frequency
- Polymer structures
 - Rigid polymer chain shows higher T_g (e.g. PS)
 - Flexible polymer chain shows lower T_g (e.g. PE)
 - Crystallization
 - Degree of crosslinking



The Effect of Test Frequency on T_g

- The glass transition is strongly influenced by the frequency of the test. The T_g is a molecular relaxation that involves cooperative segmental motion.
- Because the *RATE* of segmental motion depends on temperature, as the frequency increases, the relaxations associated with the T_g (where $\tau \sim t$) can only happen at higher temperatures (higher temperatures corresponds to higher kinetic energy that is imparted to the system).
- In general, increasing the frequency will
 - Increase the T_g

17

- Broaden the peak
- Decrease the intensity of tan δ or loss modulus
- Decrease the slope of the storage modulus curve in the region of the transition.



PET Film: Effect of Frequency on T_g



 $\tan \delta =$ $\frac{E}{E'}$



Time-period of deformation is short (high frequency) Increased elastic response Higher E', lower E'' Lower tan δ

Time-period of deformation is long (low frequency) Increased viscous response Higher E", lower E' Higher tan δ



Effect of Crystallinity in Polymer

"The major effect of the crystallite in a sample is to act as a crosslink in the polymer matrix. This makes the polymer behave as though it was a crosslinked network, but as the crystallite anchoring points are thermally labile, they disintegrate as the temperature approaches the melting temperature, and the material undergoes a progressive change in structure until beyond T_m, when it is molten"





Random Chain 100% Amorphous

Fringed Micell Crystalline

Cowie, J.M.G., Polymers: Chemistry & Physics of Modern Materials, 2nd Edition, 1991p. 330-332. ISBN 0 7514 0134 X



Effect of crystallinity on modulus and tan $\delta\text{:}$ PET



20 TA Applications note: RH100 Measurement of Glass Transition Temperatures by Dynamic Mechanical Analysis and Rheology https://www.tainstruments.com/pdf/literature/RH100.pdf

Polymer crosslinking

- Linear polymers can be chemically or physically joined at points to other chains along their length to create a crosslinked structure
- Chemically crosslinked systems are typically known as **thermosetting polymers** because the crosslinking agent is heat activated.
- Thermosetting polymers typically decompose before they can be molten.





Effect of crosslinking density on Tg



Increasing crosslinking density:

- Tg shifted to higher temperature
- Transition becomes broader and weaker (tan δ decreases)
- Rubbery plateau modulus increases





Quantification of crosslinking density



Crosslinking density:

$$M_c = \frac{3RTd}{E'_{rubbery}}$$

M_c is the molecular weight between crosslinks
R is the universal gas constant
T is the absolute temperature (in K)
d is the density of the polymer.

A lower Mc implies higher extent of crosslinking



TA Applications note RH102: Quantifying Polymer Crosslinking Density Using Rheology and DMA, https://www.tainstruments.com/pdf/literature/RH102.pdf

Applications of DMA in Polymer blends





Polymer Blend – Miscible Case

Polymer A Polymer Blend: A + B Polymer B Polymer Blend Elastic Modulus A + B 100 % Polymer A 100% **Polymer B** Temperature (°C)

Aerospace Coating



Polymer Blend – Miscible Case



26

Polymer Blend - Immiscible Case





Using Glass Transition to Evaluate Blending

Insufficient blending on the manufactured batch results in two transitions



Anisotropic Materials

- Anisotropic materials have different properties in different directions.
 - fibers, wood, fiber-filled composites
 - oriented amorphous polymers, injection molded specimens
 - crystalline polymers in which the crystalline phase is not randomly oriented
- Show different moduli at different measurement direction

Nielsen, Lawrence E., Mechanical Properties of Polymers and Composites, Marcel Dekker, Inc., New York, 1974, pp. 39-40.



DMA of Polyester/Glass Fiber Reinforced Composite





DMA of Polyester/Glass Fiber Reinforced Composite





Oriented Polymer: Shrink Wrap





- "Shrink Wrap" is widely used in packaging to conveniently label irregularly shaped containers
- When this kind of biaxially oriented film is heated, it turns to shrink and returns to lower energy random-coil state.
- PET-based Label
 - Z-direction
 - X-direction
- RSA G2: 5° C per minute, 1 Hz



PET Film Measured at X and Z Direction



33

Iso-force Temp Ramp: measure shrinkage





Iso-strain Temp Ramp: measure shrinking force





Effect of filler on modulus and glass transition



Figure 6: An epoxy sample with 10% (A) and 30 % (B) talc filler. The glass transitions are nearly identical on all signals.

TA Applications note: RH100 Measurement of Glass Transition Temperatures by Dynamic Mechanical Analysis and Rheology https://www.tainstruments.com/pdf/literature/RH100.pdf


Effect of Aging on Elastomer O Rings



Effect of Plasticizer

- *Plasticizers* are generally low molecular weight organic additives which are used to soften rigid polymers
- Typically added to a polymer for two reasons:
 - 1. To lower the T_g to make a rigid polymer become soft
 - 2. To make the polymer easier to process
- *Plasticizers* act as lubricant in polymer chains
- Therefore plasticizers will have the effect of:
 - 1. Lowering the $T_{\!_g}$ and broadening tan δ peak
 - 2. Lowering the moduli (E' and E")



Effect of Plasticizer on Vinyl Flooring





Effect of Aging on T_g of PVC Roofing Membrane





Humidity Influence: Nylon Film



ΤA

EPDM immersed in oil is subject to significant plasticizing



O-rings: Stress Relaxation

• Squeeze the O-ring to a certain strain. Hold it constant, then measure how long it takes for the force to relax



O-ring Stress Relaxation





Creep Tests on Packaging Bags

• Apply a load force, then measure how much the bag is stretched







Memory Foam: Multi-Step Creep Recovery



Multi-Specimen Tension Fatigue BR Rubber (Tire Rubber) Example





6-Specimen Fixture with 6 Force Sensors



Multi-Specimen Tension Fatigue BR Rubber (Tire Rubber) Example

- Many SN points obtained in relatively short time
 - 4 test series (24 failures) in ~1 month
- Data is fit to exponential equation (line on plot)
- Data and exponential fit is used to compare materials





Troubleshooting Experimental Issues



Factors that influence measurement results



measurements

- Force
- Displacement
- Temperature

Check list

- Instrument
 calibrations
- Thermocouple position

Geometry

- Clamp
- Sample size
- Sample loading
- Sample integrity

<u>Check list</u>

- Clamp calibrations
- Proper loading torque/force
- Sample uniformity and thermal stability
- Bad contact or sample bending and twisting

Test Parameters

Check list

- Loading force
- Proper force track
- Equilibration time
- Heating rates
- Amplitude
- Frequency



Instrument and clamp

- Perform clamp calibration every time when attaching a new clamp on the instrument
 - http://www.youtube.com/user/TATechTips
- Perform confidence verification check using polycarbonate provided by TA Instruments (see next slide)





DMA confidence check - polycarbonate

- Load Polycarbonate (L \approx 17.5, w \approx 12.85, t \approx 1.6mm)
- Use Single Cantilever Clamp
 - 20-30 micrometer amplitude
 - 1 Hz frequency
- Storage Modulus at Room Temperature
 E' = 2.35 GPa (2350 MPa) +/- 5%
- Tan Delta at Room Temperature Tan δ < 0.01
- Transition Temperature

Tan δ peak between 155-160°C @ 1Hz, 3-5°C/min E" peak will be about 5°C lower

p/n: 982165.903





Temp ramp on polycarbonate

• Available from TA for Instrument verification



PC sample



p/n: 982165.903

Clamp: single cantilever Temperature: ambient to 180°C Heating rate: 3°C/min Frequency: 1 Hz Amplitude: 20 µm



Sample preparation

- DMA measurement results are highly depending on the quality of sample preparation
- Sample flatness and uniformity have big influence to the accuracy of the test results



E' Increase in a strain sweep

Why does E' increase with increasing amplitude? How to verify linear viscoelastic region?





Instrument: Q800 Clamp: 3-p bending Temperature: ambient Amplitude: 1-500 µm Frequency: 1Hz



E' Increase in a strain sweep

The sample is not flat and not in full contact with the clamp face. Solutions: (1) Prepare a flat sample (2) Increase force track or increase static force

> Sample: ABS strain sweep File: T:...\ABS STRAIN SWEEP.003 DMA Size: 50.0000 x 12.9100 x 3.1700 mm Operator: TC Method: Strain Sweep Run Date: 23-Jan-2018 14:17 Instrument: DMA Q800 V21.3 Build 96 10000 9000 8000 7000 6000 5000 [•] Storage Modulus (MPa) [] Dynamic Force (N) 4000 3000 0.1 2000 1000 -+0.010.1 10 1Ò0 1000 Universal V4.5A TA Instruments Amplitude (µm)







E' increase in a temp ramp

• The E' of a material should not increase with temperature unless it is crystalized or crosslinked





Noisy modulus after Tg

What is the problem with this data collected after T_g ?





Instrument: Q800 Clamp: tension Temperature: -100°C to 150°C Heating rate: 3°C/min Frequency: 1Hz Amplitude: 10 μm



Noisy modulus after Tg

What is the problem with this data collected after T_g ?





Instrument: Q800 Clamp: tension Temperature: -100°C to 150°C Heating rate: 3°C/min Frequency: 1Hz Amplitude: 10 μm

Sample sagging

• Sample sagging after T_g

• Solution: use cantilever clamp instead of 3-p bending





Instrument: RSA G2 Clamp: 3-p bending Temperature: 50°C to 180°C Heating rate: 3°C/min Frequency: 1Hz Amplitude: 10 μm



Compression testing on foam/rubber

Question: How to measure correct modulus? Answer: (1) Prepare regular sample

(2) Apply appropriate static force



Time-Temperature Superposition (TTS)



Time and temperature relationship



- Some materials show a time dependence that is proportional to the temperature dependence
 - Decreasing temperature has the same effect on viscoelastic properties as increasing the frequency, and vice versa
- For such materials, changes in temperature can be used to "re-scale" time, and predict behavior over time scales not easily measured

Dealy, J., Plazek, D., Time-Temperature Superposition – A Users Guide, Rheology Bulletin, 78(20) 16 (2009)



Merits of TTS

- TTS can be used to extend the frequency beyond the instrument's range.
- Low frequency data predicts material behavior over longer timescales that cannot be practically measured on a laboratory instrument.
 - Creep or Stress Relaxation TTS can predict behavior over longer times under static load/deformation.
- High frequency data predicts material behavior at short timescales (high-speed impact, mechanical vibrations, acoustics) that are challenging to measure accurately using an analytical instrument.
 - High frequency mechanical testers almost always are influenced by inertial contributions from the instrument itself.

- Dealy J, Plazek D. Time-Temperature Superposition A Users Guide. Rheology Bulletin 2009; 78: 16
- van Gurp M, Palmen J. *Time-temperature superposition for polymeric blends.* Rheology Bulletin 1998; 67: 5.
- Cox WP and Merz EH. Correlation of dynamic and steady flow viscosities. Journal of Polymer Science 1958; 28: 619 🔍



1. TTS typically works with polymers that are thermo-rheological simple

2. TTS may or may not work if:

a.Polymer contains multiple phases: composite or blendsb.Filled polymers

- 3. TTS does not work if:
 - a. Sample properties change (e.g. crystalize, melt, cure, decompose) within temperature of interest



Temperature Step & Hold- Single /Multi-Frequency

• A step and hold temperature profile is applied. The material response is monitored over a range of frequencies, at constant amplitude of deformation





Master curve from TTS Shifting, polystyrene



Principle of time-temperature superposition (TTS)





Method development for TTS studies



Selecting a reference temperature

- Decide on the Reference Temperature: T_r.
- Higher frequencies = shorter times. Test temperatures $< T_r$.
- Lower frequencies = longer times. Test temperatures > T_r .



Importance of thermal analysis in TTS studies

70

- Sample properties change (e.g. crystalize, melt, cure, decompose) within temperature of interest
- Scan material over temperature range to get an idea of transition behavior and modulus-temperature.



Importance of thermal analysis in TTS studies

71

- Scan material over temperature range to get an idea of transition behavior and modulus-temperature.
- Allows for optimizing experimental method (axial force, force track, % strain, etc.) prior to longer TTS experiments.



TTS example: shifting of individual curves, polystyrene








Refer to http://www.tainstruments.com/tts-and-trios-software/







77







Master curve from TTS Shifting, polystyrene





How accurate is the master curve?





- 1. Van-Gurp Palmen plot should be a smooth curve
- 2. Cole-Cole plot should be a smooth curve
- 3. Shift factor (a_T) should fit the WLF or Arrhenius equation
- 4. Same a_T should be applicable to all viscoelastic parameters (E', E", tan δ)

If any one of the above criteria are not met, TTS may not be applicable for your material.



1. van Gurp-Palmen plot to Validate TTS



Angular frequency ω (rad/s)

1. van Gurp-Palmen plot to Validate TTS

• The van Gurp-Palmen plot can be directly plotted in TRIOS to validate the application of (TTS).







2. Cole-Cole plot to Validate TTS

• The Cole-Cole plot can be directly plotted in TRIOS to validate the application of (TTS).



Storage modulus G' (Pa)

3. Shift factors should fit WLF or Arrhenius Equations

• Master Curves can be generated using shift factors derived from the Williams, Landel, Ferry (WLF) model

 $\log a_T = -c_1(T-T_0)/(c_2+(T-T_0))$

 a_T = temperature shift factor T_0 = reference temperature

 $c_1 \& c_2$ = constants from curve fitting

Generally, c_1 =17.44 & c_2 =51.6 when $T_0 = T_g$

- The Arrhenius model works better if
 - $T > T_g + 100^{\circ}C$; or $T < T_g$ and polymer is not elastomeric
 - temperature range is small, then $c_1 \& c_2$ cannot be calculated precisely

 $\ln a_{T} = (E_{a}/R)(1/T-1/T_{0})$

- a_T = temperature shift factor E_a = Apparent activation energy
- T_0 = reference temperature R = gas constant

3. Shift Factors a_T vs. Temperature



88

4. Same a_{τ} should be applicable to all viscoelastic parameters (E', E", tan δ)



Limitations of TTS

- TTS typically works with polymers that are thermo-rheological simple
- TTS may or may not work if:
 - Polymer contains multiple phases: composite or blends
 - Filled polymers
- TTS does not work if:
 - Sample is partially crystalized, crosslinked, or highly filled
 - Sample properties change (e.g. crystalize, melt, cure, decompose) within temperature of interest



Limitations of TTS

- TTS typically works with polymers that are thermo-rheological simple
- TTS may or may not work if:
 - Polymer contains multiple phases: composite or blends
 - Filled polymers trans polyisoprene (PIP) with 80% inorganic filler
- TTS does not work if:
 - Sample is partially crystalized, crosslinked, or highly filled
 - Sample properties change (e.g. crystalize, melt, cure, decompose) within temperature of interest



Temperature sweep Trans PIP





Instrument: RSA-G2 DMA Fixture: 3-point bending Purge gas: Nitrogen



Van Gurp plot PIP

Van-Gurp plot shows discontinuity indicating TTS might not be applicable to this sample.





Master curve PIP

Same a_{T} cannot be applied to all viscoelastic parameters to obtain a continuous master curve.



Frequency f (Hz)

Setting up TTS procedure on DMA850

Sample: Example TTS procedure					
× Clamp: Film Clamp					
Scillation V Temperature St	weep (Multif	requency) Y	হি 🐸 🛃		
Strain Initial/preload force I Use Force Track	0.01 0.5 125.0	% N %			
Sweep from Temperature increment Soak time	140 5 00:05:00	°C to 175 °C hh:mm:ss] °C		
Sweep Mode Logarithmic O Linea Frequency Points per decade	or O Discr 0.1 5	ete Hz to 10.0	Hz		
Test Settings	Post	Test Conditions			



Setting up TTS procedure on Q800

Summary Procedure Notes	Frequency Table
Summary Procedure Notes Procedure Information Test Temp Step / Freq Sweep Notes Material is exposed to a series of increasing isothermal temperatures. At each temperature, the material is deformed at a constant amplitude (strain) over one or more frequencies and the mechanical properties	C Single C Log C Linear C Discrete Frequency 100.00 to 0.10 Hz Points per decade: 5
Method Amplitude : 15.0000 µm Strain : 0.0000 % Amplitude within the linear region	Frequency A Refresh <u>I</u> able 1 100.00 Image: Second
Start temperature:	
Final temperature: 150.00 *C Temperature increment: 2.50 *C Isothermal soak time 5.00 min Method (Frequency Table/ *C	# Running Segment Description 1 Image: Data storage Off 2 \$↑ Equilibrate at 35.00 °C 3 \$↑ Equilibrate at 35.00 °C 3 \$↑ Isothermal for 5.00 min 4 Hz Frequency sweep 5 \$↑ Increment by 2.50 °C
	 Isothermal for 5.00 min Hz Frequency sweep Repeat segment 5 until 150.00 °C



Setting up TTS procedure on RSA G2

👽 Sampl	e: PET film LN2 only				
👽 Geom	etry: Tension fixture (rec	tangle)			
Proces	dure of 2 steps				\$ 8 % #
♥ 1♦ 2	: Conditioning Options A : Oscillation Temperature	ctive, Enabled 9 Sweep			
	Environmental Control – Start temperature Soak time End temperature Temperature step Step soak time	-100 300.0 200 10 300.0]°C]s]°C]°C]s	Inherit	
	Test Parameters Strain % Logarithmic sweep	0.02	%	▼	
	Frequency Points per decade	0.1 5	to	10.0 Hz	

Data acquisition

Advanced



- Time-temperature superposition (TTS) is a valuable tool for describing the viscoelastic properties of polymeric materials over a wide range of times/frequencies.
- Thermal analysis techniques such as DSC and DMA can be valuable tools for developing the method for TTS experiments.
- The modulus vs frequency curves at various temperatures can be shifted horizontally and/or vertically depending on whether the material is thermo-rheologically simple.
- Experimental verification of part of a master curve provides additional confidence in the applicability of the TTS principles to a material.



References for TTS

- 1. Ferry JD. *Viscoelastic Properties of Polymers.* Wiley 1970, Chapter 11.
- 2. Dealy J, Plazek D. *Time-Temperature Superposition A Users Guide.* Rheology Bulletin 2009; 78: 16
- 3. van Gurp M, Palmen J. *Time-temperature superposition for polymeric blends.* Rheology Bulletin 1998; 67: 5.
- 4. Cox WP and Merz EH. *Correlation of dynamic and steady flow viscosities.* Journal of Polymer Science 1958; 28: 619

- 5. Tajvidi M, Falk RH, Hermanson JC. *Time-temperature superposition principe applied to a kenaf-fiber/high density polyethylene composite*. J. Appl. Polym. Sci. 2005; 97: 1995
- 6. Lovering EG, Wooden DC. *Transitions in Trans-1,4-Polyisoprene.* J. Polym. Sci. 1969; 7: 1639
- Mavridis H, Shroff RV. *Temperature dependence of Polyolefin Melt Rheology*. Poly. Eng & Sci 1992; 32: 1778



Summary





Summary: DMA applications

- DMA Applications and data interpretation
 - Glass transition temperature
 - Definition
 - Factors influencing Tg
 - Determination of cross-linking density
 - Evaluation of polymer blends and compatibility
 - Effect of anisotropy on mechanical properties
 - Influence of additives
 - Fillers
 - Plasticizers
 - Application of creep-recovery and stress relaxation tests



Summary: Troubleshooting DMA experiments

- Troubleshooting common experimental issues
 - Confidence check with polycarbonate
 - Sample preparation
 - Plotting the stiffness signal
 - Tips for compression testing

Summary: Time-temperature superposition (TTS)

- Time-Temperature Superposition (TTS)
 - Introduction to time-temperature superposition (TTS)
 - Performing TTS in TRIOS
 - Thermal analysis techniques to aid method development in TTS
 - Does TTS apply to my sample?
 - Van-Gurp Palmen and Cole-Cole plot
 - WLF and Arrhenius equation
 - Fitting multiple viscoelastic parameters
 - Experimental verification
 - Limitations of TTS



Getting Started Manuals on your desktop

TA



TA Instruments

Discovery DMA 850 Manuals

To view the desired manual using Acrobat Reader, click the name in the list below:

TA Manual Supplement (Contains important information applicable to all manuals.)

Instrument Documentation

Discovery DMA 850 Getting Started Guide

Accessory Documentation

Air Chiller System (ACS) Getting Started Guide DMA-RH Accessory Getting Started Guide Gas Cooling Accessory (GCS) Getting Started Guide Nitrogen Purge Cooler (NPC) Getting Started Guide

Software Documentation

What's New in TRIOS Software Installing TRIOS Software

Issued February 2020

Site Preparation Guides and Installation Requirements Discovery DMA 850 Site Preparation Guide

Additional Information DMA Clamping Factors for Compression Clamps



TA Instruments

RSA-G2 Manuals

To view the desired manual using Acrobat Reader, click the name in the list below:

TA Manual Supplement (Contains important information applicable to all manuals.)

Instrument & Accessory Documentation

RSA-G2 Getting Started Guide RSA-G2 FCO Camera Kit Installation Guide RSA-G2 LN2 Kit Installation Guide RSA-G2 Chiller Panel Kit Installation Instructions RSA-G2 Dielectric Accessory Getting Started Guide ACS Getting Started Guide - UPDATED

Software Documentation What's New in TRIOS Software Installing TRIOS Software Configuring a New Geometry in TRIOS Software

Miscellaneous Documentation RSA-G2 Site Preparation Guide





Guide

RMX File Utilities

To view the desired manual using Acrobat Reader, click on the name in the list below: TA Manual Supplement (Contains important information applicable to all manuals)

Instrument & Accessory Manuals

Software Manuals

Q Series[™] Instrument Control Getting Started Guide

Advantage Integrity[™] Getting Started

Specialty Library Getting Started Guide

Universal Analysis Getting Started DSC Q Series[™] Getting Started Guide Guide

RCS Getting Started Guide LNCS Getting Started Guide

Tzero[®] PDSC Getting Started Guide

PCA Getting Started Guide DSC Pressure Cell Getting Started Guide

DSC High Pressure Capsule Kit DSC High Volume Pan Kit

DSC Circulator-Based Cooling System TGA Q5000 IR Getting Started Guide

Miscellaneous Documents Installing/Updating Advantage[™]

TGA Q Series[™] Getting Started Guide Updating Q Series[™] Instrument

Software

DMA Q Series[™] Getting Started Guide New Features in Advantage Q Series[™] GCA Getting Started Guide DMA Humidity Accessory Getting Started New Features in

Advantage Integrity[™] SDT Q Series[™] Getting Started Guide

TA Update TMA O Series[™] Getting Started Guide MCA Getting Started Guide

MCA70 Getting Started Guide

TÀ

TGA Hi-Res[™] Option

Q5000 SA Getting Started Guide



Trios and Advantage Help











Instructional Videos

From <u>www.tainstruments.com</u> click on Videos, Support or Training

🔀 info@tainstrume	nts.com	\$ 302-427-4000		👙 💠 🤨 Language selector		Training Y News & Ever
Instruments	Product	s ~ About TA Instrume	nts ~ Videos ~ Su	pport ~ Training ~	News & Events ~ Caree	Initial Installation & Training
		Service Support	Application Support	Soft vare Downloads & Support	Support Plans	Course Schedule
		Service Support Helpline	Applications Support		Lifetime Support Plan	
	in	Installation	Helpline Took Tipo	Downloads	Premium Support Plan	Training Courses >
No. of the		Requirements & Repairs Tech Tips	recirrips	software	Plus Support Plan	
	6	The IQ/OQ Product Applications No Offering Library	Applications Notes Library	Software Sorted by	Basic Support Plan	ElectroForce Training
	Calibration with Certified Standards	Training	Instruments Report a Bug	Performance Maintenance Visit (PMV)		
		Safety Data Sheets		Request a Feature	Academic Support Plan	Strategies for Better Data
	- U	Supported Instruments				T 11 540
						Training FAQ

Select Videos for TA Tech Tips, Webinars and Quick Start Courses



See also: <u>https://www.youtube.com/user/TATechTips</u>



Instructional Video Resources

Quickstart e-Training Courses

Web based e-Training Courses

TA Instruments offers a variety of training opportunities via the Internet. e-Training opportunities include the following:

QUICKSTART e-TRAINING COURSES

QuickStart e-Training courses are designed to teach a new user how to set up and run samples on their analyzers. These 60-90 minute courses are available whenever you are. These pre-recorded courses are available to anyone at no charge. Typically these courses should be attended shortly after installation.

Contact Us for Web based e-Training Courses

https://www.tainstruments.com/videos/quick-start-guides/



- For additional questions:
- Email: rheologysupport@tainstruments.com
- Please put Online Training Questions in the subject line
- Link to download the presentation:

https://www.tainstruments.com/online-training-coursedownloads/


Thank you for your time!

