

## Newsletter At-A-Glance:

- Re-introducing ElectroForce
- New System Status Indicator Lights
- Upcoming WinTest Training Session

## ElectroForce Linear Motors: The Heart of Our Test Instruments

*New video highlights the superior performance of ElectroForce test systems*

ElectroForce test instruments integrate a proprietary, friction-free motor design that offers quiet, energy-efficient performance over wide ranges of force and frequency. Our new video re-introduces the technology at the heart of our test instruments which are used every day by some of the world's largest and most prestigious institutions.




[Click Here](#) to watch the video.

[Click Here](#) to download our latest ElectroForce brochure.

## Introducing New System Status Indicator Lights

*An easy visual cue provides the status of your system*

Many of our customers perform long-term tests that run for days, weeks or even months! To provide a visual cue so that users can more easily determine the status of a test, we are introducing System Status Indicator (SSI) Lights. SSI lights will indicate whether the test instrument is powered up, whether a test is running, or when an error state has been reached. SSI Lights can be added to any ElectroForce test instrument that is running WinTest 7.1 software or newer.




To learn more about SSI Lights, [Click Here](#).

## Complimentary WinTest Software Training

*Register for our next online training session*

Perhaps you've used WinTest in the past and are in need of a refresher or you are a new user and would like more training to gain a better understanding. As part of our Above & Beyond™ Support, we will be hosting an online training session on Monday, December 7th, that will cover the essential functions within WinTest. There will also be time to answer any questions that you might have. Sign up now as registration is limited.




For more information or to sign up, [Click Here](#).

To view our library of WinTest training videos, [Click Here](#).

## CURRENT PROMOTIONS

### Buy One, Get One



*Purchase a new ElectroForce 3200, 3330 or 3510 test instrument, and receive credit towards the purchase of system accessories as well as free Dynamic Mechanical Analysis (DMA) software!*

[More Info](#)

## ElectroForce Access Free Accessory



*Every 5500 test instrument purchased through the ElectroForce Access program now includes a free accessory for a limited time.*

[More Info](#)

## Academic Matching Grant Program



*Announcing the addition of ElectroForce test instruments to the TA Instruments Academic Matching Grant (AMG) Program.*

[More Info](#)

