



TA Instruments *New Features in* *TAM Assistant™ Software*

Notice

The material contained in this manual, and in the online help for the software used to support TA Instruments products, is believed adequate for the intended use of the instruments. If the instruments or procedures are used for purposes other than those specified herein, confirmation of their suitability must be obtained from TA Instruments. Otherwise, TA Instruments does not guarantee any results and assumes no obligation or liability. TA Instruments also reserves the right to revise this document and to make changes without notice.

TA Instruments may have patents, patent applications, trademarks, copyrights, or other intellectual property covering subject matter in this document. Except as expressly provided in written license agreement from TA Instrument, the furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property.

TA Instruments Operating Software, as well as Module, Data Analysis, and Utility Software and their associated manuals and online help, are proprietary and copyrighted by TA Instruments. Purchasers are granted a license to use these software programs on the module and controller with which they were purchased. These programs may not be duplicated by the purchaser without the prior written consent of TA Instruments. Each licensed program shall remain the exclusive property of TA Instruments, and no rights or licenses are granted to the purchaser other than as specified above.

Software Upgrade Access

The individual instrument software described in this manual is shipped with all appropriate shipments of TAM instruments. For all previous versions of the software provided with these instruments, upgrades are available via on-line downloads from the TA Instruments website. Contact TA Instruments customer service or your TA Instruments account manager for details.

Important: TA Instruments Manual Supplement

Please click the [TA Manual Supplement](#) link to access the following important information supplemental to this Getting Started Guide:

- TA Instruments Trademarks
- TA Instruments Patents
- Other Trademarks
- TA Instruments End-User License Agreement
- TA Instruments Offices

Table of Contents

Notice	1
Software Upgrade Access	1
Important: TA Instruments Manual Supplement	2
What's New for TAM Assistant v3.0.4	4

What's New for TAM Assistant v3.0.4

Note: not yet available for TAM III customers, only TAM IV and TAM Air are supported

- Communication infrastructure is improved
- TAM IV firmware is updated and communication stability significantly improved
- Hot-swap possibility of TAM IV accessories
- Simplified firmware updates
- New office inspired ribbon user interface
- Dashboard device view for better overview
- Run sequences are significantly improved
- Updated experimental wizards
- Calibration number dropdowns for better overview
- Various bug fixes and improvements

TA Instruments Offices

For information on our latest products, contact information, and more, see our web site at:
<http://www.tainstruments.com>

TA Instruments — Waters LLC
Corporate Headquarters
159 Lukens Drive
New Castle, DE 19720
USA

Telephone: 302-427-4000
Fax: 302-427-4001
Email: info@tainstruments.com